ACCENT-200 ALBUMIN

Nr kat. 7-238 (PL)

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(EN)

INTENDED USE

Diagnostic kit for determination of albumin concentration intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Albumin is the major serum protein, but is present also in other body fluids: cerebrospinal, pleural and peritoneal. Albumin regulates blood oncotic pressure and serves as amino acids reservoir. Beyond of these functions albumin is very important transport protein – binds and keeps dispersed bilirubin, hormones, vitamins, calcium, magnesium, fatty acids and medicines. Decreased albumin blood level is caused usually by liver or kidney disease, malabsorption or malnutrition.

METHOD PRINCIPLE

Bromocresol green (BCG) forms with albumin, in succinate buffer (acid medium), a coloured complex. The colour intensity of the formed complex measured at 630 nm is proportional to albumin concentration in the sample.

REAGENTS

Package

1-REAGENT 1 x 34 ml

The reagent, stored at 2-8°C is stable up to the kit expiry date printed on the package. The reagent stored on board of the analyser at 2-10°C is stable for 9 weeks (ACCENT-200 II GEN) or 12 weeks (ACCENT MC240, ACCENT S120).

Concentrations in the test

succinate buffer ≤ 108 mmol/l
bromocresol green (BCG) ≤ 0.35 mmol/l
sodium hydroxide ≤ 60 mmol/l
detergent
preservative

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not freeze the reagent.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meets the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Warning

H319 Causes serious eye irritation.
P280 Wear protective gloves, protective clothing, eye protection or face protection.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

SPECIMEN

Serum free from hemolysis.
Serum should be separated from red blood cells as soon as possible after blood collection.
Serum can be stored up to 3 days at 2-8°C or 6 months at -20°C.² Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT is ready to use.
Deionised water is recommended as a reagent blank.

Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120, there is a probability of **cross-contamination** affecting the tests results: ALBUMIN - UA, ALBUMIN - UA PLUS, ALBUMIN-MICROALBUMIN. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

REFERENCE VALUES ⁷

serum		g/dl	g/l
children	3 – 5 years	3.9 – 5.0	39 - 50
	6 – 15 years	4.1 – 5.1	41 - 51
females	16 – 54 years	3.9 – 5.0	39 – 50
males	16 – 29 years	4.6 – 5.3	46 - 53
	30 – 54 years	4.4 – 5.1	44 - 51
females and males	55 – 79 years	4.2 – 5.0	42 – 50

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples. For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, BS-120, systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended.

For the calibration of automatic analysers: ACCENT S120, ACCENT MC240, ACCENT M320 systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionized water is recommended as a calibrator 0.

The calibration curve should be prepared every 8 weeks (ACCENT-200 II GEN) or every 12 weeks (ACCENT MC240, ACCENT S120), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers: ACCENT-200 and/or ACCENT-200 II GEN and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

Sensitivity

1.14 g/dl (11.4 g/l) - ACCENT-200
0.70 g/dl (7 g/l) – ACCENT MC240

Linearity

up to 6.5 g/dl (65 g/l) - ACCENT-200
up to 6.1 g/dl (61 g/l) – ACCENT MC240

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 2.5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1200 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [g/dl]	SD [g/dl]	CV [%]
ACCENT-200 II GEN n=20	level 1	2.80	0.02	0.85
	level 2	4.40	0.03	0.57
ACCENT MC240 n=20	level 1	2.88	0.02	0.67
	level 2	4.03	0.03	0.70
Reproducibility (day to day)		Mean [g/dl]	SD [g/dl]	CV [%]
ACCENT-200 II GEN n=80	level 1	2.82	0.05	1.87
	level 2	4.39	0.07	1.60
ACCENT MC240 n=80	level 1	4.2	0.09	2.2
	level 2	3.0	0.07	2.4

Method comparison

A comparison between albumin values determined at ACCENT-200 (y) and at ADVIA 1650 (x) using 21 samples gave following results:

$y = 0.9908x + 0.1259$ g/dl;
 $R = 0.998$ (R – correlation coefficient)

A comparison between albumin values determined at ACCENT MC240 (y) and at ADVIA 1800 (x) using 68 serum samples gave following results:

$y = 0.9294x + 0.2903$ g/dl;
 $R = 0.989$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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Calibration Rule	
Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	56
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-200 II GEN

Parameters	
Test Name	ALB
Test No	4
Full Name	ALBUMIN
Reference No	4
Analy. Type	Endpoint
Pri. Wave.	630 nm
Secon. Wave.	
Trend	Ascending
Reac. Time	0 19
Incuba. Time	
Unit	g/dl
Precision	0.01

R1	250
R2	
Sample Volume	2.5
R1 Blank	
Mixed Reag. Blank	
Concentration	1.14 6.50
Linearity Limit	
Substrate Limit	
Factor	
<input type="checkbox"/> Prozone check	
q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule	
Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	56
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters	
Test	ALB
No	4
Full Name	ALBUMIN
Standard No	4
Reac. Type	Endpoint
Pri. Wave.	630 nm
Sec. Wave.	
Direction	Increase
Reac. Time	0 19
Incuba. Time	
Unit	g/dl
Precision	0.01

R1	300
R2	
Sample Volume	3
R1 Blank	
Mixed Rtg. Blank	
Concentration	1.35 6.5
Linearity Limit	
Substrate Limit	
Factor	
<input type="checkbox"/> Prozone check	
q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule	
Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	56
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• BS-120

Parameters	
Test	ALB
No	4
Full Name	Albumin
Standard No	4
Reac. Type	Endpoint
Pri. Wave.	630 nm
Sec. Wave.	
Direction	Increase
Reac. Time	0 17
Incuba. Time	16
Unit	g/dl
Precision	0.01

R1	300
R2	
Sample Volume	3
R1 Blank	
Mixed Rtg. Blank	
Concentration	1.15 6.5
Linearity Limit	
Substrate Limit	
Factor	
<input type="checkbox"/> Prozone check	
q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule	
Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	0
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	ALB	No.	004	Sample Type	SERUM							
Chemistry	ALBUMIN			Print name	ALBUMIN							
Reaction Type	Endpoint			Reaction Direction	positive							
Pri Wave	630 nm			Sec Wave								
Unit	g/dl			Decimal	0.01							
Blank Time	-3	-1		Incubation Time	0							
				Reaction Time	32 34							
Standard	2	μL	Aspirated		μL	Diluent		μL	Reagent Vol		μL	
Decreased	2	μL	20	μL	180	μL	R1	200	μL	R2		μL
Increased		μL		μL								μL
			Sample Blank	V	Auto Rerun							
Linearity range (Standard)	0.26	6.60		Linearity Limit								
Linearity Range (Decreased)				Substrate Depletion								
Linearity Range (Increased)				Mixed Blank Abs	-40000	40000						
R1 Blank Abs	-40000	40000		On-board Stability	84	Day(s)						
Blank Response	-40000	40000		Reagent Alarm Limit								
Twin Chemistry				<input type="checkbox"/> Enzyme Linear Extension								
<input type="checkbox"/> Prozone Check												
Q1		Q2		V1		Q3		Q4		V2		
Q5		Q6		V3		PC1		PC2				
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		<input type="checkbox"/> Calibrator Pretreatment								
		Pretreat Sample Vol		μL		Pretreat Sample Vol		μL				
CALIBRATION SETTINGS						AUTO CALIBRATION						
Math model	Multi-point linear											
Factor		Replicates	2									
ACCEPTANCE LIMITS						Cal Time						
Cal Time	2016	Hour										
Slope Diff				SD								
Sensitivity				Repeatability	40000							
Deter Coeff												

ACCENT-200 ALBUMIN

• ACCENT MC240

Chem	ALB	No.	004	Sample Type	SERUM
Chemistry	ALBUMIN	Print name	ALB	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave		Decimal	0.01
Pri Wave	605 nm	Incubation Time	0	Reaction Time	34 36
Unit	g/dl	Blank Time	-3 -1	Reagent Vol	
Standard	2 μL	Aspirated		R1	180 μL
Decreased	2 μL	Diluent		R2	
Increased		Sample Blank	V	Auto Rerun	

Linearity range (Standard)	0.70	6.10	Linearity Limit			
Linearity Range (Decreased)			Substrate Depletion			
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000		
R1 Blank Abs	-35000	35000	On-board Stability	84 Day(s)		
Blank Response	-35000	35000	Reagent Alarm Limit			
Twin Chemistry		Prozone Check	Enzyme Linear Extension			
Q1	Q2	V1	Q3	Q4	V2	
Q5	Q6	V3	PC1	PC2		
Sample Pretreatment	Control Pretreatment	Calibrator Pretreatment	Pretreat Sample Vol	μL	Pretreat Sample Vol	μL

CALIBRATION SETTINGS	AUTO CALIBRATION			
Math model	Multi-point linear	Bottle Changed		
Factor	Replicates	2	Lot Changed	
			Cal Time	

ACCEPTANCE LIMITS			
Cal Time	2016 Hour		
Slope Diff		SD	
Sensitivity		Repeatability	35000
Deter Coeff			

• ACCENT M320

Chem	ALB	No.	004	Sample Type	SERUM
Chemistry	ALBUMIN	Print name	ALBUMIN	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave		Decimal	0.01
Pri Wave	605 nm	Incubation Time	0	Reaction Time	28 30
Unit	g/dl	Blank Time	-3 -1	Reagent Vol	
Standard	2 μL	Aspirated		R1	180 μL
Decreased	2 μL	Diluent		R2	
Increased		Sample Blank	V	Auto Rerun	

Linearity range (Standard)	1.04	5.80	Linearity Limit			
Linearity Range (Decreased)			Substrate Depletion			
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000		
R1 Blank Abs	-35000	35000	On-board Stability	84 Day(s)		
Blank Response	-35000	35000	Reagent Alarm Limit			
Twin Chemistry		Prozone Check	Enzyme Linear Extension			
Q1	Q2	V1	Q3	Q4	V2	
Q5	Q6	V3	PC1	PC2		
Sample Pretreatment	Control Pretreatment	Calibrator Pretreatment	Pretreat Sample Vol	μL	Pretreat Sample Vol	μL

CALIBRATION SETTINGS	AUTO CALIBRATION			
Math model	Multi-point linear	Bottle Changed		
Factor	Replicates	2	Lot Changed	
			Cal Time	

ACCEPTANCE LIMITS			
Cal Time	2016 Hour		
Slope Diff		SD	
Sensitivity		Repeatability	35000
Deter Coeff			

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