ACCENT-200 LDH

Nr kat. 7-239 (PL)

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Cat. No **7-239** (EN)

INTENDED USE

Diagnostic kit for determination of lactate dehydrogenase activity intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, ACCENT 400, ACCENT Neo200 and BS-120.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Lactate dehydrogenase (LDH, LD) is intracellular enzyme occurred in all tissues. LDH catalyzes the reversible conversion of lactate to pyruvate using NAD⁺ as a cofactor. LD is a tetramer containing two possible forms of subunits: H and M. The result is five isoenzymes termed LD-1 (H₄) through LD-5 (M₄). The isoenzymes are present in different proportion in each tissue and have different electrophoresis mobility, what is very useful for diagnostic.

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The rate of absorbance changing at $\lambda=340$ nm is directly proportional to lactate dehydrogenase activity.

REAGENTS

Package

1-Reagent 1 x 30 ml
2-Reagent 1 x 8 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: 6 weeks (ACCENT-200) or 12 weeks (ACCENT MC240).

Concentrations in the test

phosphate buffer (pH 7.5) 50 mmol/l
pyruvate 0.6 mmol/l
NADH 0.25 mmol/l
preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum, heparinized plasma free from hemolysis. Do not use hemolyzed blood or serum because erythrocytes contain 150 times more LDH activity than serum.

As an anticoagulant for plasma preparation use heparin lithium or ammonium salt.

LDH activity is unstable and is rapidly lost during storage. Specimens can be stored up to 4 hours at 15-25°C or 1-2 days at 2-8°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use. Deionised water is recommended as a reagent blank.

Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120, there is a probability of **cross-contamination** affecting the tests results: ALAT - LDH. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

REFERENCE VALUES ⁴

serum / plasma	37°C	37°C
adults	225 – 450 U/l	3.75 – 7.50 μ kat/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples. For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320 BS-120, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 6 weeks (ACCENT-200) or every 12 weeks (ACCENT MC240), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers ACCENT-200 and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 25.5 U/l (0.43 μ kat/l) - ACCENT-200
- LOQ:** 21 U/l (0.35 μ kat/l) - ACCENT MC240

Linearity

up to 2020 U/l (33.7 μ kat/l) - ACCENT-200
up to 1670 U/l (27.83 μ kat/l) - ACCENT MC240

For higher activity dilute sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

Specificity / Interferences

Haemoglobin up to 5 g/dl, bilirubin up to 20 mg/dl, ascorbate up to 62 mg/l and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [U/l]	SD [U/l]	CV [%]
ACCENT-200 n=20	level 1	257.25	3.52	1.37
	level 2	762.55	3.97	0.52
ACCENT MC240 n=20	level 1	329.5	3.33	1.0
	level 2	774.3	4.05	0.5
Reproducibility (day to day)		Mean [U/l]	SD [U/l]	CV [%]
ACCENT-200 n=80	level 1	434.14	10.88	2.51
	level 2	719.59	12.44	1.73
ACCENT MC240 n=80	level 1	338.4	4.45	1.3
	level 2	805.5	7.55	0.9

Method comparison

A comparison between LDH values determined at ACCENT-200 (y) and at OLYMPUS AU 400 (x) using 30 samples gave following results:

$$y = 1.1108x - 62.097 \text{ U/l};$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

A comparison between LDH values determined at ACCENT MC240 (y) and at BS-800 (x) using 60 serum samples gave following results:

$$y = 1.007x + 1.5792 \text{ U/l};$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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ACCENT-200 LDH

ACCENT MC240

Chem <input type="text" value="LDH"/>	No. <input type="text" value="024"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="LDH"/>	Print name <input type="text" value="LDH"/>	
Reaction Type <input type="text" value="Kinetic"/>	Reaction Direction <input type="text" value="negative"/>	
Pri Wave <input type="text" value="340 nm"/>	Sec Wave <input type="text" value="450 nm"/>	
Unit <input type="text" value="U/l"/>	Decimal <input type="text" value="0.1"/>	
Blank Time <input type="text"/>	Incubation Time <input type="text" value="21"/>	Reaction Time <input type="text" value="7"/> <input type="text" value="20"/>
Standard <input type="text" value="2.4"/> <input type="text" value="2.4"/> <input type="text" value="2.4"/> <input type="text" value="2.4"/> <input type="text" value="2.4"/>	Aspirated <input type="text"/>	Diluent <input type="text"/>
Decreased <input type="text"/>	Reagent Vol <input type="text" value="160"/> <input type="text" value="160"/> <input type="text" value="160"/> <input type="text" value="160"/> <input type="text" value="160"/>	R1 <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/>
Increased <input type="text"/>	R2 <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/>	
<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun	

Linearity range (Standard) <input type="text" value="21"/> <input type="text" value="1670"/>	Linearity Limit <input type="text" value="0.2"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	
Q1 <input type="text"/> Q2 <input type="text"/> V1 <input type="text"/> Q3 <input type="text"/> Q4 <input type="text"/> V2 <input type="text"/>	
Q5 <input type="text"/> O6 <input type="text"/> V3 <input type="text"/> PC1 <input type="text"/> PC2 <input type="text"/>	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment
<input type="checkbox"/> Pretreat Sample Vol <input type="text"/> μ L	<input type="checkbox"/> Calibrator Pretreatment
<input type="checkbox"/> Pretreat Sample Vol <input type="text"/> μ L	

CALIBRATION SETTINGS	AUTO CALIBRATION
Math model <input type="text" value="Multi-point linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="checkbox"/> Cal Time

ACCEPTANCE LIMITS
Cal Time <input type="text" value="2016"/> Hour
Slope Diff <input type="text"/> SD <input type="text"/>
Sensitivity <input type="text"/> Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text"/>

ACCENT M320

Chem <input type="text" value="LDH"/>	No. <input type="text" value="024"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="LDH"/>	Print name <input type="text" value="LDH"/>	
Reaction Type <input type="text" value="Kinetic"/>	Reaction Direction <input type="text" value="negative"/>	
Pri Wave <input type="text" value="340 nm"/>	Sec Wave <input type="text" value="450 nm"/>	
Unit <input type="text" value="U/l"/>	Decimal <input type="text" value="0.1"/>	
Blank Time <input type="text"/>	Incubation Time <input type="text" value="15"/>	Reaction Time <input type="text" value="10"/> <input type="text" value="30"/>
Standard <input type="text" value="2.4"/> <input type="text" value="2.4"/> <input type="text" value="2.4"/> <input type="text" value="2.4"/> <input type="text" value="2.4"/>	Aspirated <input type="text"/>	Diluent <input type="text"/>
Decreased <input type="text"/>	Reagent Vol <input type="text" value="160"/> <input type="text" value="160"/> <input type="text" value="160"/> <input type="text" value="160"/> <input type="text" value="160"/>	R1 <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/>
Increased <input type="text"/>	R2 <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/>	
<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun	

Linearity range (Standard) <input type="text" value="21"/> <input type="text" value="1560"/>	Linearity Limit <input type="text" value="0.2"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	
Q1 <input type="text"/> Q2 <input type="text"/> V1 <input type="text"/> Q3 <input type="text"/> Q4 <input type="text"/> V2 <input type="text"/>	
Q5 <input type="text"/> O6 <input type="text"/> V3 <input type="text"/> PC1 <input type="text"/> PC2 <input type="text"/>	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment
<input type="checkbox"/> Pretreat Sample Vol <input type="text"/> μ L	<input type="checkbox"/> Calibrator Pretreatment
<input type="checkbox"/> Pretreat Sample Vol <input type="text"/> μ L	

CALIBRATION SETTINGS	AUTO CALIBRATION
Math model <input type="text" value="Multi-point linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="checkbox"/> Cal Time

ACCEPTANCE LIMITS
Cal Time <input type="text" value="2016"/> Hour
Slope Diff <input type="text"/> SD <input type="text"/>
Sensitivity <input type="text"/> Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text"/>

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