ACCENT-200 ASO

Nr kat. **7-240** (PL)

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Cat. No **7-240** (EN)

INTENDED USE

Diagnostic kit for determination of anti-streptolysin O levels intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320 ACCENT 400 and ACCENT Neo200.
 The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Most people infected with hemolytic streptococcus produce anti-streptolysin O (ASO), antibodies against streptolysin O (SLO), an exotoxin of *Streptococcus* spp. Measuring the level of ASO is effective for diagnosing, judging the progress of medical treatment, and assessing recovery from diseases caused by hemolytic streptococcus such as rheumatic fever, acute glomerulonephritis, scarlatina and tonsillitis.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between ASO in a sample and SLO which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of ASO in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

1-REAGENT 1 x 25 ml
 2-REAGENT 1 x 17.5 ml

The reagents, stored at 2-10°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 7 weeks.

Concentrations in the test

suspension of latex particles sensitized with SLO (pH 8.2) 0.17 w/v%
 glycine buffer solution (pH 8.3)
 preservative

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Reagents with different lot numbers should not be interchanged or mixed.
- Do not interchange the caps of reagent bottles.
- Reagents should be mixed before use by gentle inverting the bottle several times
- Immediately after use, recap the bottles and store at 2-10°C.
- EUH210 Safety data sheet available on request.

SPECIMEN

Serum or plasma collected on heparine (sodium or lithium salt), EDTA (sodium or potassium salt) or citric acid.
 After the blood is completely clotted, the sample should be centrifuged and separated from the cells and fibrinogen. If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Avoid repeated freezing and thawing. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.
 0.9% NaCl is recommended as a reagent blank

REFERENCE VALUES ³

serum, plasma	< 160 IU/ml
---------------	-------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL I (Cat. No 4-288) with each batch of samples.

For the calibration of automatic analysers: ACCENT-200 and ACCENT-200 II GEN the CORMAY ASO CALIBRATOR kit (Cat. No 4-278) is recommended.

For the calibration of automatic analysers: ACCENT-220S, ACCENT S120, ACCENT MC240 and ACCENT M320 the CORMAY ASO CALIBRATOR kit (Cat. No 4-278) is recommended. 0.9% NaCl should be used as a calibrator 0.
 The calibration curve should be prepared every week, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers ACCENT-200 and/or ACCENT-220S and ACCENT MC240. Results may vary if a different instrument is used.

Sensitivity

20.4 IU/ml – ACCENT-220S
 40 IU/ml – ACCENT MC240

Linearity

up to 560 IU/ml – ACCENT-220S
 up to 845 IU/ml – ACCENT MC240

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 0.5 g/dl, bilirubin up to 20 mg/dl and triglycerides up to 500 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [IU/ml]	SD [IU/ml]	CV [%]
ACCENT-200 n=20	level 1	28.39	0.92	3.22
	level 2	320.15	2.94	0.92
ACCENT MC240 n=20	level 1	142.82	6.92	4.85
	level 2	306.37	7.07	2.31
Reproducibility (day to day)		Mean [IU/ml]	SD [IU/ml]	CV [%]
ACCENT-200 n=56	level 1	35.84	1.63	4.54
	level 2	129.76	4.11	3.17
ACCENT MC240 n=80	level 1	139.4	5.16	3.7
	level 2	317.9	7.46	2.3

Method comparison

A comparison between ASO values determined at **ACCENT-220S** (y) and at **ADVIA 1650** (x) using 56 samples gave following results:

$y = 0.8814x + 4.0817$ IU/l;
 $R = 0.966$ (R – correlation coefficient)

A comparison between ASO values determined at **ACCENT MC240** (y) and at **BECKMAN COULTER AU 680** (x) using 50 serum samples gave following results:

$y = 0.9917x + 2.35$ IU/l;
 $R = 0.996$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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ACCENT-200 ASO

• ACCENT MC240

Chem	ASO	No.	030	Sample Type	SERUM
Chemistry	ASO	Print name	ASO	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave		Decimal	0.1
Pri Wave	700nm	Incubation Time	21	Reaction Time	28 30
Unit	IU/ml	Blank Time	-3 -1	Standard	Sample Vol 3.5 μL Aspirated μL Diluent μL
		Decreased	3.5 μL 20 μL 180 μL	Reagent Vol	R1 160 μL R2 110 μL
		Increased	μL μL μL	Sample Blank	<input type="checkbox"/> V Auto Rerun
Linearity range (Standard)	40 845	Linearity Limit		Linearity Range (Decreased)	
Linearity Range (Increased)		Substrate Depletion		Linearity Range (Increased)	
R1 Blank Abs	-35000 35000	Mixed Blank Abs	-35000 35000	On-board Stability	Day(s)
Blank Response	-35000 35000	Reagent Alarm Limit		Enzyme Linear Extension	<input type="checkbox"/>
Twin Chemistry		Prozone Check	<input type="checkbox"/>	Q1	Q2 V1 Q3 Q4 V2
Q5	Q6 V3 PC1 PC2	Sample Pretreatment		Control Pretreatment	
		Calibrator Pretreatment		Pretreat Sample Vol	μL
		Pretreat Sample Vol	μL		
CALIBRATION SETTINGS		AUTO CALIBRATION			
Math model	Two-point Linear	Bottle Changed			
Factor	Replicates 2	Lot Changed			
		Cal Time			
ACCEPTANCE LIMITS					
Cal Time	Hour				
Slope Diff		SD			
Sensitivity	Repeatability 35000				
Deter Coeff					

• ACCENT M320

Chem	ASO	No.	030	Sample Type	SERUM
Chemistry	ASO	Print name	ASO	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave		Decimal	0.01
Pri Wave	700nm	Incubation Time	15	Reaction Time	28 30
Unit	IU/ml	Blank Time	-3 -1	Standard	Sample Vol 3.5 μL Aspirated μL Diluent μL
		Decreased	3.5 μL 20 μL 180 μL	Reagent Vol	R1 160 μL R2 110 μL
		Increased	μL μL μL	Sample Blank	<input type="checkbox"/> V Auto Rerun
Linearity range (Standard)	31 800	Linearity Limit		Linearity Range (Decreased)	
Linearity Range (Increased)		Substrate Depletion		Linearity Range (Increased)	
R1 Blank Abs	-35000 35000	Mixed Blank Abs	-35000 35000	On-board Stability	Day(s)
Blank Response	-35000 35000	Reagent Alarm Limit		Enzyme Linear Extension	<input type="checkbox"/>
Twin Chemistry		Prozone Check	<input type="checkbox"/>	Q1	Q2 V1 Q3 Q4 V2
Q5	Q6 V3 PC1 PC2	Sample Pretreatment		Control Pretreatment	
		Calibrator Pretreatment		Pretreat Sample Vol	μL
		Pretreat Sample Vol	μL		
CALIBRATION SETTINGS		AUTO CALIBRATION			
Math model	Two-point Linear	Bottle Changed			
Factor	Replicates 2	Lot Changed			
		Cal Time			
ACCEPTANCE LIMITS					
Cal Time	Hour				
Slope Diff		SD			
Sensitivity	Repeatability 35000				
Deter Coeff					

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