ACCENT-200 GLUCOSE HEX

Nr kat. **7-252** (PL)

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LITERATURA

- Barham P., Trinder P.: Analyst 97, 142-145 (1972).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 868-869, 2006.
- McPherson Richard A., Pincus Matthew R.: Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd ed. Philadelphia, PA: WB Saunders, 2011.
- Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine, Volume 48, Issue 2, 2010, pp. 209-212.
- Zalecenia kliniczne dotyczące postępowania u chorych na cukrzycę 2019, Diabetologia Praktyczna, tom 5, nr 1, 2019.
- Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical chemistry 48.3 (2002): 436-472.
- Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clin Chem 2004;50:1704-1706.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
- Young DS., Effects of drugs on clinical laboratory tests, 5. Vol. 2. Washington DC, USA: AACC Press (2000).

Data wydania: 07.2023.

ACCENT-200 GLUCOSE HEX

Cat. No 7-252 (EN)

INTENDED USE

Diagnostic kit for determination of glucose concentration used in automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, ACCENT 400 and ACCENT Neo200.



The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Glucose is a simple six-carbon sugar. Oxidative metabolism of glucose provides the energy for most cellular processes. Glucose level in the blood is tightly controlled by several hormones. Elevated glucose level is the classic sign of diabetes mellitus. Glucose level abnormalities (hyper- or hypoglycemia) might be caused also by pancreas tumors and diseases of liver, thyroid gland or adrenal glands.

METHOD PRINCIPLE

Enzymatic method with hexokinase and glucose-6-phosphate dehydrogenase (G6P-DH).

The rate of NADH formation is directly proportional to the glucose concentration in the sample.

REAGENTS

Package
1-REAGENT 4 x 28.5 ml
2-REAGENT 4 x 6.5 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 12 weeks (ACCENT-200, ACCENT MC240).

Concentrations in the test

1-REAGENT
PIPES buffer (pH 7.5) 80 mmol/l
Mg²⁺ 10 mmol/l
ATP 4 mmol/l
NAD 3 mmol/l
2-REAGENT
hexokinase ≥ 4500 U/l
glucose-6-phosphate dehydrogenase (G6P-DH) ≥ 14000 U/l
preservative

WARNINGS AND NOTES

- Do not use after expiry date.
- Do not freeze reagents.
- Do not interchange caps.
- Protect from direct sunlight and avoid contamination!
- Reagents should be mixed before use by gentle inverting the bottles several times.
- EUH210 Safety data sheet available on request.

SPECIMEN

EDTA or heparinized plasma / serum, free from hemolysis, cerebrospinal fluid, urine.

Plasma / Serum. Serum and plasma specimens should be separated from cells within 30 minutes after collection. Plasma specimen which is not assayed immediately after collection should be kept in tubes containing sodium fluoride or sodium iodoacetate. These compounds adding prevent glycolysis and stabilize glucose level.

Serum and plasma can be stored up to 2 days at 4°C.³

Plasma is the specimen recommended for the glucose determination in the blood.⁵

Cerebrospinal fluid. Glucose concentration in cerebrospinal fluid should be measured directly after specimen collection. Cerebrospinal fluid must be analysed simultaneously with a blood sample.

After centrifuge CSF sample can be stored up to 24 hours at 4°C.⁴

Urine. Collect 24-hour sample in dark bottle and keep on ice. Preserve sample by adding 5 ml of glacial acetic acid to the container before starting the collection. The final pH of the sample should be between 4 and 5. Centrifuge samples with visible turbidity or precipitates before analysis.

Urine can be stored up to 24 hour at 4°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT and 2-REAGENT are ready to use. Deionised water is recommended as a reagent blank.

Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120, there is a probability of **cross-contamination** affecting the tests results: GLUCOSE HEX - CK-MB, GLUCOSE HEX - MG, ASAT - GLUCOSE HEX. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

REFERENCE VALUES

	mg/dl	mmol/l
plasma, serum ^{5,6,7}	70 - 99	3.9 - 5.5
urine (24h) ⁸	1 - 15	0.1 - 0.8
cerebrospinal fluid ⁸	40 - 70	2.2 - 3.9

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the following controls for each batch of samples: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended.

For the calibration of automatic analysers: ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120 the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 11 weeks (ACCENT-200) or every 12 weeks (ACCENT MC240), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers: ACCENT-200 and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

Sensitivity

8.7 mg/dl (0.48 mmol/l) - ACCENT-200

LOQ (Limit of Quantitation)

2.6 mg/dl (0.143 mmol/l) - ACCENT MC240

Linearity

up to 1000 mg/dl (55 mmol/l) - ACCENT-200

up to 1050 mg/dl (57.75 mmol/l) - ACCENT MC240

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 1.25 g/dl, bilirubin up to 40 mg/dl, ascorbate up to 62 mg/L and triglycerides up to 1000 mg/dl do not interfere with the test. Some medicines can interfere.⁹

Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=10	level 1	88.16	0.92	1.04
	level 2	293.31	3.02	1.03
ACCENT MC240 n=20	level 1	82.26	1.04	1.27
	level 2	283.02	3.45	1.22
Repeatability (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=20	level 1	87.62	1.56	1.78
	level 2	283.02	6.70	2.37
ACCENT MC240 n=80	level 1	83.10	1.15	1.40
	level 2	286.30	4.93	1.70

Method comparison

A comparison between glucose values determined at **ACCENT-200** (y) and at **COBAS INTEGRA 400** (x) using 39 samples of serum gave following results:

$$y = 0.99x + 1.2288 \text{ mg/dl;} \\ R = 0.999 \quad (R - \text{correlation coefficient})$$

A comparison between glucose values determined at **ACCENT MC240** (y) and at **BS-800** (x) using 60 samples of serum gave following results:

$$y = 0.9773x + 2.8799 \text{ mg/dl;} \\ R = 0.998 \quad (R - \text{correlation coefficient})$$

A comparison between glucose values determined at **ACCENT MC240** (y) and at **ADVIA 1800** (x) using 57 plasma samples gave following results:

$$y = 1.0324x - 1.7485 \text{ mg/dl;} \\ R = 0.998 \quad (R - \text{correlation coefficient})$$

A comparison between glucose values determined at **ACCENT MC240** (y) and at **BS-800** (x) using 57 urine samples gave following results:

$$y = 1.0349x - 3.1983 \text{ mg/dl;} \\ R = 1.000 \quad (R - \text{correlation coefficient})$$

A comparison between glucose values determined at **ACCENT MC240** (y) and at **BS-800** (x) using 64 cerebrospinal fluid samples gave following results:

$$y = 1.0058x + 0.7701 \text{ mg/dl;} \\ R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Barham P., Trinder P.: Analyst 97, 142-145 (1972).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 868-869, 2006.
- McPherson Richard A., Pincus Matthew R.: Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd ed. Philadelphia, PA: WB Saunders, 2011.
- Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine, Volume 48, Issue 2, 2010, pp. 209-212.
- Zalecenia kliniczne dotyczące postępowania u chorych na cukrzycę 2019, Diabetologia Praktyczna, tom 5, nr 1, 2019.
- Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical chemistry 48.3 (2002): 436-472.
- Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clin Chem 2004;50:1704-1706.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
- Young DS., Effects of drugs on clinical laboratory tests, 5. Vol. 2. Washington DC, USA: AACC Press (2000).

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- Barham P., Trinder P.: Analyst 97, 142-145 (1972).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 868-869, 2006.
- McPherson Richard A., Pincus Matthew R.: Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd ed. Philadelphia, PA: WB Saunders, 2011.
- Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine, Volume 48, Issue 2, 2010, pp. 209-212.
- Zalecenia kliniczne dotyczące postępowania u chorych na cukrzycę 2019, Diabetologia Praktyczna, tom 5, nr 1, 2019
- Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical chemistry 48.3 (2002): 436-472.
- Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clin Chem 2004;50:1704-1706.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
- Young DS., Effects of drugs on clinical laboratory tests, 5. Vol. 2. Washington DC, USA: AACC Press (2000).

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Calibration Rule

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	77
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-200 II GEN

Parameters			
Test Name	GLUCOSE HEX	R1	200
Test No	76	R2	40
Full Name	Glucose Hexokinase	Sample Volume	3.5
Reference No	76	R1 Blank	
Analy. Type	Endpoint	Mixed Reag. Blank	
Pri. Wave.	340 nm	Concentration	9.1 1000
Secon. Wave.		Linearity Limit	
Trend	Ascending	Substrate Limit	
Reac. Time	-1 15	Factor	
Incuba. Time	10	<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="checkbox"/>	Abs <input type="checkbox"/>

Calibration Rule

Rule	Two-point Linear
Sensitivity	0
Replicates	3
Interval (day)	77
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters			
Test	GLUC HEX	R1	200
No	76	R2	40
Full Name	Glucose Hexokinase	Sample Volume	3
Standard No	76	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	340 nm	Linearity Range	4.3 1000
Sec. Wave.		Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	-1 17	Factor	
Incuba. Time	11	<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="checkbox"/>	Abs <input type="checkbox"/>

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	77
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• BS-120

Parameters			
Test	GLUCOSE HEX	R1	180
No	1	R2	36
Full Name	Glucose Hexokinase	Sample Volume	3.5
Standard No	1	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	340 nm	Linearity Range	4.5 990
Sec. Wave.		Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	-1 15	Factor	
Incuba. Time	16	<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="checkbox"/>	Abs <input type="checkbox"/>

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	77
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	GLUC HEX	No.	076	Sample Type	SERUM/ Urine/CSF
Chemistry	GLUCOSE HEXOKINASE			Print name	GLUC HEX
Reaction Type	Endpoint			Reaction Direction	positive
Pri Wave	340 nm			Sec Wave	
Unit	mg/dl			Decimal	0.1
Blank Time	-3	-1		Incubation Time	11
Standard	2.3	Aspirated		Reaction Time	16 18
Decreased	2.3	20	180	Reagent Vol	
Increased				R1	180
				R2	36
		Sample Blank	<input checked="" type="checkbox"/>	Auto Rerun	
Linearity range (Standard)	10	1020		Linearity Limit	
Linearity Range (Decreased)				Substrate Depletion	
Linearity Range (Increased)				Mixed Blank Abs	-40000 40000
R1 Blank Abs	-40000	40000		On-board Stability	
Blank Response	-40000	40000		Reagent Alarm Limit	
Twin Chemistry				<input type="checkbox"/> Enzyme Linear Extension	
<input type="checkbox"/> Prozone Check					
Q1		Q2		V1	
Q5		Q6		V3	
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		PC1	
		<input type="checkbox"/> Calibrator Pretreatment		PC2	
		Pretreat Sample Vol		uL	
		Pretreat Sample Vol		uL	
CALIBRATION SETTINGS		AUTO CALIBRATION			
Math model	Multi-point Linear	<input type="checkbox"/> Bottle Changed			
Factor		<input type="checkbox"/> Lot Changed			
		<input type="checkbox"/> Cal Time			
ACCEPTANCE LIMITS					
Cal Time		Hour			
Slope Diff		SD			
Sensitivity		Repeatability	40000		
Deter Coeff					