ACCENT-200 LACTATE

Nr kat. 7-266

(PL)

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LITERATURA

- Tietz Textbook of Clinical Chemistry (Edited by Burtis CA and Ashwood ER Eds): Third Edition WB Saunders Company 787-8, (1999).
- Alan H. B. Wu, Tietz Clinical Guide to Laboratory Tests, W.B. Saunders Company, 4th edition, 650-652, (2006).

Data wydania: 05. 2022.

ACCENT-200 LACTATE

Cat. No 7-266

(EN)

INTENDED USE

Diagnostic kit for determination of lactate concentration intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240 ACCENT M320, BS-120, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Lactate is produced in Cori cycle, by anaerobic conversion of glucose, mainly in skeletal muscle. Its determination, frequently done together with pyruvate, is useful in discovering lactic acidosis due to i.a. reduced tissue oxygenation, enzymatic deficiencies, diabetes mellitus, liver and kidneys diseases.

METHOD PRINCIPLE

Lactate is oxidized by lactate oxidase to pyruvate and hydrogen peroxide, which, in presence of peroxidase (POD), reacts with 4-aminoantipyrine and phenol forming a compound, which colour intensity is proportional to the concentration of lactate in the examined sample.

REAGENTS

Package

1-Reagent 2 x 18.5 ml

Unopened reagent is stable up to the kit expiry date printed on the package when stored at 2-8°C. The reagents stored on board of the analyser at 2-10°C are stable for 11 weeks.

Concentrations in the test

Tris buffer (pH 7.5) ≥ 50 mmol/l
lactate oxidase ≥ 0.2 kU/l
peroxidase ≥ 2 kU/l
4-aminoantipyrine ≥ 0.4 mmol/l

Warnings and notes

- Do not freeze the reagent!
- Protect from direct sunlight and avoid contamination!
- Do not use after expiry date.
- The appearance of turbidity or control sera values outside the manufacturer's acceptable range may indicate of the reagents instability.
- Lactate concentration rapidly increases during physical activities. Normal levels are reached again after usually 30 minutes but it may vary according to individuals.
- Draw blood with lowest venous stasis as possible (max. 30 seconds) from fasting and completely resting patient and avoid using a tourniquet.

SPECIMEN

Plasma. free from haemolysis.

Collect samples in tubes containing sodium fluoride and potassium oxalate. Keep samples on ice. Centrifuge within 15 minutes after collection and separate from cells. Analyze promptly. Note whether sample is venous or arterial. It is recommended to follow NCCLS procedures regarding specimen collecting and handling.

Lactate in plasma is stable up to 8 hours at room temperature or up to 14 days at 2-8°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent is ready to use.

Deionised water is recommended as a reagent blank.

Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120, there is a probability of **cross-contamination** affecting the tests results: LACTATE – CRP ULTRA, LACTATE - FERRITIN. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

REFERENCE VALUES ?

plasma (venous)	4.5 – 19.8 mg/dl	0.5 – 2.2 mmol/l
plasma (arterial)	4.5 – 14.4 mg/dl	0.5 – 1.6 mmol/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320 the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

For the calibration of automatic analyser BS-120 the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 11 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers ACCENT-200 and/or ACCENT-200 II GEN and/or ACCENT-220S. Results may vary if a different instrument or a manual procedure is used.

▪ Sensitivity

1.7 mg/dl (0.19 mmol/l) - ACCENT-200
2.7 mg/dl (0.3 mmol/l) – ACCENT MC240

▪ Linearity

up to 90 mg/dl (9.99 mmol/l) - ACCENT-200
up to 116 mg/dl (12.88 mmol/l) – ACCENT MC240

For higher concentration dilute the sample 1:1 with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

▪ Specificity / Interferences

a) In plasma samples containing approximately 12 mg lactate/dl, there is no interference up to: 0.23 g/dl haemoglobin, 8 mg/dl bilirubin, 337 mg/dl triglycerides, 15.5 mg/l ascorbic acid.

b) In plasma samples containing approximately 40 mg lactate/dl, there is no interference up to: 1.25 g/dl haemoglobin, 10 mg/dl bilirubin, 1000 mg/dl triglycerides, 62 mg/l ascorbic acid.

▪ Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=10	level 1	9.2	0.13	1.49
	level 2	41.59	0.29	0.69
ACCENT MC240 n=20	level 1	8.59	0.07	0.83
	level 2	41.81	0.46	1.11
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=10	level 1	42.84	1.01	2.35
	level 2	9.12	0.11	1.20
ACCENT MC240 n=80	level 1	40.9	1.35	3.3
	level 2	9.2	0.31	3.4

▪ Method comparison

A comparison between CORMAY reagent (y) and commercially available assay (x) using 40 samples gave following results:

$$y = 1.0219x + 0.2743 \text{ mg/dl;}$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

A comparison between lactate values determined at ACCENT MC240 (y) and at BS-400 (x) using 60 samples gave following results:

$$y = 0.9901x + 1.5281 \text{ mg/dl;}$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Tietz Textbook of Clinical Chemistry (Edited by Burtis CA and Ashwood ER Eds): Third Edition WB Saunders Company 787-8, (1999).
- Alan H. B. Wu, Tietz Clinical Guide to Laboratory Tests, W.B. Saunders Company, 4th edition, 650-652, (2006).

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