## ACCENT-200 TG mono

Nr kat. 7-273

(PL)

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lipoprotein lipase (LPL) ~2000 U/l  
 detergents, preservatives

### INTENDED USE

Diagnostic kit for determination of triglycerides concentration intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, ACCENT 400 and ACCENT Neo200.



The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Triglycerides are built of glycerol molecule esterified with three fatty acids molecules. Triglycerides are delivered with food or are synthesized endogenously in liver. Triglycerides stored in adipose tissue constitute a reserve of energy. Elevated triglycerides serum level is a risk factor of atherosclerosis. Triglycerides measurement is useful for hyperlipidemia diagnosis and treatment or for estimation of atherosclerosis progression.

### METHOD PRINCIPLE

Colorimetric, enzymatic method with glycerophosphate oxidase.

The colour intensity is proportional to the triglycerides concentration.

### REAGENTS

**Package**  
 1-Reagent 6 x 31 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 7 weeks (ACCENT-200) or 12 weeks (ACCENT M320).

### Concentrations in the test

buffer TRIS (PH 8.0)	200 mmol/l
4-aminoantipyrine (4-AA)	< 0.4 mmol/l
ATP	< 1.5 mmol/l
Mg <sup>2+</sup>	< 1.6 mmol/l
4-chlorophenol	< 2.5 mmol/l
chlorophenicol	1.6 mmol/l
potassium hexacyanoferrate (II)	< 1 mmol/l
FAD-2Na	< 1 mmol/l
glycerol kinase (GK)	~2500 U/l
glycerol phosphate oxidase (GPO)	~2500 U/l
peroxidase (POD)	~1900 U/l

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum, EDTA or heparinized plasma (recommended: heparine lithium, sodium or ammonium salt) free from hemolysis. Blood should be collected only if the patient has been fasting for minimum of 12 hours. Before blood collection patient should stay in rest position for about 30 minutes. Venous blood is recommended for triglycerides measurement. Plasma triglycerides values have been reported to be 2% to 4% lower than serum triglycerides values. Serum and plasma can be stored up to 3 days at 2-8°C or up to 3 months at -20°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use. Deionised water is recommended as a reagent blank. For analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 it is recommended to determine the reagent blank during each calibration. Deionized water should be used as reagent blank. When performing calibration, the task type **Calib+Rgt.Blk** should be selected.

### Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 there is a probability of **cross-contamination** affecting the tests results: TG mono - LIPASE, HDL DIRECT II GEN – TG mono, LDL DIRECT II GEN – TG mono, URINE PROTEINS II GEN – TG mono. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_001\_ACCENT-200\_CARRYOVER.

### REFERENCE VALUES

serum, plasma	< 150 mg/dl < 1.7 mmol/l
---------------	-----------------------------

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples. For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, BS-120, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

For the calibration of automatic analysers ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320 the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended.

Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 10 weeks (ACCENT-200) or 12 weeks (ACCENT M320), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers ACCENT-200 and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

#### Sensitivity

1.36 mg/dl (0.02 mmol/l) - ACCENT-200  
 4.00 mg/dl (0.05 mmol/l) - ACCENT MC240

#### Linearity:

up to 1300 mg/dl (14.7 mmol/l) - ACCENT-200  
 up to 1460 mg/dl (16.5 mmol/l) - ACCENT MC240

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Haemoglobin up to 0.31 g/dl, bilirubin up to 8.6 mg/dl and ascorbate up to 31 mg/l do not interfere with the test.

#### Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=20	level 1	82.77	0.47	0.57
	level 2	157.62	1.28	0.81
ACCENT MC240 n= 20	level 1	83.55	0.95	1.14
	level 2	178.94	2.39	1.34
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=20	level 1	84.27	2.10	2.49
	level 2	168.60	3.70	2.20
ACCENT MC240 n= 80	level 1	86.8	1.77	2.0
	level 2	189.3	3.82	2.0

#### Method comparison

A comparison between triglycerides values determined at ACCENT-200 (y) and at COBAS INTEGRA 400 (x) using 92 samples gave following results:  
 $y = 0.9833x + 4.4927$  mg/dl;  
 $R = 0.997$  (R – correlation coefficient)

A comparison between triglycerides values determined at ACCENT MC240 (y) and at ADVIA SIEMENS 1800 (x) using 60 serum samples gave following results:  
 $y = 0.9242x + 4.9586$  mg/dl;  
 $R = 0.999$  (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

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### Calibration Rule

Rule	One-point Linear
Sensitivity	1
Replicates	3
Interval (day)	49
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

### • ACCENT-200 II GEN

<b>Parameters</b>	
Test Name	TG mono
Test No	75
Full Name	Triglycerides mono
Reference No	75
Analy. Type	Endpoint
Pri. Wave.	546 nm
Secon. Wave.	670 nm
Trend	Ascending
Reac. Time	0   30
Incuba. Time	
Unit	mg/dl
Precision	0.1
R1	300
R2	
Sample Volume	3
R1 Blank	
Mixed Reag. Blank Concentration	1.38   1100
Linearity Limit	
Substrate Limit	
Factor	
<input type="checkbox"/> Prozone check	
q1	<input type="checkbox"/>
q2	<input type="checkbox"/>
q3	<input type="checkbox"/>
q4	<input type="checkbox"/>
PC	<input type="checkbox"/>
Abs	<input type="checkbox"/>

### Calibration Rule

Rule	One-point Linear
Sensitivity	1
Replicates	3
Interval (day)	49
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

### • ACCENT-220S

<b>Parameters</b>	
Test	TG mono
No	75
Full Name	Triglycerides mono
Standard No	75
Reac. Type	Endpoint
Pri. Wave.	546 nm
Sec. Wave.	670 nm
Direction	Increase
Reac. Time	0   30
Incuba. Time	
Unit	mg/dl
Precision	Integer
R1	300
R2	
Sample Volume	3
R1 Blank	
Mixed Rtg. Blank Concentration	2.6   1500
Linearity Limit	
Substrate Limit	
Factor	
<input type="checkbox"/> Prozone check	
q1	<input type="checkbox"/>
q2	<input type="checkbox"/>
q3	<input type="checkbox"/>
q4	<input type="checkbox"/>
PC	<input type="checkbox"/>
Abs	<input type="checkbox"/>

### Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	49
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

### • BS-120

<b>Parameters</b>	
Test	TG mono
No	3
Full Name	Triglycerides mono
Standard No	3
Reac. Type	Endpoint
Pri. Wave.	546 nm
Sec. Wave.	670 nm
Direction	Increase
Reac. Time	0   16
Incuba. Time	
Unit	mg/dl
Precision	0.1
R1	300
R2	
Sample Volume	3
R1 Blank	
Mixed Rtg. Blank Concentration	0.7   1400
Linearity Limit	
Substrate Limit	
Factor	
<input type="checkbox"/> Prozone check	
q1	<input type="checkbox"/>
q2	<input type="checkbox"/>
q3	<input type="checkbox"/>
q4	<input type="checkbox"/>
PC	<input type="checkbox"/>
Abs	<input type="checkbox"/>

### Calibration Rule

Rule	One-point Linear
Sensitivity	1
Replicates	3
Interval (day)	0
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

### • ACCENT S120

Chem	TG MONO	No.	075	Sample Type	SERUM
Chemistry	TG MONO	Print name	TG MONO	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave	670nm	Decimal	0.1
Pri Wave	546nm	Incubation Time	0	Reaction Time	28   30
Unit	mg/dL	Blank Time	-3   -1	Standard	Sample Vol
Sample Vol	3 μL	Aspirated	μL	Diluent	μL
Decreased	3 μL	20 μL	180 μL	Reagent Vol	R1 200 μL
Increased	μL	μL	μL	R2	μL
<input type="checkbox"/> Sample Blank		<input type="checkbox"/> V		<input type="checkbox"/> Auto Rerun	
Linearity range (Standard)	4,1	1170	Linearity Limit		
Linearity Range (Decreased)			Substrate Depletion		
Linearity Range (Increased)			Mixed Blank Abs	-40000	40000
R1 Blank Abs	-40000	40000	On-board Stability		
Blank Response	-40000	40000	Reagent Alarm Limit		
Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension		
<input type="checkbox"/> Prozone Check					
Q1		Q2		V1	
Q5		Q6		V3	
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		<input type="checkbox"/> Calibrator Pretreatment	
Pretreat Sample Vol			μL	Pretreat Sample Vol	
μL			μL		
<b>CALIBRATION SETTINGS</b>					
Math model	Multi-point Linear		<input type="checkbox"/> Bottle Changed		
Factor		Replicates	2	<input type="checkbox"/> Lot Changed	
<b>ACCEPTANCE LIMITS</b>					
Cal Time		Hour			
Slope Diff		SD			
Sensitivity		Repeatability	40000		
Deter Coeff					