ACCENT-300 IgG

Nr kat. 7-303

(PL)

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Cat. No **7-303** (EN)

INTENDED USE

Diagnostic kit for determination of IgG concentration, used in automatic analyser ACCENT-300. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Immunoglobulins (Igs) are the instrumental proteins of immunity. Immunity is a property of the lymphoid system which is made of organs (spleen, thymus, bone marrow) and of cells (lymphocytes). Circulating immunoglobulins are secreted in the blood by B lymphocytes and they thereby export far-away the specific biological functions of humoral immunity.

Immunoglobulin G (IgG) - the predominant serum immunoglobulin (75 % of Igs)- is of particular importance in the body's long-term defense against infection; IgG deficiency is associated with recurrent and occasionally severe pyogenic infections. IgG serum levels are increased in response to chronic or recurrent infections or autoimmune diseases.

METHOD PRINCIPLE

The IgG present in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum is proportional to IgG concentration in the sample.

REAGENTS

Package
 1-Reagent 1 x 40 ml
 2-Reagent 1 x 13.5 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: tricine buffer, polymer, inorganic salt, preservative.

2-Reagent: anti-human IgG antiserum, buffer, inorganic salt, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found

to be non-reactive. However this material should be handled as potentially infectious.

- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN ⁴

Serum. Nonhemolyzed, nonlipaemic, fresh serum is recommended.

Serum can be stored up to 3 days at 2-8°C or up to 6 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ⁴

adults	5.65 – 17.65 g/l
children (> 6 years)	6.50 – 16.00 g/l
children (1 month – 6 years)	2.00 – 12.40 g/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers systems the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using an automatic analyser Hitachi. Results may vary if a different instrument is used.

- Measurement range:** 0.125 g/l to 35 g/l.

- Specificity / Interferences**

Hemoglobin up to 0.32 g/dl, bilirubin up to 22 mg/dl, triglycerides up to 940 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	6.83	0.01	1.43
level 2	14.03	0.16	1.12
level 3	19.86	0.36	1.80

Reproducibility (day to day) n = 66	Mean [g/l]	SD [g/l]	CV [%]
level 1	6.57	0.29	4.46
level 2	13.66	0.43	3.15
level 3	18.87	0.56	2.95

- Method comparison**

A comparison between CORMAY reagent (y) and another commercially available assay (x) using minimum 30 samples gave following results:

$$y = 1.02x - 0.00325 \text{ g/l;}$$

$$R = 0.993 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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