## ACCENT-300 UA

Nr kat. 7-308 (PL)

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Cat No **7-308** (EN)

### INTENDED USE

Diagnostic kit for determination of uric acid concentration, used in automatic analyser ACCENT-300.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Uric acid is a product of purine catabolism. It is produced in the liver and excreted in the urine. Both, the amount of uric acid production and the efficiency of renal excretion, affect serum urate level. Elevated serum uric acid level is caused usually by gout, leukemia, diabetes mellitus, hyperfunction of parathyroid and thyroid, renal failure, renal calculus. Urate concentration in serum and in urine depends on glomerular filtration, thus is useful for renal function monitoring.

### METHOD PRINCIPLE

Enzymatic, colorimetric method with uricase and peroxidase.

The colour intensity is proportional to the uric acid concentration.

### REAGENTS

#### Package

1-Reagent 2 x 40 ml  
2-Reagent 1 x 20.5 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 11 weeks.

#### Concentrations in the test

buffer PIPES (pH 7.0) 100 mmol/l  
4-aminoantipyrine 0.78 mmol/l  
ADPS 0.67 mmol/l  
ferricyanide potassium 3.8 µmol/l  
peroxidase (POD) > 38.34 µkat/l  
uricase > 1.65 µkat/l  
preservative

#### Warnings and notes

- Protect from direct sunlight and contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

#### Warning

H315 Causes skin irritation.  
H319 Causes serious eye irritation.  
P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 IF ON SKIN: Wash with plenty of soap and water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### SPECIMEN

Serum, heparinized plasma free from hemolysis. Do not use EDTA and fluoride as anticoagulants. Specimen can be stored 3-5 days at 2-8°C or 6 months at -20°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use. Deionized water is recommended as a reagent blank.

#### Actions required:

When performing assays in analyser ACCENT-300, there is a probability of **cross-contamination** affecting the tests results: UA – URINE PROTEINS, FERRUM – UA, HDL DIRECT – UA, LDL DIRECT – UA, URINE PROTEINS – UA, ALBUMIN - UA. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_007\_ACCENT-300\_CARRYOVER.

### REFERENCE VALUES <sup>5</sup>

serum / plasma	mg/dl	µmol/l
females	2.5 – 6.8	149 – 405
males	3.6 – 7.7	214 – 458

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 11 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 1.15 mg/dl (68.4 µmol/l).
- Linearity:** up to 32.4 mg/dl (1927.15 µmol/l).
- Specificity / Interferences**

Haemoglobin up to 1.25 g/dl, ascorbate up to 31 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run) n = 10	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	4.60	0.07	1.55
level 2	9.49	0.09	0.94
Reproducibility (day to day) n = 40	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	6.03	0.21	3.50
level 2	9.45	0.21	2.19

### Method comparison

A comparison between uric acid values for samples obtained on **ACCENT-300** (y) and obtained on **ADVIA 1650** (x) using 80 samples gave following results:  
y = 1.0121 x + 0.065 mg/dl;  
R = 0.9968 (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

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- Henry R.J.: Clinical Chemistry, Harper & Row Publishers Inc., New York (1974).