ACCENT-300 LIPASE

Nr kat. **7-309** (PL)

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LITERATURA

- Tietz NW et al. Lipase in serum-the elusive enzyme: An overview. Clin Chem 1993;39:746-756.
- Steinberg WM, Goldstein SS, Davies ND et al. Diagnostic assays in acute pancreatitis. (Review). Ann Intern Med 1985; 102:576-580.
- Leybold A, Junge W. Importance of colipase for the measurement of serum lipase activity. Adv Clin Enzymol 1986;4:60-67.
- Alan H. B. Wu, Tietz Clinical Guide to Laboratory Tests, W.B. Saunders Company, 4th edition, 676 (2006).

Data wydania: 10.2020.

ACCENT-300 LIPASE

Cat. No **7-309** (EN)

INTENDED USE

Diagnostic kit for determination of lipase activity, used in automatic analyser ACCENT-300.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Lipase is a digestive enzyme released into the intestine from the pancreas where it breaks down triglycerides into fatty acids and glycerol prior to absorption. Lipase measurements are used in the diagnosis and treatment of diseases of the pancreas such as acute pancreatitis, obstruction of the pancreatic duct and pancreatic tumours.

METHOD PRINCIPLE

The colorimetric method is based on a lipase specific degradation of a chromogenic substrate. The specific lipase substrate-DGGMR [1,2-o-dilauryl-racglycero-3-glutaric acid-(6'-methylresorufin) ester] is cleaved by the catalytic action of lipase to form 1,2-o-dilauryl-racglycerol and an unstable intermediate, glutaric acid-(6-methyl resorufin) ester. This decomposes spontaneously in alkaline solution to form glutaric acid and methylresorufin. The lipase activity in the specimen is proportional to the production of methylresorufin in the reaction and can be determined photometrically.

REAGENTS

Package	
1-Reagent	1 x 40 ml
2-Reagent	1 x 14 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents, stored on board of the analyzer at 2-10°C are stable for 8 weeks.

Concentrations in the test

1-Reagent	
TAPS [N-Tris(hydroxymethyl)methyl-3-aminopropanesulfonic acid]	100 mM
sodium hydroxide	40 mM
sodium deoxycholate preservative	34 mM
2-Reagent	
tartaric acid	9.5 mM
sodium hydroxide	19 mM
colipase	460 IU/ml
2-propanol	0.65 M
DGGMR [1,2-o-dilauryl-rac-glycero-3-glutaric acid-(6'-methylresorufin)-ester]	0.4 mM

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum, heparinized plasma free from hemolysis. Sample may be stored for up to 5 days at 2-8°C or 24 hours at 20-25°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use. Deionised water is recommended as a reagent blank.

Actions required:

When performing assays in analyser ACCENT-300, there is a probability of **cross-contamination** affecting the tests results: CHOL – LIPASE, LDL DIRECT – LIPASE, TG – LIPASE. To avoid this effect follow the recommendations contained in the instruction 51_03_24_007_ACCENT-300_CARRYOVER.

REFERENCE VALUES ⁴

Normal range	13 – 60 U/l	0.22 – 1.00 μ kat/l
--------------	-------------	-------------------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 8 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser ACCENT-300. Results may vary if a different instrument is used.

- Sensitivity:** 12 U/l (0.2 μ kat/l).
- Linearity:** up to 420 U/l (7.0 μ kat/l).

For higher activity dilute sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

Specificity / Interferences

Haemoglobin up to 0.16 g/dl, ascorbate up to 62 mg/l, bilirubin up to 15 mg/dl and triglycerides up to 750 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 10	Mean [U/l]	SD [U/l]	CV [%]
level 1	33.34	0.83	2.49
level 2	93.60	1.12	1.19
Reproducibility (day to day) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	35.00	0.69	1.96
level 2	94.76	1.24	1.31

Method comparison

A comparison between lipase values determined at ACCENT-300 (y) and at ADVIA 1650 (x) using 29 samples gave following results:

$$y = 0.9768x - 0.8521 \text{ U/l};$$

$$R = 0.9974 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Tietz NW et al. Lipase in serum-the elusive enzyme: An overview. Clin Chem 1993;39:746-756.
- Steinberg WM, Goldstein SS, Davies ND et al. Diagnostic assays in acute pancreatitis. (Review). Ann Intern Med 1985; 102:576-580.
- Leybold A, Junge W. Importance of colipase for the measurement of serum lipase activity. Adv clin Enzymol 1986;4:60-67.
- Alan H. B. Wu, Tietz Clinical Guide to Laboratory Tests, W.B. Saunders Company, 4th edition, 676 (2006).

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