## ACCENT-300 COMPLEMENT C4

Nr kat. 7-313 (PL)

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### LITERATURA

- Tietz Textbook of Clinical Chemistry, W.B. Saunders, Philadelphia, (1994).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 2262, 2006

Data wydania: 10. 2020.

## ACCENT-300 COMPLEMENT C4

Cat. No 7-313 (EN)

### INTENDED USE

Diagnostic kit for determination of complement C4 concentration, used in automatic analyser ACCENT-300. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Complement is a group of 20 immunologically distinct proteins present in blood and tissues. They are able to interact with antigen-antibody complexes, with each other and with cell membranes, in a complex way intended at destroying viruses and bacteria. They are synthesised in liver and are present in serum as functionally inactive molecules. They are activated by antigen-antibody complexes. C4 complement is a  $\alpha$ -glycoprotein of 3 subunits. It is an acute phase reactant whose levels are increased during the acute phase. Low levels are found in immune complex diseases and in inherited angioedema, while C3 complement levels are normal.

### METHOD PRINCIPLE

The complement C4 present in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at  $\lambda=340$  nm is proportional to complement C4 concentration in the sample.

### REAGENTS

**Package**  
 1-Reagent 1 x 40 ml  
 2-Reagent 1 x 9 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

### Reagent components

**1-Reagent:** imidazole buffer, polymer, inorganic salt, buffer, preservative.

**2-Reagent:** anti-human complement C4 antiserum, inorganic salt, buffer, preservative.

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum or plasma.  
 Serum should be separated from red blood cells as soon as possible after blood collection. If the test cannot be done immediately, the sample should be stored at -70°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.  
 For reagent blank 0.9% NaCl is recommended.

### REFERENCE VALUES <sup>2</sup>

adults	0.1 – 0.4 g/l
--------	---------------

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL III (Cat. No 4-291) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using an automatic analyser Hitachi. Results may vary if a different instrument is used.

- LOQ:** 0.002 g/l
- Linearity:** up to 0.825 g/l.
- Specificity / Interferences**

Hemoglobin up to 0.32 g/dl, bilirubin up to 22 mg/dl triglycerides up to 155.5 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.147	0.001	0.77
level 2	0.321	0.003	0.76
level 3	0.462	0.003	0.64
Reproducibility (day to day) n = 60	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.137	0.005	3.48
level 2	0.321	0.008	2.52
level 3	0.472	0.010	2.08

### Method comparison

A comparison between CORMAY reagent (y) and commercially available assay (x) using at least 30 samples gave following results:

$$y = 0.9216x + 0.2927 \text{ g/l};$$

$$R = 0.906 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- Tietz Textbook of Clinical Chemistry, W.B. Saunders, Philadelphia, (1994).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 2262, 2006

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