ACCENT-300 ASAT

Nr kat. 7-314 (PL)

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Cat. No **7-314** (EN)

INTENDED USE

Diagnostic kit for determination of aspartate aminotransferase activity intended to use in automatic analyser ACCENT-300. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Aspartate aminotransferase (ASAT, AST, GOT) is an enzyme participated in amino acids metabolism. ASAT is found in all tissues but particularly high level of ASAT is observed in heart muscle, skeletal muscle, liver and kidney. This is why elevated ASAT serum level is marker of myocardial infarction and kidney, liver or skeletal muscle injury.

METHOD PRINCIPLE

Optimized, modified method according to International Federation of Clinical Chemistry (IFCC), without pyridoxal phosphate.

The rate of absorbance changing at $\lambda=340$ nm is directly proportional to aspartate aminotransferase activity.

REAGENTS

Package
 1-Reagent 4 x 48 ml
 2-Reagent 1 x 48 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 12 weeks on board the analyser at 2-10°C.

Concentrations in the test

Tris (pH 7.8)	80 mmol/l
L-aspartate	240 mmol/l
MDH	> 10 μ kat/l
LDH	> 20 μ kat/l
2-oxoglutarate	15 mmol/l
NADH	0.18 mmol/l
sodium hydroxide	< 1 %

Warnings and notes

- Protect from light and avoid contamination!
- The reagents contain (< 0.1%) sodium azide as a preservative. Avoid contact with skin and mucous membranes.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Warning.

 H315 Causes skin irritation.
 H319 Causes serious eye irritation.
 P280 Wear protective gloves/protective clothing/eye protection/face protection.
 P302+P352 IF ON SKIN: Wash with plenty of soap and water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

SPECIMEN

Serum, heparinized or EDTA plasma free from hemolysis. Do not use heparine ammonium salt. Hemolysis should be avoided, since ASAT activity in erythrocytes is 10 times higher than in normal serum. Do not freeze the samples. ASAT activity remains stable in specimen up to 1 day at 15-25°C or up to 4 days at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
 For reagent blank deionized water is recommended.

REFERENCE VALUES ⁶

serum / plasma	37°C	
female	up to 31 U/l	up to 0.518 μ kat/l
male	up to 37 U/l	up to 0.618 μ kat/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples. For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. The calibration curve should be prepared every 12 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyzer ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

Sensitivity

9 U/l (0.150 μ kat/l).

Linearity

up to 1100 U/l (18.4 μ kat/l).

Specificity / Interferences

Haemoglobin up to 0.16 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	44.30	1.87	4.23
level 2	178.10	4.15	2.33
Reproducibility (day to day) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	35.13	1.63	4.63
level 2	141.78	4.01	2.83

Method comparison

A comparison between ASAT values determined at ACCENT-300 (y) and at COBAS INTEGRA 400 (x) using 76 samples gave following results:

$$y = 1.0953 x - 1.3024 \text{ U/l};$$

$$R = 0.9998 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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