ACCENT-300 CERULOPLASMIN

Nr kat. 7-318 (PL)

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(EN)

INTENDED USE

Diagnostic kit for determination of ceruloplasmin concentration, used in automatic analyser ACCENT-300.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Ceruloplasmin is an α_2 -glycoprotein containing 6-7 atoms of copper per molecule. It has been considered for a long time as a copper transporter but since recently, it has been shown to be a serum ferroxidase playing a major role in oxidising iron (II) to iron (III) in serum and at the cell surface, thereby regulating its binding by transferrin.

Ceruloplasmin is a late acute phase reactant. Low levels are found in malnutrition, nephrosis, severe liver diseases such as primary biliary cirrhosis and in Wilson's disease, an autosomal recessive defect in the regulation of copper metabolism.

METHOD PRINCIPLE

The ceruloplasmin presents in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to ceruloplasmin concentration in the sample.

REAGENTS

Package

1-Reagent 1 x 40 ml
 2-Reagent 1 x 9 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: MES buffer, polymer, inorganic salt, preservative.

2-Reagent: anti-human ceruloplasmin antiserum, buffer, inorganic salt, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store reagents closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN ⁷

Serum.

Specimen without lipemia or hemolysis is recommended. Samples remain stable for 3 days at 2-8°C or 4 weeks at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ⁷

age	g/l	mg/dl
1 day – 3 months	0.05 – 0.18	5 – 18
6 – 12 months	0.33 – 0.43	33 – 43
1 – 7 years	0.24 – 0.56	24 – 56
> 7 years	0.18 – 0.45	18 – 45

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers systems the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended.

The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using an automatic analyser Hitachi. Results may vary if a different instrument or a manual procedure is used.

LOQ (Limit of Quantitation):

0.948 mg/dl (0.00948 g/l)

Linearity: up to 105.12 mg/dl (1.0512 g/l).

Specificity / Interferences

Hemoglobin up to 0.19 g/dl, bilirubin up to 22 mg/dl triglycerides up to 1669,5 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

Precision

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.25	0.002	0.6
level 2	0.40	0.003	0.6
level 3	0.52	0.004	0.7

Reproducibility (day to day) n = 60	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.19	0.006	3.0
level 2	0.36	0.006	1.6
level 3	0.49	0.006	1.3

Method comparison

A comparison between CORMAY reagent (y) and commercially available assay (x) using at least 30 samples gave following results:

$y = x + 26.28$ mg/dl;

$R = 0.986$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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Parameters

No.	42	Prim.Wave.	340
Test	CERUL	Sec.Wave.	700
Method	Endpoint	Sample Vol.	5
Direction	Ascend	R1 Vol.	250
Unit	g/l	R2 Vol.	50
Decimals	2	Line. Limit	
Incubation	25	Antigen Check	
Reaction	-1 25	Substrat	0
R1 Blank		Mix. R Blank	
Lower	0	Lower	0
Upper	0	Upper	0
Response		Linearity	
Lower	-2.5	Lower	
Upper	2.5	Upper	
Sample Vol.	45	Full Name	CERUL
Dilution	5	Print No.	42

Calibration

Rule	Logistic 5P
K Factor	0
Replicates	1
Interval	0
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0 2.5
Coefficient	0
Difference	
Non-linear SD	0

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