ACCENT-300 MYOGLOBIN

Nr kat. 7-326 (PL)

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Cat.No 7-326 (EN)

INTENDED USE

Diagnostic kit for determination of myoglobin concentration intended to use in automatic analyzer ACCENT-300.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Myoglobin (Mb) is a hemo-protein present in cardiac and skeletal muscle cells and is released into blood circulation when these cells are damaged. The determination of serum Mb level is useful in the diagnosis of myocardial infarction, muscular dystrophy, myositis and myopathy, and also for the assessment of treatment and disease prognosis.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between Mb in a sample and anti-Mb antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of Mb in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package	
1-Reagent	1 x 31 ml
2-Reagent	1 x 11 ml

The reagents when stored at 2-10°C are stable up to expiry date printed on the package. The reagents stored on board of the analyzer at 2-10°C are stable for 8 weeks (Hitachi 911).

Concentrations in the test

suspension of latex particles sensitized with anti-Mb (rabbit) antibodies (pH 7.3) glycine buffer solution (pH 9.0) preservative 0.12 w/v%

Warnings and notes

- Protect from light and avoid contamination!
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Care should be taken not to interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum or plasma (Na-EDTA, K-EDTA, Na-Heparin, Li-Heparin, citric acid).

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
 For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ³

serum, plasma	< 70 ng/ml
---------------	------------

It is recommended for each laboratory to establish its own reference ranges for local population. Diagnosis should only be made after taking clinical symptoms and the results of other tests into consideration.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL II (Cat. No 4-290) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY MYOGLOBIN CALBRATORS kit (Cat. No 4-279) is recommended.

The calibration curve should be prepared every 8 weeks (Hitachi 911), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using an automatic analyzer Hitachi 917. Results may vary if a different instrument is used.

- Analytical range:** 20 – 1000 ng/ml.

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

- Specificity / Interferences**

Haemoglobin up to 0.96 g/dl, ascorbate up to 500 mg/l, bilirubin up to 62 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

- Precision**

	Repeatability (run to run) n = 21		
	Mean [ng/ml]	SD [ng/ml]	CV [%]
level 1	89.38	1.02	1.1
level 2	323.97	1.97	0.6
	Reproducibility (day to day) n = 21		
	Mean [ng/ml]	SD [ng/ml]	CV [%]
level 1	77.0	1.703	2.2
level 2	364.4	6.004	1.6

- Method comparison**

A comparison between CORMAY reagent (y) and another commercially available assay (x) using 50 samples gave following results:

$$y = 1.47x - 19.10 \text{ ng/ml;}$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Parameters

No.	46	Prim.Wave.	578
Test	MYO	Sec.Wave.	0
Method	Kinetic	Sample Vol.	6
Direction	Ascend	R1 Vol.	180
Unit	ng/ml	R2 Vol.	60
Decimals	0	Line. Limit	
Incubation	15	Antigen Check	
Reaction	10 20	Substrat	0
R1 Blank		Mix. R Blank	
Lower	0	Lower	0
Upper	0	Upper	0
Response		Linearity	
Lower	-2.5	Lower	
Upper	2.5	Upper	
Sample Vol.	45	Full Name	Myoglobin
Dilution	5	Print No.	46

Calibration

Rule	Logistic 5P
K Factor	0
Replicates	1
Interval	0
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0 2.5
Coefficient Difference	0
Non-linear SD	0

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