## ACCENT-300 TOTAL PROTEIN

Nr kat. 7-336

(PL)

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Cat. No. **7-336** (EN)

### INTENDED USE

Diagnostic kit for determination of total protein concentration, intended to use in automatic analyser ACCENT-300.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Most serum proteins except gamma globulins and hemoglobin are synthesized in the liver. Proteins participate in transport, catalysis and coagulation, act as hormones and receptors, antigens and antibodies, regulate osmotic pressure and play structural functions. Correct serum level of total protein depends mainly on balance between synthesis and degradation of albumin and immunoglobulins. Total protein level abnormalities are caused usually by dehydration, liver or kidney disease and starvation.

### METHOD PRINCIPLE

The method is based on the biuret reaction. Protein forms the coloured complex with cupric ions in alkaline medium. The colour intensity is proportional to the protein concentration.

### REAGENTS

#### Package

1-Reagent 2 x 40 ml

The reagent when stored at 2-8°C is stable up to expiry date printed on the package. The reagents are stable for 5 weeks on board the analyser at 2-10°C.

### Concentrations in the test

copper sulfate ≤ 14 mmol/l  
sodium-potassium tartrate ≤ 36 mmol/l  
potassium iodide ≤ 36 mmol/l  
sodium hydroxide ≤ 570 mmol/l

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

#### Ingredients:

1-Reagent contains sodium hydroxide and copper (II) sulfate pentahydrate.

#### Danger

H314 Causes severe skin burns and eye damage.  
H411 Toxic to aquatic life with long lasting effects.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P273 Avoid release to the environment.

### SPECIMEN

Serum or plasma. Nonhemolyzed, nonlipaemic, fresh serum is recommended.

Recommended anticoagulants: EDTA, heparine lithium, sodium or ammonium salt.

Serum should be separated from red blood cells as soon as possible after blood collection.

Serum and plasma can be stored up to 3 days at 2-8°C or 6 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent is ready to use.

For reagent blank deionized water is recommended.

### REFERENCE VALUES <sup>7</sup>

serum	g/dl	g/l
premature 1 day	3.4 – 5.0	34 – 50
children 1day – 4 week	4.6 – 6.8	46 – 68
2 – 12 month	4.8 – 7.6	48 – 76
≥1 year	6.0 – 8.0	60 – 80
adults	6.6 – 8.7	66 – 87

Results obtained from plasma analysis might be slightly elevated due to fibrinogen presence<sup>8</sup>:

Sample	Protein concentration increase from serum to plasma
Blood donors	+ 0,25 g/dl
Nonhospitalized patients	+ 0,36 g/dl
Hospitalized patients	+ 0,46 g/dl
Hospitalized patients with CRP<50 mg/dl	+ 0,66 g/dl

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176). Deionized water should be used as a calibrator 0.

The calibration curve should be prepared every 4 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

▪ **Sensitivity:** 0.27 g/dl (2.7 g/l).

▪ **Linearity:** up to 16.50 g/dl (165 g/l)

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

### Specificity / Interferences

Haemoglobin up to 0.31 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 950 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run) n = 10	Mean [g/dl]	SD [g/dl]	CV [%]
level 1	4.36	0.04	0.88
level 2	6.92	0.03	0.49
Reproducibility (day to day) n = 20	Mean [g/dl]	SD [g/dl]	CV [%]
level 1	4.38	0.04	0.93
level 2	6.91	0.04	0.64

### Method comparison

A comparison between total protein values determined at ACCENT-300 (y) and at ADVIA 1650 (x) using 21 samples gave following results:

$$y = 0.910 x + 0.573 \text{ g/dl};$$

$$R = 0.957 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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PROGRAM NA ANALIZATOR / APPLICATION for / АДАПТАЦИЯ для:

### Parameters

No.	4	Prim.Wave.	546
Test	TP	Sec.Wave.	700
Method	Endpoint	Sample Vol.	3
Direction	Ascend	R1 Vol.	300
Unit	g/dl	R2 Vol.	
Decimals	2	Line. Limit	
Incubation	0	Antigen Check	
Reaction	0   50	Substrat	0
<b>R1 Blank</b>		<b>Mix. R Blank</b>	
Lower	0	Lower	0
Upper	0	Upper	0
<b>Response</b>		<b>Linearity</b>	
Lower	-2.5	Lower	0.27
Upper	2.5	Upper	16.50
Sample Vol.	45	Full Name	Total Protein
Dilution	5	Print No.	4

### Calibration

Rule	Two-Point Linear
K Factor	0
Replicates	3
Interval	28
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0   2.5
Coefficient	0
Difference	
Non-linear SD	0

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