ACCENT-300 TG

Nr kat. 7-353	(PL)	peroksydaza (POD) > 20,00 µkat/l lipaza lipoproteinowa (LPL) > 16,67 µkat/l konserwant
----------------------	-------------	--

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Cat. No **7-353** (EN)

INTENDED USE

Diagnostic kit for determination of triglycerides concentration, used in automatic analyser ACCENT-300.



The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Triglycerides are built of glycerol molecule esterified with three fatty acids molecules. Triglycerides are delivered with food or are synthesized endogenously in liver. Triglycerides stored in adipose tissue constitute a reserve of energy. Elevated triglycerides serum level is a risk factor of atherosclerosis. Triglycerides measurement is useful for hyperlipidemia diagnosis and treatment or for estimation of atherosclerosis progression.

METHOD PRINCIPLE

Colorimetric, enzymatic method with glycerophosphate oxidase.

The colour intensity is proportional to the triglycerides concentration.

REAGENTS

Package

1-Reagent 4 x 40 ml
 2-Reagent 1 x 40 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 11 weeks.

Concentrations in the test

buffer PIPES (pH 7.0)	40 mmol/l
4-aminoantipyrine (4-AA)	0.4 mmol/l
ATP	1.5 mmol/l
Mg ²⁺	1.6 mmol/l
ADPS	0.6 mmol/l
glycerol kinase (GK)	> 66.67 µkat/l
glycerol-3-phosphate oxidase (GPO)	> 60.00 µkat/l
peroxidase (POD)	> 20.00 µkat/l
lipoprotein lipase (LPL)	> 16.67 µkat/l
preservative	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum, EDTA or heparinized plasma (recommended: heparine lithium, sodium or ammonium salt) free from hemolysis.

Blood should be collected only if the patient has been fasting for minimum of 12 hours. Before blood collection patient should stay in rest position for about 30 minutes. Venous blood is recommended for triglycerides measurement.

Plasma triglycerides values have been reported to be 2% to 4% lower than serum triglycerides values.

Serum and plasma can be stored up to 3 days at 2-8°C or 3 months at -20°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

Deionised water is recommended as a reagent blank.

Actions required:

When performing assays in analyser ACCENT-300, there is a probability of **cross-contamination** affecting the tests results: TG – LIPASE, TG – LIPASE II GEN, HDL DIRECT – TG, LDL DIRECT - TG. To avoid this effect follow the recommendations contained in the instruction: 51_03_24_007_ACCENT-300_CARRYOVER.

REFERENCE VALUES²

serum, plasma	< 150 mg/dl < 1.7 mmol/l
---------------	-----------------------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, CORMAY LIPID CONTROL 1 (Cat. No 5-179) and CORMAY LIPID CONTROL 2 (Cat. No 5-180) or CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and/or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 11 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 2.5 mg/dl (0.028 mmol/l).
- Linearity:** up to 2540 mg/dl (28.7 mmol/l).

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 2.5 g/dl, bilirubin up to 20 mg/dl and ascorbate up to 62 mg/l do not interfere with the test.

Precision

Repeatability (run to run) n = 10	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	82.61	0.97	1.17
level 2	182.55	1.55	0.85
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	100.90	3.90	3.86
level 2	273.83	6.84	2.50

Method comparison

A comparison between triglycerides values determined at ACCENT-300 (y) and at ADVIA 1650 (x) using 30 samples gave following results:

$$y = 0.9754x - 0.9606 \text{ mg/dl;}$$

$$R = 0.9991 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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PROGRAM NA ANALIZATOR / APPLICATION / АДАПТАЦИЯ:

Parameters

No.	7	Prim.Wave.	546
Test	TG	Sec.Wave.	700
Method	Endpoint	Sample Vol.	3
Direction	Ascend	R1 Vol.	200
Unit	mg/dl	R2 Vol.	50
Decimals	1	Line. Limit	
Incubation	10	Antigen Check	
Reaction	-2 35	Substrat	0
R1 Blank		Mix. R Blank	
Lower	0	Lower	0
Upper	0	Upper	0
Response		Linearity	
Lower	-2.5	Lower	2.5
Upper	2.5	Upper	2540
Sample Vol.	45	Full Name	Triglycerides
Dilution	5	Print No.	7

Calibration

Rule	Two Point Linear / One Point Linear
K Factor	0
Replicates	3
Interval	77
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0 2.5
Coefficient	0
Difference	
Non-linear SD	0