ACCENT-300 CREA ENZYMATYC

Nr kat. 7-377 (PL)

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LITERATURA

- Newman DJ, Pnce CP, Renal function and nitrogen metabolites. In: Burtis CA, Ashwood ER, Bruns DE, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th ed, St. Louis: W.B Saunders Company; 2006. p. 797-801.
- Alan H.B. Wu. editor. Tietz Clinical Guide to Laboratory Tests, 4th ed. St. Louis: W.B Saunders Company; 2006. p.316.
- Mazzachi BC, Peake M, Erhardt V, Reference range and method comparison for enzymatic and Jaffe Creatinine assays in plasma and serum and early morning urine. Clin Lab 2000; 46: 53-5.
- Susumu Osawa, Medical Technology 1982, Vol. 10, No 7, 575-579.
- Minoru Konno, Medical Technology 1984, Vol. 12, No 3, 270-276.
- Schlebusch H, Liappis N, Klein G. Ultrasensitive CRP and Creatinine: Reference intervals from infancy to childhood. Clin Chem Lab Med. 2001; 39 Special supplement pp S1-S448; May 2001. PO-T042.

Data wydania: 10.2020

ACCENT-300 CREA ENZYMATIC

Cat. No 7-377

(EN)

INTENDED USE

Diagnostic kit for determination of creatinine concentration, used in automatic analyser ACCENT-300.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Creatinine is a product of creatine nonenzymatic dehydration in skeletal muscle. The amount of creatinine generated and excreted by kidney is proportional to muscle mass and usually is higher in men than women. Daily creatinine generation remains fairly constant, with the exception of crushing injury or degenerative diseases that cause massive damage to muscle. Creatinine blood and urine level depends on glomerular filtration so creatinine clearance is excellent index of renal function.

METHOD PRINCIPLE

Enzymatic, colorimetric method.

The colour intensity measured at 546 nm is proportional to the creatinine concentration.

REAGENTS

Package

1-Reagent 2 x 40 ml
2-Reagent 1 x 27 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 12 weeks.

Components and concentrations

1-Reagent

Good's buffer ≤ 5%
creatininase ≤ 5%
N-ethyl-N-(3-sulfoethyl)-3-methylaniline (ESPMT) ≤ 5%
sarcosine oxidase ≤ 0.01%
ascorbate oxidase ≤ 1%
detergents, stabilizers and preservatives

2-Reagent

Good's buffer ≤ 5%
creatininase ≤ 1%
peroxidase ≤ 5%
4-aminoantipyrine (4-AA) ≤ 0.01%
stabilizers and preservatives

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum and urine.

Urine preparation: before analysis urine sample should be diluted with 0,9% NaCl 2-10 times. Multiply the result by dilution factor.

Serum can be stored up to 1 day at 2-8°C. For longer storage samples should be frozen at -20°C.

Urine can be stored up to 1 day at 20-25°C, 4 days at 2-8°C. For longer storage samples should be frozen at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

Deionised water is recommended as a reagent blank.

Actions required:

When performing assays in analyser ACCENT-300, there is a probability of **cross-contamination** affecting the tests results: CREA ENZYMATIC – URINE PROTEINS, CK-MB – CREA ENZYMATIC, GLUCOSE – CREA ENZYMATIC. To avoid this effect follow the recommendations contained in the instruction 51_03_24_007_ACCENT-300_CARRYOVER.

REFERENCE VALUES ^{2,3}

serum / plasma	mg/dl	µmol/l
newborns	0.3 – 1.0	26.5 – 88.4
infants	0.2 – 0.4	17.7 – 35.4
children	0.2 – 0.8	17.7 – 70.7
female	0.5 – 1.0	44.2 – 88.4
male	0.7 – 1.2	61.9 – 106.1
urine (morning)	mg/dl	mmol/l
female	29 – 226	2.56 – 20.0
male	40 – 278	3.54 – 24.6

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) or LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 7 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

▪ **Sensitivity:** 0.08 mg/dl (7.07 µmol/l).

▪ **Linearity:** up to 25 mg/dl (2210 µmol/l).

Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbic acid up to 62 mg/l, bilirubin up to 20 mg/dl, triglycerides up to 1000 mg/dl and creatine up to 20 mg/dl do not interfere with the test both at serum and urine.

Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	1.29	0.01	1.11
level 2	7.88	0.03	0.37
Repeatability (run to run) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	1.30	0.08	6.07
level 2	7.94	0.46	5.75

Method comparison

A comparison between CORMAY reagent (y) and another commercially available assay (x) using 47 serum samples gave following results:

$y = 1.0673x - 0.0752$ mg/dl;

$R = 0.9986$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Newman DJ, Pnce CP, Renal function and nitrogen metabolites. In: Burtis CA, Ashwood ER, Brunts DE, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th ed, St. Louis: W.B Saunders Company; 2006. p. 797-801.
- Alan H.B. Wu. editor. Tietz Clinical Guide to Laboratory Tests, 4th ed. St. Louis: W.B Saunders Company; 2006, p.316.
- Mazzachi BC, Peake M, Erhardt V, Reference range and method comparison for enzymatic and Jaffe Creatinine assays in plasma and serum and early morning urine. Clin Lab 2000; 46: 53-5.
- Susumu Osawa, Medical Technology 1982, Vol. 10, No 7, 575-579.
- Minoru Konno, Medical Technology 1984, Vol. 12, No 3, 270-276.
- Schlebusch H, Liappis N, Klein G. Ultrasensitive CRP and Creatinine: Reference intervals from infancy to childhood. Clin Chem Lab Med. 2001; 39 Special supplement pp S1-S448; May 2001. PO-T042.

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- Newman DJ, Pnce CP, Renal function and nitrogen metabolites. In: Burtis CA, Ashwood ER, Bruns DE, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th ed, St. Louis: W.B Saunders Company; 2006. p. 797-801.
- Alan H.B. Wu. editor. Tietz Clinical Guide to Laboratory Tests, 4th ed. St. Louis: W.B Saunders Company; 2006, p.316.
- Mazzachi BC, Peake M, Erhardt V, Reference range and method comparison for enzymatic and Jaffe Creatinine assays in plasma and serum and early morning urine. Clin Lab 2000; 46: 53-5.
- Susumu Osawa, Medical Technology 1982, Vol. 10, No 7, 575-579.
- Minoru Konno, Medical Technology 1984, Vol. 12, No 3, 270-276.
- Schlebusch H, Liappis N, Klein G. Ultrasensitive CRP and Creatinine: Reference intervals from infancy to childhood. Clin Chem Lab Med. 2001; 39 Special supplement pp S1-S448; May 2001. PO-T042.

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PROGRAM NA ANALIZATORY / APPLICATION for / АДАПТАЦИЯ для:

Parameters

No.	77	Prim.Wave.	546
Test	CREA ENZ	Sec.Wave.	
Method	Endpoint	Sample Vol.	7
Direction	Ascend	R1 Vol.	180
Unit	mg/dl	R2 Vol.	60
Decimals	2	Line. Limit	20
Incubation	15	Antigen Check	
Reaction	-2 25	Substrat	0
R1 Blank		Mix. R Blank	
Lower	0	Lower	0
Upper	0	Upper	0
Response		Linearity	
Lower	-2.5	Lower	0.08
Upper	2.5	Upper	25
Sample Vol.	45	Full Name	Crea Enzymatic
Dilution	5	Print No.	77

Calibration

Rule	One Point Linear
K Factor	0
Replicates	3
Interval	49
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0 2.5
Coefficient	0
Difference	
Non-linear SD	0