## ACCENT-300 LDL DIRECT

Nr kat. **7-380** (PL)

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Cat. No **7-380** (EN)

### INTENDED USE

Diagnostic kit for determination of  $\alpha$ -amylase activity, used in automatic analyser ACCENT-300.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Plasma lipoproteins are spherical particles containing varying amounts of cholesterol, triglycerides, phospholipids and proteins. The relative protein and lipid determine the density of these lipoproteins and provide the basis on which to begin their classification. The classes are: chylomicron, very-low-density lipoprotein (VLDL), low-density-lipoprotein (LDL) and high-density lipoprotein (HDL).

LDL are synthesized in the liver by the action of various lipolytic enzymes on triglyceride rich VLDL. LDL-cholesterol concentration is considered to be the most important clinical predictor, of all single parameters, with respect to coronary atherosclerosis.

Accurate measurement of LDL-cholesterol is of vital importance in therapies which focus on lipid reduction to prevent atherosclerosis or reduce its progress and to avoid plaque rupture.

### METHOD PRINCIPLE

The assay is a homogenous method for directly measuring LDL-cholesterol concentrations in serum or plasma, without the need for any off-line pretreatment or centrifugation steps.

Liquid selective detergent method.

The method is in a two reagents format and depends on the properties of a unique detergent. This detergent (Reagent 1) solubilizes only the non LDL particles (HDL, VLDL, CM). The cholesterol released is consumed by cholesterol esterase and cholesterol oxidase in a non color forming reaction. A second detergent (Reagent 2) solubilizes the remaining LDL particles and chromogenic coupler allows for color formulation. The enzyme reaction with LDL-cholesterol in the presence of the coupler produces color that is proportional to the amount of LDL cholesterol present in the sample.

### REAGENTS

#### Package

1-Reagent 2 x 39 ml  
2-Reagent 1 x 26 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 12 weeks.

### Concentrations in the test

#### 1-Reagent

Buffer < 1.0 %  
Detergent 1 < 1.0 %  
Cholesterol esterase (*Pseudomonas* spp.) < 1500 U/l  
Cholesterol oxidase (*Cellulomonas* spp.) < 1500 U/l  
Peroxidase (horseradish) < 1300 ppG U/l  
4-aminoantipyrine < 0.1 %  
Ascorbic Acid Oxidase (*Curcubita* spp.) < 3000 U/l  
Preservative

### 2-Reagent

Buffer < 1.0 %  
Detergent 2 < 1.0 %  
N,N-bis(sulfobutyl)-toluidine, disodium (DSBmT) < 1.0 mM  
Preservative

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze reagents.

### SPECIMEN

Serum, sodium heparinized or EDTA plasma.

Anticoagulants containing citrate should not be used.

Serum: Collect whole blood by venipuncture and allow to clot. Centrifuge and remove the serum as soon as possible after collection (within 3 hours).

Plasma: Specimens may be collected in EDTA or sodium heparin. Centrifuge and remove the plasma as soon as possible after collection (within 3 hours).

If not analysed promptly, specimens may be stored at 2-8°C for up to 5 days. If specimens need to be stored for more than 5 days, they may be preserved at -80°C. Samples may be frozen once.

Nevertheless it is recommended to perform the assay with freshly collected samples.

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

Deionized water is recommended as a reagent blank.

#### Actions required:

When performing assays in analyser ACCENT-300, there is a probability of **cross-contamination** affecting the tests results: LDL DIRECT – TG, LDL DIRECT – LIPASE, LDL DIRECT – LIPASE II GEN, LDL DIRECT – UA, LDL DIRECT – UA PLUS. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_007\_ACCENT-300\_CARRYOVER.

### REFERENCE VALUES <sup>7</sup>

NCEP* Classification	Adults	
	mg/dl	mmol/l
Optional	< 100	< 2.59
Near optional	< 130	< 3.37
Bordeline high	130-159	3.37-4.12
High	160-189	4.14-4.90
Very high	≥ 190	≥ 4.92

\* National Cholesterol Education Program

As LDL cholesterol is affected by a number of factors such as smoking, exercise, hormones, age and sex, each laboratory should establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, CORMAY LIPID CONTROL 1 (Cat. No 5-179) and CORMAY LIPID CONTROL 2

(Cat. No 5-180) or CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers the CORMAY HDL/LDL CALIBRATOR (Cat. No 5-178) is recommended. The calibration curve should be prepared every 12 weeks with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity** : 5.2 mg/dl (0.135 mmol/l).
- Linearity**: up to 700 mg/dl (18.13 mmol/l).

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Triglycerides up to 1293 mg/dl, bilirubin conjugated up to 20 mg/dl, bilirubin total up to 20 mg/dl, haemoglobin up to 0.5 g/dl, ascorbic acid up to 500 mg/l and gamma-globulins up to 5000 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run) n = 10	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	142.72	1.74	1.22
level 2	175.31	2.73	1.56
Reproducibility (day to day) n = 20	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	150.68	2.83	1.88
level 2	180.97	6.84	3.78

#### Method comparison

A comparison between LDL cholesterol values for samples obtained on **ACCENT-300** (y) and obtained on **ADVIA 1650** (x) using 63 samples gave following results:

$y = 1.1182x - 9.305$  mg/dl;

$R = 0.984$  (R - correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

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