## A-400 IgA

Nr kat. 7-402 (PL)

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Cat. No **7-402** (EN)

### INTENDED USE

Diagnostic kit for determination of IgA concentration used in automatic analysers BS-400 and BS-480.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Immunoglobulins (Igs) are the instrumental proteins of immunity. Immunity is a property of the lymphoid system which is made of organs (spleen, thymus, bone marrow) and of cells (lymphocytes). Circulating immunoglobulins are secreted in the blood by B lymphocytes and they thereby export far-away the specific biological functions of humoral immunity. Immunoglobulin A (IgA) is the major Ig found in secretions, playing a major role in the protection of the respiratory, genitourinary, and gastrointestinal tracts against infection.

### METHOD PRINCIPLE

The IgA present in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum is proportional to IgA concentration in the sample.

### REAGENTS

**Package**  
1-Reagent 1 x 40 ml  
2-Reagent 1 x 9 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

### Reagent components

**1-Reagent:** tricine buffer, polymer, inorganic salt, preservative.

**2-Reagent:** anti-human IgA antiserum, buffer, inorganic salt, preservative.

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN <sup>4</sup>

Serum.  
Specimen without lipemia or hemolysis is recommended.

Specimen can be stored up to 3 days at 2-8°C or up to 6 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.  
0.9% NaCl is recommended as a reagent blank.

### REFERENCE VALUES <sup>4</sup>

adults	0.40 – 3.50 g/l
children (1 year – 12 years)	0.15 – 2.50 g/l
children (1 month – 12 months)	0.20 – 0.90 g/l

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended. 0.9% should be used as a calibrator 0.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers: Hitachi or BS-480. Results may vary if a different instrument is used.

- Measurement range:** 0.02 g/l to 8 g/l.

- Specificity / Interferences**

Haemoglobin up to 0.19 g/dl, bilirubin up to 22 mg/dl, triglycerides up to 155.5 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 0.8 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	1.30	0.02	1.13
level 2	2.60	0.07	2.67
level 3	3.90	0.17	4.29
Reproducibility (day to day) n = 66	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.95	0.08	8.75
level 2	2.27	0.09	4.17
level 3	3.54	0.10	3.00

- Method comparison**

A comparison between IgA values determined at **BS-480** (y) and at **ADVIA 1650** (x) using 46 samples gave following results:

$$y = 0.9738x - 0.0683 \text{ g/l}$$

$$R = 0.994 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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#### • BS-480

Chem <b>IgA</b>		No. <b>030</b>		Sample Type <b>SERUM</b>	
Chemistry <b>Immunoglobulin A</b>		Print name <b>IgA</b>			
Reaction Type <b>Endpoint</b>		Reaction Direction <b>Increase</b>			
Pri Wave <b>570</b>		Sec Wave <b>800</b>			
Unit <b>g/L</b>		Decimal <b>0.01</b>			
Blank Time <b>48 49</b>		Reaction Time <b>81 82</b>			
Standard <b>2.5</b> μL		Reagent Vol <b>200</b> μL		Diluent	
Decreased <b>2.5</b> μL		R1 <b>200</b> μL		Diluent	
Increased		R2 <b>40</b> μL		μL	
Sample Blank <input type="checkbox"/>		R3		μL	
Auto Retun <input type="checkbox"/>		R4		μL	
Linearity Range (Standard) <b>0.08 5.0</b>		Linearity Limit			
Linearity Range (Decreased)		Substrate Depletion			
Linearity Range (Increased)		Mixed Blank Abs <b>-33000 33000</b>			
R1 Blank Abs <b>-33000 33000</b>		Uncapping Time <b>56</b> Day(s)			
Blank Response <b>-33000 33000</b>		Reagent Alarm Limit			
Twin Chemistry		<input type="checkbox"/> Enzyme Linear Extension			
<input type="checkbox"/> Prozone Check		<input type="checkbox"/> Rate Check		<input type="checkbox"/> Antigen Addition	
Q1 0		Q2 0		Q3 0 Q4 0	
PC 0		ABS 0			
<b>Calibration Settings</b>		<b>Auto Calibration</b>			
Math Model <b>Logit-Log 5P</b>		<input type="checkbox"/> Bottle Changed			
Factor		<input type="checkbox"/> Lot Changed			
Replicates <b>1</b>		<input type="checkbox"/> Cal Time			
<b>Acceptance Limits</b>					
Cal Time <b>672</b> Hour		SD			
Slope Diff		Repeatability			
Sensitivity					
Deter Coeff					

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