A-400 IgG

Nr kat. 7-403 (PL)

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Cat. No **7-403** (EN)

INTENDED USE

Diagnostic kit for determination of IgG concentration used in automatic analysers BS-400 and BS-480.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Immunoglobulins (Igs) are the instrumental proteins of immunity. Immunity is a property of the lymphoid system which is made of organs (spleen, thymus, bone marrow) and of cells (lymphocytes). Circulating immunoglobulins are secreted in the blood by B lymphocytes and they thereby export far-away the specific biological functions of humoral immunity.

Immunoglobulin G (IgG) - the predominant serum immunoglobulin (75 % of Igs)- is of particular importance in the body's long-term defense against infection; IgG deficiency is associated with recurrent and occasionally severe pyogenic infections. IgG serum levels are increased in response to chronic or recurrent infections or autoimmune diseases.

METHOD PRINCIPLE

The IgG present in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum is proportional to IgG concentration in the sample.

REAGENTS

Package

1-Reagent 1 x 40 ml
 2-Reagent 1 x 16 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: tricine buffer, polymer, inorganic salt, preservative.

2-Reagent: anti-human IgG antiserum, buffer, inorganic salt, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.

- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN ⁴

Serum. Nonhemolyzed, nonlipaemic, fresh serum is recommended.

Serum can be stored up to 3 days at 2-8°C or up to 6 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ⁴

adults	5.65 – 17.65 g/l
children (> 6 years)	6.50 – 16.00 g/l
children (1 month – 6 years)	2.00 – 12.40 g/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL III (Cat. No 4-291) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended.

0.9% NaCl should be used as calibrator 0.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using an automatic analyser Hitachi or BS-480. Results may vary if a different instrument is used.

- Measurement range:** 0.125 g/l to 35 g/l.

- Specificity / Interferences**

Hemoglobin up to 0.32 g/dl, bilirubin up to 22 mg/dl, triglycerides up to 940 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	6.83	0.01	1.43
level 2	14.03	0.16	1.12
level 3	19.86	0.36	1.80

Reproducibility (day to day) n = 66	Mean [g/l]	SD [g/l]	CV [%]
level 1	6.57	0.29	4.46
level 2	13.66	0.43	3.15
level 3	18.87	0.56	2.95

- Method comparison**

A comparison between IgG values determined at BS-480 (y) and at Advia 1650 (x) using 50 samples gave following results:

$$y = 0.9635 x + 0.6834 \text{ g/l;}$$

$$R = 0.983 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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• Calibration	
Calibration	Judgment Criteria
Rule <input type="text" value="Logit-Log 5P"/>	Sensitivity <input type="text"/>
Replicate <input type="text" value="1"/>	Factor Diff. <input type="text"/>
K <input type="text"/>	SD <input type="text"/>
• QC	Blank Abs. <input type="text"/>
Rules	Error Limit <input type="text"/>
Westgard Multi-rule	Corr. Coeff. <input type="text"/>
<input type="checkbox"/> v <input type="checkbox"/> R-4S	
<input type="checkbox"/> v <input type="checkbox"/> 4-1S	
<input type="checkbox"/> v <input type="checkbox"/> 10-X	
	Cum. Sum Check <input type="text"/>
	Interval <input type="text"/>
	• 1.0 - 2.7
	• 1.0 - 3.0
	• 0.5 - 5.1

• BS-480

Chem <input type="text" value="IgG"/>	No. <input type="text" value="031"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="Immunoglobulin G"/>	Print name <input type="text" value="IgG"/>	
Reaction Type <input type="text" value="Endpoint"/>	Reaction Direction <input type="text" value="Increase"/>	
Pri Wave <input type="text" value="570"/>	Sec Wave <input type="text" value="800"/>	
Unit <input type="text" value="g/L"/>	Decimal <input type="text" value="0.01"/>	
Blank Time <input type="text" value="48"/> <input type="text" value="49"/>	Reaction Time <input type="text" value="68"/> <input type="text" value="69"/>	
Standard <input type="text" value="3"/> <input type="text" value="3"/> <input type="text" value="3"/>	Reagent Vol <input type="text" value="250"/> <input type="text" value="95"/> <input type="text" value="95"/>	Diluent <input type="text" value="3"/> <input type="text" value="20"/> <input type="text" value="180"/>
Decreased <input type="text" value="3"/> <input type="text" value="3"/> <input type="text" value="3"/>	R1 <input type="text" value="250"/> <input type="text" value="95"/> <input type="text" value="95"/>	R2 <input type="text" value="250"/> <input type="text" value="95"/> <input type="text" value="95"/>
Increased <input type="text" value="3"/> <input type="text" value="3"/> <input type="text" value="3"/>	R3 <input type="text" value="250"/> <input type="text" value="95"/> <input type="text" value="95"/>	R4 <input type="text" value="250"/> <input type="text" value="95"/> <input type="text" value="95"/>
<input type="checkbox"/> Sample Blank <input type="checkbox"/> Auto Retun		
Linearity Range <input type="text" value="0.39"/> <input type="text" value="21"/>	Linearity Limit <input type="text"/>	
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>	
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-33000"/> <input type="text" value="33000"/>	
R1 Blank Abs <input type="text" value="-33000"/> <input type="text" value="33000"/>	Uncapping Time <input type="text" value="56"/> Dav(s)	
Blank Response <input type="text" value="-33000"/> <input type="text" value="33000"/>	Reagent Alarm Limit <input type="text"/>	
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension	
<input type="checkbox"/> Prozone Check <input type="checkbox"/> Rate Check	• Antigen Addition	
Q1 <input type="text" value="0"/> Q2 <input type="text" value="0"/>	Q3 <input type="text" value="0"/> Q4 <input type="text" value="0"/>	
PC <input type="text" value="0"/> ABS <input type="text" value="0"/>		

Calibration Settings	Auto Calibration
Math Model <input type="text" value="Logi-Log 4P"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="1"/>	<input type="checkbox"/> Cal Time
Acceptance Limits	
Cal Time <input type="text" value="672"/> Hour	
Slope Diff <input type="text"/>	SD <input type="text"/>
Sensitivity <input type="text"/>	Repeatability <input type="text"/>
Deter Coeff <input type="text"/>	

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