A-400 TRANSFERRIN

Nr kat. **7-410** (PL)

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Cat. No **7-410** (EN)

INTENDED USE

Diagnostic kit for determination of transferrin concentration intended to use in automatic analyser BS-400.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Transferrin (siderophilin) is a glycoprotein synthesised in the liver. Transferrin is the major plasma transport protein for iron. Its concentration correlates with the total iron-binding capacity (TIBC).

Evaluation plasma/serum transferrin levels is useful for the differential diagnosis of anaemia and for monitoring its treatment. Causes of decreased synthesis and low plasma levels may be chronic liver disease and malnutrition. High levels of transferrin occur in pregnancy and during estrogen administration.

METHOD PRINCIPLE

Transferrin, present in a sample, form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum is proportional to transferrin concentration in the sample.

REAGENTS

Package

1-Reagent 1 x 40 ml
 2-Reagent 1 x 12 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: imidazole buffer, polymer, inorganic salt, buffer, preservative.

2-Reagent: anti-human transferrin antiserum, inorganic salt, buffer, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN ⁶

Serum or plasma. Nonhemolyzed, nonlipaemic, fresh serum is recommended.

Serum can be stored up to 3 days at 4°C or up to 6 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

0.9% NaCl is recommended as a reagent blank.

REFERENCE VALUES ⁶

newborns	1.30 – 2.75 g/l
children (3 months – 16 years)	2.03 – 3.60 g/l
adults F	2.50 – 3.80 g/l
M	2.15 – 3.65 g/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended. 0.9% NaCl should be used as a calibrator 0.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Hitachi. Results may vary if a different instrument is used.

- LOQ:** 0.076 g/l

- Linearity:** up to 8 g/l

- Specificity / Interferences**

Haemoglobin up to 0.32 g/dl, bilirubin up to 29.5 mg/dl, triglycerides up to 312 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	1.44	0.02	1.12
level 2	2.52	0.04	1.52
level 3	3.47	0.06	1.83
Reproducibility (day to day) n = 66	Mean [g/l]	SD [g/l]	CV [%]
level 1	1.49	0.11	4.56
level 2	2.51	0.12	4.06
level 3	3.51	0.12	1.64

- Method comparison**

A comparison between CORMAY reagent (y) and another commercially available assay (x) using at least 30 samples gave following results:

$y = 0.91x - 0.0007$ g/l;

$R = 0.925$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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