A-400 COMPLEMENT C3

Nr kat. 7-411

(PL)

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Cat. No **7-411** (EN)

INTENDED USE

Diagnostic kit for determination of complement C3 concentration intended to use in automatic analyser BS-400. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Complement is a group of 20 immunologically distinct proteins present in blood and tissues. They are able to interact with antigen-antibody complexes, with each other and with cell membranes, in a complex way intended at destroying viruses and bacteria. They are synthesised in liver and are present in serum as functionally inactive molecules. They are activated by antigen-antibody complexes. C3 complement is a α -glycoprotein of 2 subunits. It is an acute phase reactant whose levels are increased during the acute phase. Low levels are found in immune complex diseases and in inherited deficiency which results in recurrent infections.

METHOD PRINCIPLE

The complement C3, present in a sample, form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to complement C3 concentration in the sample.

REAGENTS

Package

1-Reagent 1 x 40 ml
 2-Reagent 1 x 9 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: Imidazole buffer, polymer, inorganic salt, buffer, preservative.

2-Reagent: anti-human complement C3 antiserum, inorganic salt, buffer, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.

- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum or plasma.

Serum should be separated from red blood cells as soon as possible after blood collection. If the test cannot be done immediately, the sample should be stored at -70°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
 For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ²

adults	0.9 – 1.8 g/l
--------	---------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended. As a 0 calibrator 0.9% NaCl should be used.

The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyzer Hitachi. Results may vary if a different instrument is used.

- LOQ:** 0.015 g/l
- Linearity:** up to 4.31 g/l
- Specificity / Interferences**
 Haemoglobin up to 0.32 g/dl, bilirubin up to 22 mg/dl, triglycerides up to 312 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.759	0.005	0.66
level 2	1.615	0.013	0.80
level 3	2.367	0.021	0.87

Reproducibility (day to day) n = 59	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.778	0.019	2.42
level 2	1.634	0.026	1.59
level 3	2.428	0.045	1.87

- Method comparison**

A comparison between CORMAY reagent (y) and commercially available assay (x) using at least 30 samples gave following results:

$$y = 1.0577x - 0.0004 \text{ g/l;}$$

$$R = 0.930 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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A-400 COMPLEMENT C3

PROGRAM NA ANALIZATOR / APPLICATION / АДАПТАЦИЯ:

• Basic																									
Test information	<table border="1"> <tr> <td>No.</td> <td>35</td> </tr> <tr> <td>Test</td> <td>C3</td> </tr> <tr> <td>Full Name</td> <td>Complement C3</td> </tr> <tr> <td>Std. No.</td> <td>35</td> </tr> </table>	No.	35	Test	C3	Full Name	Complement C3	Std. No.	35																
No.	35																								
Test	C3																								
Full Name	Complement C3																								
Std. No.	35																								
Reagent Volume	<table border="1"> <tr> <td>R1</td> <td>200</td> </tr> <tr> <td>R2</td> <td>40</td> </tr> <tr> <td>R3</td> <td></td> </tr> <tr> <td>R4</td> <td></td> </tr> </table>	R1	200	R2	40	R3		R4																	
R1	200																								
R2	40																								
R3																									
R4																									
Sample Volume	<table border="1"> <tr> <td>Standard</td> <td>3.5</td> <td>15</td> <td>10</td> </tr> <tr> <td>Increased</td> <td>7</td> <td>15</td> <td>10</td> </tr> <tr> <td>Decreased</td> <td></td> <td></td> <td></td> </tr> </table>	Standard	3.5	15	10	Increased	7	15	10	Decreased															
Standard	3.5	15	10																						
Increased	7	15	10																						
Decreased																									
Reaction Parameters	<table border="1"> <tr> <td>Reac. Type</td> <td>Endpoint</td> <td>Direction</td> <td>Increase</td> </tr> <tr> <td>Pri. Wave</td> <td>340</td> <td>Rtg. Blank</td> <td>41 42</td> </tr> <tr> <td>Sec. Wave</td> <td>700</td> <td>Reac. Time</td> <td>76 77</td> </tr> </table>	Reac. Type	Endpoint	Direction	Increase	Pri. Wave	340	Rtg. Blank	41 42	Sec. Wave	700	Reac. Time	76 77												
Reac. Type	Endpoint	Direction	Increase																						
Pri. Wave	340	Rtg. Blank	41 42																						
Sec. Wave	700	Reac. Time	76 77																						
Judgment Criteria	<table border="1"> <tr> <td>Absorbance</td> <td>0 0</td> <td>Lin. Range</td> <td></td> </tr> <tr> <td>Incre. Test</td> <td>0</td> <td>Lin. Limit</td> <td></td> </tr> <tr> <td>Decre. Test</td> <td>0</td> <td>Subs. Limit</td> <td></td> </tr> </table>	Absorbance	0 0	Lin. Range		Incre. Test	0	Lin. Limit		Decre. Test	0	Subs. Limit													
Absorbance	0 0	Lin. Range																							
Incre. Test	0	Lin. Limit																							
Decre. Test	0	Subs. Limit																							
Result Setup	<table border="1"> <tr> <td>Decimal</td> <td>0.001</td> <td>Slope</td> <td>1</td> </tr> <tr> <td>Unit</td> <td>g/l</td> <td>Inter</td> <td>0</td> </tr> </table>	Decimal	0.001	Slope	1	Unit	g/l	Inter	0																
Decimal	0.001	Slope	1																						
Unit	g/l	Inter	0																						
QC	<table border="1"> <tr> <td>Q1</td> <td>0</td> <td>Prozone</td> <td>0</td> <td>Rate</td> <td>0</td> <td>Antigen</td> <td>0</td> </tr> <tr> <td>PC</td> <td>0</td> <td>Q2</td> <td>0</td> <td>Q3</td> <td>0</td> <td>Q4</td> <td>0</td> </tr> <tr> <td>ABS</td> <td>0</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> </table>	Q1	0	Prozone	0	Rate	0	Antigen	0	PC	0	Q2	0	Q3	0	Q4	0	ABS	0						
Q1	0	Prozone	0	Rate	0	Antigen	0																		
PC	0	Q2	0	Q3	0	Q4	0																		
ABS	0																								

• Calibration																					
Calibration	<table border="1"> <tr> <td>Rule</td> <td>Logit-Log 5P</td> </tr> <tr> <td>Replicate</td> <td>1</td> </tr> <tr> <td>K</td> <td></td> </tr> </table>	Rule	Logit-Log 5P	Replicate	1	K															
Rule	Logit-Log 5P																				
Replicate	1																				
K																					
Judgment Criteria	<table border="1"> <tr> <td>Sensitivity</td> <td></td> <td>Blank Abs.</td> <td></td> </tr> <tr> <td>Factor Diff.</td> <td></td> <td>Error Limit</td> <td></td> </tr> <tr> <td>SD</td> <td></td> <td>Corr. Coeff.</td> <td></td> </tr> </table>	Sensitivity		Blank Abs.		Factor Diff.		Error Limit		SD		Corr. Coeff.									
Sensitivity		Blank Abs.																			
Factor Diff.		Error Limit																			
SD		Corr. Coeff.																			
• QC	<table border="1"> <tr> <td>Rules</td> <td></td> <td>Auto QC</td> <td></td> </tr> <tr> <td>Westgard Multi-rule</td> <td></td> <td>Interval</td> <td></td> </tr> <tr> <td>v 1-2S</td> <td>v R-4S</td> <td></td> <td></td> </tr> <tr> <td>v 1-3S</td> <td>v 4-1S</td> <td></td> <td></td> </tr> <tr> <td>v 2-2S</td> <td>v 10-X</td> <td></td> <td></td> </tr> </table>	Rules		Auto QC		Westgard Multi-rule		Interval		v 1-2S	v R-4S			v 1-3S	v 4-1S			v 2-2S	v 10-X		
Rules		Auto QC																			
Westgard Multi-rule		Interval																			
v 1-2S	v R-4S																				
v 1-3S	v 4-1S																				
v 2-2S	v 10-X																				
Cum. Sum Check	<ul style="list-style-type: none"> 1.0 - 2.7 1.0 - 3.0 0.5 - 5.1 																				

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