## A-400 ASAT

Nr kat. **7-414** (PL)

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Cat. No **7-414** (EN)

### INTENDED USE

Diagnostic kit for determination of aspartate aminotransferase activity intended to use in automatic analyzers BS-400 and BS-480.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Aspartate aminotransferase (ASAT, AST, GOT) is an enzyme participated in amino acids metabolism. ASAT is found in all tissues but particularly high level of ASAT is observed in heart muscle, skeletal muscle, liver and kidney. This is why elevated ASAT serum level is marker of myocardial infarction and kidney, liver or skeletal muscle injury.

### METHOD PRINCIPLE

Optimized, modified method according to International Federation of Clinical Chemistry (IFCC), without pyridoxal phosphate.

The rate of absorbance changing at  $\lambda=340$  nm is directly proportional to aspartate aminotransferase activity.

### REAGENTS

#### Package

1-Reagent 6 x 30,5 ml  
 2-Reagent 3 x 16,5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 12 weeks on board the analyser at 2-10°C.

### Concentrations in the test

#### 1-Reagent

Tris buffer, pH 7,7 ≤ 120 mmol/l  
 L-aspartate ≤ 360 mmol/l  
 MDH ≤ 1.4 U/ml  
 LDH ≤ 2.3 U/ml

stabilizer, preservative

#### 2-Reagent

2-oxoglutarate ≤ 74 mmol/l  
 NADH ≤ 1.7 mmol/l  
 buffer  
 preservatives

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Warning

H315 Causes skin irritation.  
 H319 Causes serious eye irritation.  
 P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 IF ON SKIN: Wash with plenty of soap and water.  
 P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### SPECIMEN

Serum, heparinized or EDTA plasma free from hemolysis. Do not use heparine ammonium salt. Hemolysis should be avoided, since ASAT activity in erythrocytes is 10 times higher than in normal serum. Do not freeze the samples. ASAT activity remains stable in specimen up to 1 day at 15-25°C or up to 4 days at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.  
 For reagent blank deionized water is recommended.

#### Actions required:

When performing assays at analyser **BS-400**, there is a probability of **cross-contamination** affecting the tests results: : ASAT – LDH, LACTATE - ASAT. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_002\_BS-400\_CARRYOVER.

### REFERENCE VALUES <sup>6</sup>

| serum / plasma | 37°C         |                    |
|----------------|--------------|--------------------|
| female         | up to 31 U/l | up to 0.518 µkat/l |
| male           | up to 37 U/l | up to 0.618 µkat/l |

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

▪ **LOQ (Limit of Quantitation):** 7.0 U/l (0.12 µkat/l)

▪ **Linearity:** up to 1000 U/l (16.7 µkat/l)

#### Specificity / Interferences

Haemoglobin up to 0.63 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

#### Precision

| Repeatability<br>(run to run)   |         | Mean<br>[U/l] | SD<br>[U/l] | CV<br>[%] |
|---------------------------------|---------|---------------|-------------|-----------|
| <b>BS-400</b><br>(n = 20)       | level 1 | 42.0          | 0.85        | 2.0       |
|                                 | level 2 | 206.0         | 1.0         | 0.5       |
| <b>BS-480</b><br>(n = 20)       | level 1 | 44.0          | 0.69        | 1.6       |
|                                 | level 2 | 206.0         | 1.2         | 0.6       |
| Reproducibility<br>(day to day) |         | Mean<br>[U/l] | SD<br>[U/l] | CV<br>[%] |
| <b>BS-400</b><br>(n = 80)       | level 1 | 42.7          | 0.96        | 2.2       |
|                                 | level 2 | 211.0         | 3.1         | 1.5       |
| <b>BS-480</b><br>(n = 80)       | level 1 | 44.7          | 0.92        | 2.1       |
|                                 | level 2 | 211.0         | 1.7         | 0.8       |

#### Method comparison

A comparison between ASAT values determined at **BS-400** (y) and at **COBAS 6000** (x) using 155 serum samples gave following results:

$$y = 1.0405x - 0.048 \text{ U/l};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

A comparison between ASAT values determined at **BS-480** (y) and at **BECKMAN COULTER AU680** (x) using 97 serum samples gave following results:

$$y = 0.9175x + 1.6552 \text{ U/l};$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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• **Calibration**

|                     |                    |                          |             |                |  |
|---------------------|--------------------|--------------------------|-------------|----------------|--|
| <b>Calibration</b>  |                    | <b>Judgment Criteria</b> |             |                |  |
| Rule                | Multi-point Linear | Sensitivity              |             | Blank Abs.     |  |
| Replicate           | 2                  | Factor Diff.             |             | Error Limit    |  |
| K                   |                    | SD                       |             | Corr. Coeff.   |  |
| <b>• QC</b>         |                    |                          |             | <b>Auto QC</b> |  |
| <b>Rules</b>        |                    | Cum. Sum Check           |             | Interval       |  |
| Westgard Multi-rule |                    |                          | 1.0 - 2.7   |                |  |
| v                   | 1-2S               | v                        | 4-S         |                |  |
| v                   | 1-3S               | v                        | 4-1S        |                |  |
| v                   | 2-2S               | v                        | 10-X        |                |  |
|                     |                    |                          | • 1.0 - 3.0 |                |  |
|                     |                    |                          | 0.5 - 5.1   |                |  |

• **BS-480**

|  |                                     |  |                |                    |            |
|--|-------------------------------------|--|----------------|--------------------|------------|
| Chem                                   | ASAT II GEN                         | No.  | 018            | Sample Type        | SERUM      |
| Chemistry                              | ASAT                                | Print name                                       | ASAT           |                    |            |
| Reaction Type                          | Kinetic                             | Reaction Direction                               | Decrease       |                    |            |
| Pri Wave                               | 340                                 | Sec Wave   | 450            |                    |            |
| Unit                                   | U/L                                 | Decimal  | 0.1            |                    |            |
| Blank Time                             | 38 48                               | Reaction Time                                    | 58 78          |                    |            |
| Standard                               | Sample Vol 10 μL                    | Aspirated 20 μL                                  | Diluent 180 μL | Reagent Vol 200 μL | Diluent μL |
| Decreased                              | 10 μL                               | 20 μL  | 180 μL         | R2 50 μL           | μL         |
| Increased                              | μL                                  | μL   | μL             | R3 μL              | μL         |
|  | μL                                  | μL   | μL             | R4 μL              | μL         |
|  | Sample Blank                        | <input checked="" type="checkbox"/> Auto Rerun   |                |                    |            |
| Linearity Range (Standard)             | 7 1000                              | Linearity Limit                                  | 0.2            |                    |            |
| Linearity Range (Decreased)            |                                     | Substrate Depletion                              | 5600           |                    |            |
| Linearity Range (Increased)            |                                     | Mixed Blank Abs                                  | -33000 33000   |                    |            |
| R1 Blank Abs                           | -33000 33000                        | Uncapping Time                                   | 84 Day(s)      |                    |            |
| Blank Response                         | -33000 33000                        | Reagent Alarm Limit                              |                |                    |            |
| Twin Chemistry                         |                                     | <input type="checkbox"/> Enzyme Linear Extension |                |                    |            |
| <input type="checkbox"/> Prozone Check | <input type="checkbox"/> Rate Check | • Antigen Addition                               |                |                    |            |
| O1 0                                   | O2 0                                | O3 0   | O4 0           |                    |            |
| PC 0                                   | ABS 0                               |  |                |                    |            |

|                             |                    |   |  |
|-----------------------------|--------------------|---|--|
| <b>Calibration Settings</b> |                    | <b>Auto Calibration</b>                 |  |
| Math Model                  | Multi-point Linear | <input type="checkbox"/> Bottle Changed |  |
| Factor                      |                    | <input type="checkbox"/> Lot Changed    |  |
| Replicates                  | 2                  | <input type="checkbox"/> Cal Time       |  |
| <b>Acceptance Limits</b>    |                    |   |  |
| Cal Time                    | 2016 Hour          | SD                                      |  |
| Slope Diff                  |                    | Repeatability                           |  |
| Sensitivity                 |                    |   |  |
| Deter Coeff                 |                    |   |  |

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