A-400 ALAT

Nr kat. 7-416

(PL)

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Cat. No. **7-416** (EN)

INTENDED USE


Diagnostic kit for determination of alanine aminotransferase activity used in automatic analysers BS-400 and BS-480. The reagents must be used only for in vitro diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Alanine aminotransferase (ALAT, ALT, GPT) is an enzyme participated in amino acids metabolism. ALAT is present in all tissues but the highest level is found in liver and kidney cells. Damage of hepatocytes or kidney cells causes significant release of ALAT into the circulation. Measurement of ALT activity in serum is valuable in the diagnosis of liver diseases: jaundice, mononucleosis or hepatic cirrhosis.

METHOD PRINCIPLE

Optimized, modified method according to International Federation of Clinical Chemistry (IFCC), without pyridoxal phosphate.

The rate of absorbance changing at $\lambda=340$ nm is directly proportional to alanine aminotransferase activity.

REAGENTS

Package

1-Reagent 6 x 30.5 ml
 2-Reagent 3 x 16.5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 12 weeks on board the analyser at 2-10°C.

Concentrations in the test

1-Reagent

Tris (pH 7,4) ≤ 150 mmol/l
 L-alanine ≤ 750 mmol/l
 LDH ≤ 4 U/ml

stabilizers, preservative

2-Reagent

2-oxoglutarate ≤ 74 mmol/l
 NADH ≤ 1.7 mmol/l

buffer

preservatives

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum, heparinized or EDTA plasma free from hemolysis. Hemolysis should be avoided, since ALAT activity in erythrocytes is 3 to 5 times higher than in normal serum. Do not freeze the samples. ALAT activity remains stable in specimen up to 3 days at 15-25°C or up to 7 days at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank deionized water is recommended.

Actions required:

When performing assays at analyser **BS-400**, there is a probability of **cross-contamination** affecting the tests results: ALAT – LDH, LACTATE - ALAT. To avoid this effect follow the recommendations contained in the instruction 51_03_24_002_BS-400_CARRYOVER

REFERENCE VALUES ⁶

serum / plasma	37°C	
women	up to 31 U/l	up to 0.517 μ kat/l
men	up to 41 U/l	up to 0.683 μ kat/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0. The calibration curve should be prepared every 12 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

▪ **LOQ (Limit of Quantitation):** 6 U/l (0.1 μ kat/l)

▪ **Linearity:** up to 1000 U/l (16.7 μ kat/l)

▪ **Specificity / Interferences**

Haemoglobin up to 0.31 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [U/l]	SD [U/l]	CV [%]
BS-400 (n = 20)	level 1	31.4	0.65	2.1
	level 2	103	0.65	0.63
BS-480 (n = 20)	level 1	31.5	0.77	2.4
	level 2	102	0.80	0.78
Reproducibility (day to day)		Mean [U/l]	SD [U/l]	CV [%]
BS-400 (n = 80)	level 1	32.8	1.1	3.4
	level 2	106	1.4	1.4
BS-480 (n = 80)	level 1	33.3	1.0	2.9
	level 2	105	1.2	1.1

Method comparison

A comparison between ALAT values determined at **BS-400** (y) and at **COBAS 6000** (x) using 128 serum samples gave following results:

$$y = 1.0266x + 0.9243 \text{ U/l}$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

A comparison between ALAT values determined at **BS-480** (y) and at **BECKMAN COULTER AU680** (x) using 97 serum samples gave following results:

$$y = 0.9772x + 1.6984 \text{ U/l}$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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• Calibration

Calibration		Judgment Criteria			
Rule	Multi-point Linear	Sensitivity		Blank Abs.	
Replicate	2	Factor Diff.		Error Limit	
K		SD		Corr. Coeff.	
QC				Auto QC	
Rules					
Westgard Multi-rule		Cum. Sum Check		Interval	
v	1-2S	v	R-4S		1.0 - 2.7
v	1-3S	v	4-1S		• 1.0 - 3.0
v	2-2S	v	10-X		0.5 - 5.1

• BS-480

Chem	ALAT II GEN	No.	017	Sample Type	SERUM
Chemistry	ALAT	Print name	ALAT		
Reaction Type	Kinetic	Reaction Direction	Decrease		
Pri Wave	340	Sec Wave	450		
Unit	U/L	Decimal	0.1		
Blank Time	38 48	Reaction Time	58 78		
Standard	Sample Vol 10 μL Aspirated 10 μL Diluent 180 μL	Reagent Vol	R1 200 μL R2 50 μL R3 μL R4 μL	Diluent	μL μL μL μL
Decreased	10 μL 20 μL 180 μL				
Increased	μL μL μL				
	Sample Blank Auto Rerun				
Linearity Range (Standard)	6 1000	Linearity Limit	0.2		
Linearity Range (Decreased)		Substrate Depletion	5900		
Linearity Range (Increased)		Mixed Blank Abs	-33000 33000		
R1 Blank Abs	-33000 33000	Uncapping Time	84 Dav(s)		
Blank Response	-33000 33000	Reagent Alarm Limit			
Twin Chemistry		Enzyme Linear Extension			
<input type="checkbox"/> Prozone Check	<input type="checkbox"/> Rate Check	Antigen Addition			
Q1 0	Q2 0	Q3 0	Q4 0		
PC 0	ABS 0				

Calibration Settings		Auto Calibration	
Math Model	Multi-point Linear	<input type="checkbox"/> Bottle Changed	
Factor		Replicates	2
		<input type="checkbox"/> Lot Changed	
		<input type="checkbox"/> Cal Time	
Acceptance Limits			
Cal Time	2016 Hour	SD	
Slope Diff		Repeatability	
Sensitivity			
Deter Coeff			

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