A-400 PROCALCITONIN

Nr kat. 7-417

(PL)

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LITERATURA

1. Müller B, et al., Ubiquitous expression of the calcitonin-i gene in multiple tissues in response to sepsis. *J Clin Endocrinol Metab.* 2001, 86(1): 396-404.
2. Meisner M. Procalcitonin (PCT) – A new, innovative infection parameter. *Biochemical and clinical aspects.* Thieme; Stuttgart, New York, 2000; ISBN 3-13-105503-0.
3. Christ-Crain M, et al., Procalcitonin in bacterial infections – hype, hope or more or less? *Swiss Med Wkly.* 2005; 135 (31-32): 451-60.
4. Bone RC, Balk RA, Cerra FB, Dellinger RP, Fein AM, Knaus WA, Schein RM, Sibbald WJ. Definitions for sepsis and organ failure and guidelines for the use of innovative therapies in sepsis. The ACCP/SCCM Consensus Conference Committee. American College of Chest Physicians/Society of Critical Care Medicine. *Chest* 1992; 101: 1644-55.
5. Christ-Crain M, Jaccard-Stolz D, Bingisser R, Gencay MM, Huber PR, Tamm M, Müller B. Effect of procalcitonin-guided treatment on antibiotic use and outcome in lower respiratory tract infections: cluster-randomised, single-blinded intervention trial. *Lancet.* 2004;363: 600-7.
6. Christ-Crain M, Stolz D, Bingisser R, Müller C, Miedinger D, Huber PR, Zimmerli W, Harbarth S, Tamm M, Müller B. Procalcitonin Guidance of Antibiotic Therapy in Community-acquired Pneumonia: A Randomized Trial. *Am J Respir Crit Care Med.* 2006; 174: 84-93.
7. Stolz D, Christ-Crain M, Bingisser R, Leuppi J, Miedinger D, Müller C, Huber P, Müller B, Tamm M. Antibiotic treatment of exacerbations of COPD: a randomized, controlled trial comparing procalcitonin-guidance with standard therapy. *Chest.* 2007;131: 9-19.
8. Schuetz P, et al., Long-term stability of procalcitonin in frozen samples and comparison of Kryptor and VIDAS automated immunoassays. *Clin Biochem.* 2010; 43(3):341-4.
9. Müller B, et al., Calcitonin precursors are reliable markers of sepsis in a medical intensive care unit. *Crit Care Med.* 2000; 28(4): 977-983.
10. Harbarth S, et al., Diagnostic value of procalcitonin, interleukin-6, and interleukin-8 in critically ill patients admitted with suspected sepsis. *Am J Respir Crit Care Med.* 2001; 164: 396-402.
11. Dymicka-Piekarska V., Wasiluk A. Prokalcytonina (PCT), współczesny wskaźnik infekcji i stanów zapalnych. *Postepy Hig Med Dosw (online).* 2015; 69: 723-728.
12. Herzum I, Renz H.: Inflammatory markers in SIRS, sepsis and septic shock. *Curr. Med. Chem.,* 2008; 15: 581-587.
13. Boeken U, Feindt P, Micek M, Petzold T, Schulte HD, Gams E., Procalcitonin (PCT) in cardiac surgery: diagnostic value in systemic inflammatory response syndrome (SIRS), sepsis and after heart transplantation (HTX). *Cardiovasc Surg.* 2000 Dec; 8(7): 550-4.
14. Vijayan et al., Procalcitonin: a promising diagnostic marker for sepsis and antibiotic Therapy. *Journal of Intensive Care* (2017) 5:51.

Data wydania: 12. 2019.

A-400 PROCALCITONIN

Cat.No 7-417

(EN)

INTENDED USE

Diagnostic kit for determination of procalcitonin concentration intended to use in automatic analyzers BS-400 and BS-480.
The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Procalcitonin (PCT) is a 116 amino acid polypeptide, the precursor of calcitonin (CT). Under physiological conditions, calcitonin is produced exclusively in the C-cells of the thyroid gland through specific intracellular proteolysis of PCT. In a healthy individual, the blood PCT level is very low (< 0.05 ng/ml). In response to bacterial infection/ sepsis and inflammatory response, bacterial toxins and pro-inflammatory mediators stimulate release of PCT directly to bloodstream by many body tissues (PCT is not transfer to calcitonin). PCT level increases about 3 hours after bacterial infection, reaching maximum values after 6-12 hours. PCT half life is from 25 to 30 hours, therefore determination of procalcitonin concentration in order to monitor the progress of infection can be done once per day. [1-8]

Elevated PCT level may not always be caused by systemic bacterial infection [2]. If there is a disagreement between the laboratory findings and the clinical signs, additional tests should be performed.
Measurements of PCT in conjunction with other laboratory findings and clinical assessment aids in the risk assessment of critically ill patients of their first day of ICU admission for progression to severe sepsis and septic shock.

PCT determination allows to differentiate the systemic inflammatory response syndrome (SIRS) from sepsis, differentiate bacterial infection from a viral one (PCT concentration is low during a viral infection), early detection of sepsis and monitoring the progress of bacterial infection, as well as fungal and parasitic infections. [11-13]

METHOD PRINCIPLE

A-400 PROCALCITONIN is a latex immunoturbidimetric assay. PCT proteins in the sample bind to the specific anti-PCT antibody, which is coated on latex particles, resulting in agglutination (forming PCT-antibody complexes). The degree of the turbidity caused by agglutination is measured optically and is proportional to the amount of PCT in the sample. The analyser calculates the PCT concentration of a sample by interpolation of the obtained signal with regard to a 6-point calibration curve.

REAGENTS

Package
1-Reagent 1 x 21 ml
2-Reagent 1 x 9 ml

The reagents, if stored at 2-8°C, are stable up to expiry date printed on the package. The reagents stored on board of the analyser are stable for 17 days.

Concentrations in the test

Tris buffer 100 mM
anti-human PCT antibody coated latex particles 0.2 %
Sodium azide < 0.1 %

Warnings and notes

- Do not use after the expiration date.
- Protect from light and avoid contamination and evaporation!
- Do not freeze the reagent!
- In order to maximize the stability of the reagents, recap the bottles after use and store them in the fridge at 2-8°C.**
- Do not swap the caps between the bottles, because contamination can influence the results.
- Specimens containing human sourced materials should be handled as if potentially infectious
- The reagent contains < 0.1% sodium azide as a preservative. Avoid contact with skin and mucous membranes.
- Do not mix reagent of different lots.
- As with any latex turbidimetric immunoassays, PCT assay runs should be followed with appropriate and thorough wash steps.

A-400 PROCALCITONIN

51_03_01_057_05

- It is recommended not to use samples which appear to be lipemic or hemolyzed and, if possible, to collect new sample.
- PCT assay is not indicated to be used as a stand-alone diagnostic assay and should be used in conjunction with clinical signs and symptoms of infection and other diagnostic evidence.
- Severity of renal failure or insufficiency, may influence procalcitonin values and should be considered as potentially confounding clinical factors when interpreting PCT values.
- Increased PCT levels may be observed in severe illness such as polytrauma, burns, major surgery, prolonged or cardiogenic shock.

SPECIMEN [2, 3, 8]

Serum, EDTA or heparinized plasma (recommended: lithium salt)

It should be used only dedicated test tube and container for taking and preparing samples.

Do not use preservatives.

PCT values measured in arterial blood are about 4% higher than from venous blood.

For monitoring PCT concentration it is recommended to use the same type of sample.

PCT is relatively stable in both plasma and serum samples, and there are no special requirements for preanalytical sample handling.

Turbid samples must be clarified by centrifuging.

Stability of retained sample:

- when stored at 2-8°C - **3 days**
- when stored at room temperature - **24 hours**
- when stored at -20°C:
 - ✓ **2 months** (serum samples)
 - ✓ **6 months** (plasma samples).

4 cycles of freeze/thaw show a <10% loss of PCT in sample. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

The CORMAY PROCALCITONIN CALIBRATORS Level 0 is recommended as a reagent blank.

When it is possible, the tests for determination of **procalcitonin** concentration, should be performed in the **separate order**.

REFERENCE VALUES [2-4]

PCT	< 0.5 ng/ml
-----	-------------

According to the literature [2-4, 9-10, 14], diagnosis of systemic bacterial infection / sepsis is categorized as follows:

- 0.05 ng/ml < PCT < 0.5 ng/mL** – low risk of severe sepsis and/or septic shock, do not exclude local infection or a systemic infection in its initial stages (< 6 hours);
- 0.5 ng/mL ≤ PCT < 2 ng/mL** – should be interpreted considering the patient's history. When PCT < **2 ng/ml** it is recommended to retest PCT within 6-24 hours;
- PCT > 2 ng/mL** – high risk of severe sepsis and / or septic shock.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY PROCALCITONIN CONTROLS (Cat. No 5-166) with each batch of samples.

For the calibration of automatic analysers systems CORMAY PROCALCITONIN CALIBRATORS (Cat. No 5-136) is recommended. The CORMAY PROCALCITONIN CALIBRATORS Level 0 should be used as a calibrator 0.

The calibration curve should be prepared every 3 days for every new lot of reagents or as required e.g. quality control findings fall outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers: BS-400 and/or BS-480 and/or Hitachi 917. Results may vary if a different instrument is used.

Limit of Quantitation (LOQ)

Lowest concentration, at which inaccuracy of the test does not exceed 20%

0.35 ng/ml - **BS-400**

0.37 ng/ml - **BS-480**

Measurement range

0.35 ng/ml up to 52 ng/ml – **BS-400**

0.37 ng/ml up to 52 ng/ml – **BS-480**

For higher concentration of procalcitonin dilute the sample with calibrator Level 0 and repeat the assay. Multiply the result by the dilution factor.

Specificity / Interferences

The substances listed below did not show significant (<10%) interference with the PCT assay in the concentration given in the tables:

SUBSTANCES	CONCENTRATION
ascorbic acid	129 mg/dl
free bilirubin	30 mg/dl
bilirubin conjugated	30 mg/dl
hemoglobin	750 mg/dl
triglyceride	750 mg/dl
rheumatoid factor	75 IU/ml
albumin	4 g/dl
human calcitonin	60 ng/ml
human katalcalcin	10 ng/ml
human alpha-CGRP	10 µg/ml
human beta-CGRP	10 µg/ml
human anti-mouse IgG (HAMA)	350 ng/ml

DRUGS	CONCENTRATION
imipenem	0.5 mg/ml
cefotaxime	180 mg/dl
noradrenalin	4 µg/ml
dobutamine	22.4 µg/ml
unfractionated heparin	16000 U/l
furosemide	4 mg/dl
vankomycin	3 mg/ml
dopamine	26 g/dl

There exists a possibility that other substances not mentioned above may interfere with the determination of procalcitonin concentration.

Precision

Repeatability (run to run)		Mean [ng/ml]	SD [ng/ml]	CV [%]
BS-400 (n = 10)	level 1	1.00	0.04	3.99
	level 2	13.33	0.17	1.30
BS-480 (n = 10)	level 1	1.07	0.06	5.50
	level 2	14.00	0.29	2.04
Reproducibility (day to day)		Mean [ng/ml]	SD [ng/ml]	CV [%]
BS-400 (n = 10)	level 1	1.30	0.09	6.84
	level 2	16.09	0.88	5.47
BS-480 (n = 10)	level 1	1.04	0.06	5.83
	level 2	11.77	0.54	4.62

Method comparison

A comparison between procalcitonin values determined at **BS-400** (y) and at **Advia CENTAUR** (x), using 45 samples gave following results:

$y = 0.8211 x + 0.0565$ ng/ml;

$R = 0.992$ (R – correlation coefficient)

A comparison between procalcitonin values determined at **BS-480** (y) and at **Advia CENTAUR** (x), using 44 samples gave following results:

$y = 0.8446 x + 0.0779$ ng/ml;

$R = 0.983$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

A-400 PROCALCITONIN

LITERATURE

- Müller B, et al., Ubiquitous expression of the calcitonin-i gene in multiple tissues in response to sepsis. *J Clin Endocrinol Metab.* 2001, 86(1): 396-404.
- Meisner M. Procalcitonin (PCT) – A new, innovative infection parameter. *Biochemical and clinical aspects.* Thieme; Stuttgart, New York, 2000; ISBN 3-13-105503-0.
- Christ-Crain M, et al., Procalcitonin in bacterial infections – hype, hope or more or less? *Swiss Med Wkly.* 2005; 135 (31-32): 451-60.
- Bone RC, Balk RA, Cerra FB, Dellinger RP, Fein AM, Knaus WA, Schein RM, Sibbald WJ. Definitions for sepsis and organ failure and guidelines for the use of innovative therapies in sepsis. The ACCP/SCCM Consensus Conference Committee. American College of Chest Physicians/Society of Critical Care Medicine. *Chest* 1992; 101: 1644-55.
- Christ-Crain M, Jaccard-Stolz D, Bingisser R, Gencay MM, Huber PR, Tamm M, Muller B. Effect of procalcitonin-guided treatment on antibiotic use and outcome in lower respiratory tract infections: cluster-randomised, single-blinded intervention trial. *Lancet.* 2004;363: 600-7.
- Christ-Crain M, Stolz D, Bingisser R, Muller C, Miedinger D, Huber PR, Zimmerli W, Harbarth S, Tamm M, Muller B. Procalcitonin Guidance of Antibiotic Therapy in Community-acquired Pneumonia: A Randomized Trial. *Am J Respir Crit Care Med.* 2006; 174: 84-93.
- Stolz D, Christ-Crain M, Bingisser R, Leuppi J, Miedinger D, Muller C, Huber P, Muller B, Tamm M. Antibiotic treatment of exacerbations of COPD: a randomized, controlled trial comparing procalcitonin-guidance with standard therapy. *Chest.* 2007;131: 9-19.
- Schuetz P, et al., Long-term stability of procalcitonin in frozen samples and comparison of Kryptor and VIDAS automated immunoassays. *Clin Biochem.* 2010; 43(3):341-4.
- Müller B, et al., Calcitonin precursors are reliable markers of sepsis in a medical intensive care unit. *Crit Care Med.* 2000; 28(4): 977-983.
- Harbarth S, et al., Diagnostic value of procalcitonin, interleukin-6, and interleukin-8 in critically ill patients admitted with suspected sepsis. *Am J Respir Crit Care Med.* 2001; 164: 396-402.
- Dymicka-Piekarska V., Wasiluk A. Prokalcitonina (PCT), wspólny wskaźnik infekcji i stanów zapalnych. *Postepy Hig Med Dosw (online).* 2015; 69: 723-728.
- Herzum I, Renz H.: Inflammatory markers in SIRS, sepsis and septic shock. *Curr. Med. Chem.*, 2008; 15: 581-587.
- Boeken U, Feindt P, Micek M, Petzold T, Schulte HD, Gams E., Procalcitonin (PCT) in cardiac surgery: diagnostic value in systemic inflammatory response syndrome (SIRS), sepsis and after heart transplantation (HTX). *Cardiovasc Surg.* 2000 Dec; 8(7): 550-4.
- Vijayan et al., Procalcitonin: a promising diagnostic marker for sepsis and antibiotic therapy. *Journal of Intensive Care* (2017) 5:51.

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- Müller B, et al., Ubiquitous expression of the calcitonin-i gene in multiple tissues in response to sepsis. J Clin Endocrinol Metab. 2001, 86(1): 396-404.
- Meisner M. Procalcitonin (PCT) – A new, innovative infection parameter. Biochemical and clinical aspects. Thieme; Stuttgart, New York, 2000; ISBN 3-13-105503-0.
- Christ-Crain M, et al., Procalcitonin in bacterial infections – hype, hope or more or less? Swiss Med Wkly. 2005; 135 (31-32): 451-60.
- Bone RC, Balk RA, Cerra FB, Dellinger RP, Fein AM, Knaus WA, Schein RM, Sibbald WJ. Definitions for sepsis and organ failure and guidelines for the use of innovative therapies in sepsis. The ACCP/SCCM Consensus Conference Committee. American College of Chest Physicians/Society of Critical Care Medicine. Chest 1992; 101: 1644-55.
- Christ-Crain M, Jaccard-Stolz D, Bingisser R, Gencay MM, Huber PR, Tamm M, Müller B. Effect of procalcitonin-guided treatment on antibiotic use and outcome in lower respiratory tract infections: cluster-randomised, single-blinded intervention trial. Lancet. 2004;363: 600-7.
- Christ-Crain M, Stolz D, Bingisser R, Müller C, Miedinger D, Huber PR, Zimmerli W, Harbarth S, Tamm M, Müller B. Procalcitonin Guidance of Antibiotic Therapy in Community-acquired Pneumonia: A Randomized Trial. Am J Respir Crit Care Med. 2006; 174: 84-93.
- Stolz D, Christ-Crain M, Bingisser R, Leuppi J, Miedinger D, Müller C, Huber P, Müller B, Tamm M. Antibiotic treatment of exacerbations of COPD: a randomized, controlled trial comparing procalcitonin-guidance with standard therapy. Chest. 2007;131: 9-19.
- Schuetz P, et al., Long-term stability of procalcitonin in frozen samples and comparison of Kryptor and VIDAS automated immunoassays. Clin Biochem. 2010; 43(3):341-4.
- Müller B, et al., Calcitonin precursors are reliable markers of sepsis in a medical intensive care unit. Crit Care Med. 2000; 28(4): 977-983.
- Harbarth S, et al., Diagnostic value of procalcitonin, interleukin-6, and interleukin-8 in critically ill patients admitted with suspected sepsis. Am J Respir Crit Care Med. 2001; 164: 396-402.
- Dymicka-Piekarska V., Wasiluk A. Prokalcitonina (PCT), wspólniecznie wskaźnik infekcji i stanów zapalnych. Postępy Hig Med Dosw (online), 2015; 69: 723-728.
- Herzum I, Renz H.: Inflammatory markers in SIRS, sepsis and septic shock. Curr. Med. Chem., 2008; 15: 581-587.
- Boeken U, Feindt P, Micek M, Petzold T, Schulte HD, Gams E., Procalcitonin (PCT) in cardiac surgery: diagnostic value in systemic inflammatory response syndrome (SIRS), sepsis and after heart transplantation (HTX). Cardiovasc Surg. 2000 Dec; 8(7): 550-4.
- Vijayan et al., Procalcitonin: a promising diagnostic marker for sepsis and antibiotic Therapy. Journal of Intensive Care (2017) 5:51.

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• Calibration

Calibration		Judgment Criteria	
Rule	Logit-Log 5P	Sensitivity	
Replicate	2	Factor Diff.	
K		SD	
• QC		Auto QC	
Rules		Interval	
Westgard Multi-rule		Cum. Sum Check	
v	1-2S	v	R-4S
v	1-3S	v	4-1S
v	2-2S	v	10-X
		•	1.0 - 3.0
			0.5 - 5.1

• BS-480

Chem	PCT	No.	086	Sample Type	SERUM
Chemistry	Procalcitonin	Print name	PCT		
Reaction Type	Fixed-time	Reaction Direction	Increase		
Pri Wave	605	Sec Wave			
Unit	ng/ml	Decimal	0.01		
Blank Time		Reaction Time	55	82	
Sample Vol	35 μL	Aspirated		Diluent	
Standard		Reagent Vol	150 μL	Diluent	

Decreased	<input type="text"/> μL	<input type="text"/> μL	<input type="text"/> μL	R2	50	μL	<input type="text"/> μL
Increased	<input type="text"/> μL	<input type="text"/> μL	<input type="text"/> μL	R3	<input type="text"/>	μL	<input type="text"/> μL
	<input type="text"/> Sample Blank	<input type="text"/> Auto Return		R4	<input type="text"/>	μL	<input type="text"/> μL
Linearity Range (Standard)	0.37	52		Linearity Limit	<input type="text"/>		
Linearity Range (Decreased)	<input type="text"/>	<input type="text"/>		Substrate Depletion	<input type="text"/>		
Linearity Range (Increased)	<input type="text"/>	<input type="text"/>		Mixed Blank Abs	-33000	33000	
R1 Blank Abs	-33000	33000		Uncapping Time	17	Day(s)	
Blank Response	-33000	33000		Reagent Alarm Limit	<input type="text"/>		
Twin Chemistry	<input type="text"/>			<input type="checkbox"/> Enzyme Linear Extension			
<input type="checkbox"/> Prozone Check	<input type="checkbox"/> Rate Check			• Antigen Addition			
Q1	0	Q2	0	Q3	0	Q4	0
PC	0	ABS	0				

Calibration Settings		Auto Calibration	
Math Model	Logit-Log 5P	<input type="checkbox"/> Bottle Changed	
Factor	<input type="text"/>	Replicates	2
		<input type="checkbox"/> Lot Changed	
		<input type="checkbox"/> Cal Time	
Acceptance Limits			
Cal Time	72	Hour	
Slope Diff	<input type="text"/>	SD	<input type="text"/>
Sensitivity	<input type="text"/>	Repeatability	<input type="text"/>
Deter Coeff	<input type="text"/>		

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