A-400 TOTAL IgE

Nr kat. 7-431 (PL)

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(EN)

INTENDED USE

Diagnostic kit for determination of IgE level, used in automatic analysers BS-400 and BS-480.

The reagents must be used only for in vitro diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

IgE is an immunoglobulin with a molecular weight of approximately 190 kD normally present in the blood in trace amounts. Continual production of IgE antibodies in response to common naturally occurring allergens, however, often results in elevated serum levels and in the development of such clinically important Type I allergic reactions as asthma, hay fever, dermatitis and food allergies. Elevated IgE levels are also seen in parasitic (helminth) diseases, IgE myeloma, and in hepatitis. The measurement of IgE in human serum is thus considered to be useful in the diagnosis, treatment, assessment of disease progression, or postoperative prognosis for such conditions.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between IgE in a sample and anti-IgE antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of IgE in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

1-Reagent 1 x 33 ml
 2-Reagent 1 x 18 ml

The reagents when stored at 2-10°C are stable up to expiry date printed on the package. The reagents are stable for 11 weeks on board the analyser at 2-10°C.

Concentrations in the test

suspension of latex particles sensitized with 0.125
 (mouse) anti-IgE antibodies (pH 7.3) w/v%
 glycine buffer solution (pH 8.3)
 preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Care should be taken not to interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum or plasma (Na-EDTA, K-EDTA, Na-Heparin, Li-Heparin, citric acid).

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
 For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ³

serum, plasma	< 358 IU/ml
---------------	-------------

It is recommended for each laboratory to establish its own reference ranges for local population. Diagnosis should only be made after taking clinical symptoms and the results of other tests into consideration.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL II (Cat. No 4-290) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY IgE CALIBRATORS kit (Cat. No 4-280) is recommended.

As a 0 calibrator 0.9% NaCl should be used.

The calibration curve should be prepared every 2 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using automatic analysers BS-400 and/or BS-480. Results may vary if a different instrument is used.

- Sensitivity:**
25 IU/ml- BS-400, BS-480
- Linearity:**
Up to 1000 IU/ml- BS-400, BS-480

For higher concentration of IgE dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

- Specificity / Interferences**

Haemoglobin up to 0.5 g/dl, bilirubin up to 30 mg/dl, triglycerides up to 1500 mg/dl, RF up to 500 IU/ml do not interfere with the test.

- Precision**

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 (n = 10)	level 1	40,5	2,7	6,57
	level 2	427,4	7,7	1,80

Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 (n = 10)	level 1	68,45	2,63	3,84
	level 2	402,05	2,69	0,67

- Method comparison**

A comparison between IgE values determined at **BS-400** (y) and at **BS-800** (x) using 23 samples gave following results:

$$y = 1,063x - 3,7893$$

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A comparison between IgE values determined at **BS-480** (y) and at **BS-800** (x) using 21 samples gave following results:

$$y = 1,0507x - 1,8765$$

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WASTE MANAGEMENT

Please refer to local legal requirements.

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PROGRAM NA ANALIZATOR: / APPLICATION for: / АДАПТАЦИЯ для:

• **BS-400**

• Basic		Reagent Volume		Sample Volume	
Test information		R1	150	Standard	10 15 10
No.	39	R2	75	Increased	20 15 10
Test	IgE	R3		Decreased	5 15 10
Full Name	IgE	R4			
Std. No.	39				
Reaction Parameters		Direction		Result Setup	
Reac. Type	Kinetic	Rtg. Blank	Increase	Decimal	1 Slope 1
Pri. Wave	570	Reac. Time	0 0	Unit	IU/ml Inter 0
Sec. Wave	800		45 52		
Judgment Criteria		Judgment Criteria			
Absorbance	0 0	Lin. Range	25 1000	Prozone	<input type="checkbox"/>
Incr. Test	0	Lin. Limit	0.20	Rate	<input type="checkbox"/>
Decr. Test	0	Subs. Limit		Antigen	<input type="checkbox"/>
		Q1	0	Q2	0
		PC	0	Q3	0
				Q4	0
				ABS	0
• Calibration		Judgment Criteria			
Calibration		Sensitivity		Blank Abs.	
Rule	Logit-Log 5P	Factor Diff.		Error Limit	
Replicate	1	SD		Corr. Coeff.	
• QC		Auto QC			
Rules		Westgard Multi-rule		Cum. Sum Check	
v	1-2S	v	R-4S	1.0 - 2.7	
v	1-3S	v	4-1S	• 1.0 - 3.0	
v	2-2S	v	10-X	0.5 - 5.1	
• BS-480					
Chem	IgE	No.	039	Sample Type	SERUM
Chemistry	Total IgE			Print name	IgE
Reaction Type	Kinetic	Reaction Direction	Increase	Sec Wave	800
Pri Wave	570	Decimal	0.1	Reaction Time	53
Unit	IU/mL	Reaction Time	61	Reagent Vol	
Blank Time	0	Sample Vol		Diluent	
Standard	10 μL	Aspirated		Diluent	
Decreased	10 μL	20 μL		R1	150 μL
Increased	10 μL	180 μL		R2	75 μL
				R3	
				R4	
		Sample Blank	v	Auto Retun	
Linearity Range (Standard)	25	1000		Linearity Limit	0.6
Linearity Range (Decreased)				Substrate Depletion	
Linearity Range (Increased)				Mixed Blank Abs	-33000 33000
R1 Blank Abs	-33000	33000		Uncapping Time	77 Day(s)
Blank Response	-33000	33000		Reagent Alarm Limit	
Twin Chemistry				Enzyme Linear Extension	
Prozone Check	<input type="checkbox"/>	Rate Check	<input type="checkbox"/>	Antigen Addition	
Q1	0	Q2	0	Q3	0
PC	0	ABS	0	Q4	0
Calibration Settings		Auto Calibration			
Math Model	Logit-Log 5P	Bottle Changed			
Factor		Lot Changed			
Replicates	2	Cal Time			
Acceptance Limits					
Cal Time	336	Hour			
Slope Diff		SD			
Sensitivity		Repeatability			
Deter Coeff					

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