## A-400 ALBUMIN

Nr kat. 7-438

(PL)

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Cat. No 7-438 (EN)

### INTENDED USE

Diagnostic kit for determination of albumin concentration, intended to use in automatic analyzers: BS-400 and BS-480. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Albumin is the major serum protein, but is present also in other body fluids: cerebrospinal, pleural and peritoneal. Albumin regulates blood oncotic pressure and serves as amino acids reservoir. Beyond of these functions albumin is very important transport protein – binds and keeps dispersed bilirubin, hormones, vitamins, calcium, magnesium, fatty acids and medicines. Decreased albumin blood level is caused usually by liver or kidney disease, malabsorption or malnutrition.

### METHOD PRINCIPLE

Bromocresol green (BCG) forms with albumin, in succinate buffer (acid medium), a coloured complex. The absorbance of this complex is proportional to the albumin concentration in the sample. The colour intensity of the formed complex measured at 630 nm is proportional to albumin concentration in the sample.

### REAGENTS

#### Package

1-Reagent 2 x 33 ml

The reagent is stable up to the kit expiry date printed on the package when stored at 2-8°C. The reagents are stable for 12 weeks on board the analyser at 2-10°C.

### Concentrations in the test

succinate buffer ≤ 108 mmol/l  
bromocresol green (BCG) ≤ 0.35 mmol/l  
sodium hydroxide ≤ 60 mmol/l  
detergent  
preservative

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze the reagent!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-Reagent meets the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Warning

H319 Causes serious eye irritation.  
P280 Wear protective gloves, eye protection and face protection.  
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### SPECIMEN

Serum free from hemolysis.  
Serum should be separated from red blood cells as soon as possible after blood collection.  
Serum can be stored up to 3 days at 2-8°C or 6 months at -20°C.<sup>2</sup>  
Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent is ready to use.  
For reagent blank deionized water is recommended.

### REFERENCE VALUES<sup>7</sup>

serum		g/dl	g/l
children	3 – 5 years	3.9 – 5.0	39 - 50
	6 – 15 years	4.1 – 5.1	41 - 51
females	16 – 54 years	3.9 – 5.0	39 – 50
males	16 – 29 years	4.6 – 5.3	46 - 53
	30 – 54 years	4.4 – 5.1	44 - 51
females and males	55 – 79 years	4.2 – 5.0	42 – 50

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples. For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionized water should be used as calibrator 0.

The calibration curve should be prepared every 12 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

#### ▪ Sensitivity:

1.11 g/dl (11.1 g/l) – BS-400  
0.7 g/dl (7.0 g/l) – BS-480

#### ▪ Linearity:

up to 6.5 g/dl (65.0 g/l) – BS-400; BS-480

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### ▪ Specificity / Interferences

Haemoglobin up to 2.5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1200 mg/dl do not interfere with the test.

#### ▪ Precision

Repeatability (run to run)		Mean [g/dl]	SD [g/dl]	CV [%]
BS-400 (n = 20)	level 1	2.86	0.02	0.53
	level 2	4.47	0.01	0.21
BS-480 (n = 10)	level 1	2.89	0.01	0.50
	level 2	4.19	0.01	0.32
Reproducibility (day to day)		Mean [g/dl]	SD [g/dl]	CV [%]
BS-400 (n = 80)	level 1	2.87	0.03	1.08
	level 2	4.44	0.04	0.88
BS-480 (n = 20)	level 1	2.91	0.02	0.79
	level 2	4.23	0.03	0.67

#### ▪ Method comparison

A comparison between albumin values determined at BS-400 (y) and at ADVIA 1650 (x) using 50 samples gave following results:

$$y = 1.0231x + 0.0496 \text{ g/dl;} \\ R = 0.998 \quad (R - \text{correlation coefficient})$$

A comparison between albumin values determined at BS-480 (y) and at Cobas Integra 400 Plus (x) using 39 samples gave following results:

$$y = 0.8979x + 0.4783 \text{ g/dl;} \\ R = 0.990 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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• **Calibration**

<b>Calibration</b>	<b>Judgment Criteria</b>	<b>Blank Abs.</b>	
Rule	Sensitivity	Error Limit	
Replicate	Factor Diff.	Corr. Coeff.	
K	SD		

• **QC**

**Rules**

Westgard Multi-rule	Cum. Sum Check	<b>Auto QC</b>
v 1-2S	1.0 - 2.7	Interval
v 1-3S	• 1.0 - 3.0	
v 2-2S	0.5 - 5.1	

• **BS-480**

Chem	ALB	No.	003	Sample Type	SERUM
Chemistry	ALBUMIN	Print name	ALB		
Reaction Type	Endpoint	Reaction Direction	Increase		
Pri Wave	605	Sec Wave	740		
Unit	g/dl	Decimal	0.01		
Blank Time	11 12	Reaction Time	46 48		
Standard	Sample Vol 2.5 μL Aspirated 30 μL Diluent 120 μL	Reagent Vol	R1 250 μL R2 R3 R4	Diluent	
Decreased	2.5 μL				
Increased					
	Sample Blank v Auto Retun				
Linearity Range (Standard)	0.7 6.5	Linearity Limit			
Linearity Range (Decreased)		Substrate Depletion			
Linearity Range (Increased)		Mixed Blank Abs	-33000 33000		
R1 Blank Abs	-33000 33000	Uncapping Time	84 Day(s)		
Blank Response	-33000 33000	Reagent Alarm Limit			
Twin Chemistry		Enzyme Linear Extension			
	Prozone Check <input type="checkbox"/>	Rate Check <input type="checkbox"/>	Antigen Addition <input type="checkbox"/>		
Q1	0	Q2	0	Q3	0
PC	0	ABS	0	Q4	0

**Calibration Settings**

Math Model	Multi-point Linear	Bottle Changed	<input type="checkbox"/>
Factor		Lot Changed	<input type="checkbox"/>
	Replicates 3	Cal Time	<input type="checkbox"/>

**Acceptance Limits**

Cal Time	2016 Hour	SD	
Slope Diff		Repeatability	
Sensitivity			
Deter Coeff			

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