DILUENT PLUS

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INTENDED USE

DILUENT PLUS is intended to dilute blood samples before assay and to maintain appropriate environment during assay on hematology analyzers SYSMEX KX-21 and SYSMEX KX-21N. It is only for *in vitro* diagnostics, for healthcare professional users.

OPERATION PRINCIPLE

DILUENT PLUS is an isotonic reagent connected to the analyzer.

By diluting the blood sample, it allows it to flow through the measuring aperture of the analyzer and to make a measurement.

CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

preservatives	0,024 %
sodium chloride	
borate buffer	
metal ion chelator	

REAGENT STABILITY

The reagent must be used before the expiry date indicated on the label.

Once opened, the reagent is stable 90 days at room temperature.

WARNINGS AND NOTES

- Protect from direct sunlight!
- Do not use reagent beyond expiry date printed on the label.
- Do not combine the remaining reagent with the reagent from a new package.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Avoid discharge to sewage system or to environment.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

DEVICES INTENDED TO BE USED IN COMBINATION WITH THE REAGENT

DILUENT PLUS reagent is intended to be used in combination with the lysing reagent, FLUSH / FLUSH SET cleaning solution and appropriate hematology analyzer.

Name	Cat. No
LYSING REAGENT KX21 CN FREE	8-877
FLUSH (1 l)	8-832
FLUSH (500 ml)	8-891
FLUSH SET	8-896
SYSMEX KX-21	N/A
SYSMEX KX-21N	N/A

ASSAY PROCEDURE

DILENT PLUS is ready for use. Prior the first use, unopened reagent must be stored for minimum 24 hours at room temperature. Follow the analyzer's user manual.

WASTE MANAGEMENT¹

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Reagents residues: 18 01 07
- Empty packages: 15 01 02
- Wastewater from the analyzer: 18 01 03*

INCIDENTS²

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer (website address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

LITERATURE

- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on *in vitro* diagnostic medical devices.

LIST OF CHANGES

Previous version:	Current version:
51_03_13_006_01	51_03_13_019_01

Creation the new format of the document.

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