ENZYMATIC CLEANER FORTE

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INTENDED USE

ENZYMATIC CLEANER FORTE is intended to remove biological material's residues from the measurement system of the hematology analyzer MINDRAY BC-3000Plus. It is only for *in vitro* diagnostics, for healthcare professional users.

OPERATION PRINCIPLE

ENZYMATIC CLEANER FORTE is a washing reagent on board the analyzer. Removal of the biological material residues from the measurement system of the analyzer takes place due to the cleaning properties of the reagent. ENZYMATYC CLEANER FORTE contains surfactants and an enzyme that dissolve and remove large organic particles.

CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

SAVINASE	0,05 %
LIQUANASE	0,08 %
preservative	0,08 %
sodium sulfate	
sodium chloride	
sodium hydroxide	
hydrochloric acid	
detergent	
BRILLIANT BLUE R	

REAGENT STABILITY

The reagent must be used before the expiry date indicated on the label.

WARNING AND NOTES

- Protect from direct sunlight!
- Do not use reagent beyond expiry date printed on the label.
- Do not combine the remaining reagent with the reagent from a new package.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Avoid discharge to sewage system or to environment.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- EUH208 Contains subtilisin. May produce an allergic reaction.
- EUH210 Safety data sheet available on request.

DEVICES INTENDED TO BE USED IN COMBINATION WITH THE REAGENT

ENZYMATIC CLEANER FORTE is intended to be used in combination with the lysing reagent, diluent, cleaner and appropriate hematology analyzer.

Name	Cat. No
DILUENT	8-892
LYSING REAGENT CN FREE	8-893
CLEANING REAGENT	8-894
MINDRAY BC-3000Plus	N/A

ASSAY PROCEDURE

ENZYMATIC CLEANER FORTE is ready for use. Prior the first use, unopened reagent must be stored for minimum 24 hours at room temperature. Follow the analyzer's operator manual.

WASTE MANAGEMENT¹

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Reagents residues: 18 01 07
- Empty packages: 15 01 02
- Wastewater from the analyzer: 18 01 03*

INCIDENTS²

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer (websites address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

LITERATURE

- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on *in vitro* diagnostic medical devices.

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