Cystatin C

Diagnostic reagent for the quantitative in vitro determination of Cystatin C (CysC) in human serum or plasma by turbidimetric assay

REF	Content
906840	2x 25 mL Buffer + 1x 10 mL Latex Reagent
906850	1x 25 mL Buffer + 1x 5 mL Latex Reagent
909860	1x 12.5 mL Buffer + 1x 2.5 mL Latex Reagent
9A1840	1x 50 mL Buffer + 1x 12.5 mL Latex Reagent
9T1340	1x 50 mL Buffer + 1x 12.5 mL Latex Reagent

Additionally offered:

909890	1x 2 mL Cystatin C Calibrator High
906880	6x 2 mL Cystatin C Control Set (2 Levels)

GENERAL INFORMATION

Method	Immunoturbidimetric
Reaction	2 Point End
Wavelength	600 nm
Assay Temperature	37 °C
Sample	Serum, Plasma
Measuring Range	0.0 – 10.0 mg/L
Sensitivity	0.08 mg/L
Hook Effect	> 10 mg/L

Test Procedure	Tests/Kit
906840	200 (Hitachi 911)
906850	100 (Hitachi 911)

Automated Test Procedure

Instrument dependent – please ask for applications

REAGENT COMPOSITION

COMPONENTS	FINAL CONCENTRATION
Cystatin C Latex Reagent	
Goat anti human CysC antibody sensitized latex	0.5 %
Sodium azide	0.95 g/L
Glycine buffer, pH 7.3	
Cystatin C Buffer	
Phosphate buffered saline, pH 7.43	
Sodium azide	0.95 g/L
Detergent	0.5 %

REAGENT PREPARATION

The reagents are liquid and ready to use.

REAGENT STABILITY AND STORAGE

Conditions: Protect from light. Close immediately after use.

Stability: at 2 – 8 °C up to the expiration date

Do not freeze!

SAMPLE STABILITY AND STORAGE

Stability: at 2 – 8 °C 48 hours
 at – 20 °C 3 months

Freeze only once!

TEST PROCEDURE ON HITACHI 911

TEST NAME	CysC
ASSAY CODE	[2-POINT END];[10]-[0]
ASSAY POINT	[15] – [31] – [0] – [0]
SAMPLE VOLUME	[2]
R1 VOLUME	[230]
R2 VOLUME	[40]
WAVELENGTH (2 nd /Prim.)	[] [600]
STD (1) CONC-POS	[0.0*] [1] [2] [0] [0]
STD (2) CONC-POS	[**] [2] [5] [5] [195]
STD (3) CONC-POS	[**] [3] [10] [5] [190]
STD (4) CONC-POS	[**] [4] [20] [5] [180]
STD (5) CONC-POS	[**] [5] [20] [10] [180]
STD (6) CONC-POS	[**] [6] [2] [0] [0]
SD LIMIT	[999]
DUPLICATE LIMIT	[1000]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[32000] [INCREASE]
PROZONE LIMIT	[] [-]
EXPECTED VALUE	[**] [**]
PANIC VALUE	[] [-]
INSTRUMENT FACTOR	[1.0]

* use saline for the 0.0 calibrator

** input the values of the calibrator being used

CALCULATION

The Cystatin C concentration of unknown samples is derived from the calibration curve using an appropriate mathematical model such as spline. The calibration curve is obtained with calibrator set and NaCl solution (9 g/l) as zero value.

REFERENCE RANGE

0.4 – 1.2 mg/L

It is recommended that each laboratory establishes its own normal range.

TEST PRINCIPLE

The assay of Cystatin C is based on turbidimetric measurement. Turbidity is caused by the formation of antigen – antibody insoluble immuno complexes.

DIAGNOSTIC IMPLICATIONS

Cystatin C is a cysteine proteinase inhibitor with a relative molecular weight of 13,250 daltons and is formed by all nucleated cells investigated. Since it is formed at a constant rate and freely filtered by the healthy kidney, this protein is a good marker of renal function. Serum concentrations of Cystatin C are almost totally dependent on the glomerular filtration rate. A reduction in the glomerular filtration rate (GFR) causes a rise in the concentration of Cystatin C. Cystatin C has not been shown to be affected by factors such as muscle mass and nutrition, factors which have been demonstrated to affect creatinine values. In addition, a raise in creatinine does not become evident until the GFR has fallen by approximately 50%.

PERFORMANCE CHARACTERISTICS

MEASURING RANGE

0.0 – 10.0 mg/L

DETECTION LIMIT

0.08 mg/L

HOOK EFFECT

> 10 mg/L

SPECIFICITY

Monospecific

SENSITIVITY

0.0306 ABS units/concentration units

PRECISION [%CV]

	Low	Medium	High
Intra-Run	4,86	3,35	3,10
Inter-Run		2,48	3,8

ACCURACY [mg/L]

Control	Assigned	Measured
Behring L1	0,97 (0,82 – 1,11)	1,08
Behring L2	1,91 (1,62 – 2,20)	2,19

METHOD COMPARISON

A comparison with Behring Nephelometry gave the following results:
 $y = 0,9009x + 0,1108 / r = 0,9626$

INTERFERING SUBSTANCES

No interference up to:

Triglycerides	5000 mg/dL	Hemoglobin	500 mg/dL
Bilirubin	60 mg/dL	Sodium Citrate	2000 mg/dL
Heparin	100 mg/dL	EDTA	10 mg/mL

QUALITY CONTROL

All Control sera with Cystatin C values measured by this method may be used. We recommend the Dialab Cystatin C Control Set.

CALIBRATION

The assay requires the use of Cystatin C Calibrators. We recommend the Dialab Cystatin C Calibrator High.

AUTOMATION

Applications for automated systems are available upon request.

WARNINGS AND PRECAUTIONS

- The Cystatin C reagents are intended for in vitro diagnostic use only.
- Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion.
- Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV antibodies, as well as for Hepatitis B surface antigen, using a method approved by the FDA

WASTE MANAGEMENT

Please refer to local requirements.

REFERENCES

- Simonsen O, Grubb A, Thysell H. The blood serum concentration of Cystatin C (gamma-trace) as a measure of the glomerular filtration rate. Scand J Clin Lab Invest 1985; 45: 97-101.
- Randers E, erlandsen EJ. Serum cystatin C as an endogenous marker of the renal function – a review. Clin chem. Lab Med 1999; 37: 389-95.
- Grubb A. Diagnostic value of analysis of cystatin C and protein HC in biological fluids. Clin Nephrol 1992; 38 Suppl 1: S20-7.
- Coll E, Botey A, Alvarez L, et al. Serum cystatin C as a new marker of glomerular filtration rate and as a marker for early renal impairment. Am J Kidney Dis 2000; 36: 39-34

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BIBLIOGRAPHIE

- Simonsen O, Grubb A, Thysell H. The blood serum concentration of Cystatin C (gamma-trace) as a measure of the glomerular filtration rate. Scand J Clin Lab Invest 1985; 45: 97-101.
- Randers E, erlandsen EJ. Serum cystatin C as an endogenous marker of the renal function – a review. Clin chem. Lab Med 1999; 37: 389–95.
- Grubb A. Diagnostic value of analysis of cystatin C and protein HC in biological fluids. Clin Nephrol 1992; 38 Suppl 1: S20-7.
- Coll E, Botey A, Alvarez L, et al. Serum cystatin C as a new marker of glomerular filtration rate and as a marker for early renal impairment. Am J Kidney Dis 2000; 36: 39-34