Lyophilised Control Serum

LITHIUM CONTROL SET (Assayed Control Serum) 2 levels

Human based serum for quality control of quantitative in vitro determination of lithium in human serum on photometric systems.

REF

Cont.

910290 **2 x 3 ml** control set, 2 levels

Additionally offered:

910210 3 x 20 ml Lithium reagent kit

CONTROL COMPOSITION

The controls consist of lithium in human serum.

CONTROL PREPARATION

Open vial carefully, avoiding any loss of material, and reconstitute with exactly 3.0 mL of distilled water. Close vial and allow standing at room temperature for 30 minutes, dissolving contents completely by swirling or rotating gently.

CONTROL STABILITY AND STORAGE

	Temperature:	Stability:
Lyophilized:	2 – 8 °C	up to expiration date
After reconstitution:	2 – 8 °C	1 week

CLOSE IMMEDIATELY AFTER USE
PROTECT FROM LIGHT!

WARNINGS AND PRECAUTIONS

Each donor unit of serum in the preparation of this material was tested by FDA-approved methods and found negative for the Human Immunodeficiency Virus Antibody (HIV I/II Ab), Hepatitis B Surface Antigen (HBsAg), and Hepatitis C Virus Antibody (HCV). Because no method can offer complete assurance as to the absence of infectious agents, however, this material and all patient samples should be handled as though capable of transmitting infectious disease and disposed of accordingly.

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