α-1 Microglobulin

Diagnostic reagent for the quantitative in vitro determination of α-1 Microglobulin in human urine by turbidimetric assay

REF	Content
A00503	1x 5 mL α-1 Microglobulin Latex Reagent 2x 25 mL α-1 Microglobulin Buffer
Additionally offered:	
A00737	1x 1 mL α-1 Microglobulin Calibrator
A00738	1x 5 mL α-1 Microglobulin Calibrator
A01810	1x 1 mL α-1 Microglobulin Control

GENERAL INFORMATION

Method	Immunoturbidimetric
Reaction	Nonlinear, endpoint
Wavelength	600 nm
Assay Temperature	18 – 37 °C
Sample	Urine
Measuring Range	approx. 0 – 50 mg/L
Sensitivity	1 mg/L (Cobas Mira)
Hook Effect	without sample dilution: > 500 mg/L (Cobas Mira) with sample dilution: > 500 mg/L (Cobas Mira)

Manual Test Procedure	Tests/Kit*
without sample dilution	25
with sample dilution	50

Automated Test Procedure
Instrument dependent – please ask for applications

* calculated on amount of antibody reagent; additional buffer on request

REAGENT COMPOSITION

COMPONENTS	FINAL CONCENTRATION
α-1 Microglobulin Latex Reagent	
Polystyrene latex particles sensitized with α-1 Microglobulin antibodies	0.5 %
Sodium azide	0.095 %
α-1 Microglobulin Buffer	
Phosphate buffered saline	
Sodium azide	0.095 %
Detergent	0.1 %

REAGENT PREPARATION

The reagents are liquid and ready to use.

REAGENT STABILITY AND STORAGE

Conditions:	Protect from light. Close immediately after use.	
Stability:	at 2 – 8 °C	up to the expiration date
	at 18 – 25 °C	1 month

Do not freeze!

SAMPLE STABILITY AND STORAGE

Use fresh, centrifuged urine. Freeze only once!

Stability:	at 2 – 8 °C	48 hours
	at – 20 °C	3 months

MANUAL TEST PROCEDURE

Test Procedure without Sample Dilution:

Samples/Controls: ready to use
 Calibration curve: Use α-1 Microglobulin Calibrator to generate a calibration curve by making 1:2 serial dilutions of the calibrator with 0.9% saline as diluent. Use 0.9% saline as zero point.

Pipette into test tubes:	Calibrators	Samples/Controls
Buffer	800 µL	800 µL
Cal./Ctrls/Samples	8 µL	8 µL
Mix. Read A1 of calibrators and samples/controls at 600 nm. Then add:		
Latex Reagent	200 µL	200 µL
Mix. Incubate 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 600 nm. Calculate: $\Delta A = (A2 - A1)$		

Test Procedure with Sample Dilution:

Sample/ Control: dilute 1:10 in saline 0.9%
 Calibration curve: Use α-1 Microglobulin Calibrator to generate a calibration curve by making 1:10, 1:20, 1:40, 1:80 dilutions with 0.9% saline as diluent. Use 0.9% saline as zero point.

Pipette into test tubes	Calibrators	Samples/Controls
Buffer	900 µL	900 µL
Cal./Ctrls/Samples	50 µL	50 µL
Mix. Read A1 of calibrators and samples/controls at 600 nm. Then add:		
Latex Reagent	100 µL	100 µL
Mix. Incubate 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 600 nm. Calculate: $\Delta A = (A2 - A1)$		

CALCULATION

Calculate and plot $\Delta A = (A2 - A1)$ of the calibrators versus assigned concentration values on a linear-linear graph paper. Calculate ΔA optical densities of samples and control(s) and read values in mg/L on the reference curve. Samples yielding absorbances above highest calibrator should be retested after further dilution.

REFERENCE RANGE

< 12 mg/L

It is recommended that each laboratory establishes its own normal range.

TEST PRINCIPLE

The assay of α-1 Microglobulin is based on turbidimetric measurement. Turbidity is caused by the formation of antigen-antibody insoluble immuno complexes.

DIAGNOSTIC IMPLICATIONS

α-1 Microglobulin is a low molecular weight glycoprotein of approximately 33,000 daltons. It is mainly synthesized in the liver and widely distributed in various body fluids. The clinical relevance of urine α-1 Microglobulin relates to the identification of tubular proteinurias. Indeed, AMI is filtered through the glomeruli; reabsorption and catabolism take place in the proximal tubule. Elevated urinary concentrations of α-1 Microglobulin are indicative of tubular damage as may occur in nephritides, advanced diabetic nephropathy, after exposure of heavy metals or after administration of nephrotoxic drugs.

PERFORMANCE CHARACTERISTICS

SENSITIVITY

1 mg/L (Cobas Mira)

ACCURACY

A control from Behring was measured on the Cobas Mira to verify proper assay recovery.

Control	Assigned Value (mg/L)	Measured Value (mg/L)
Behring	33.0 (28.0 – 38.0)	33.4

PRECISION

Intra-Assay Precision

Two urine samples were consecutively measured 20 times on the Cobas Mira.

Expected Value	n	Mean	S.D.	C.V
Normal	20	6.74	0.202	2.99
Pathological	20	17.10	0.30	1.78

Inter-Assay Precision

Two samples were measured at regularly time intervals on the Cobas Mira. The samples were stored at 4 °C.

Sample	n	Mean	S.D.	C.V.
1	22	11.72	0.684	5.84
2	22	30.47	0.848	2.78

METHOD COMPARISON

A comparison with Nephelometry gave the following results:
 $y = 0.9072x + 0.1001$; $r = 0.997$

INTERFERING SUBSTANCES

No interference up to:

NH ₄ Cl	400 mg/dL	Hemoglobin	540 mg/dL
Bilirubin	30 mg/dL	Ascorbic Acid	50 mg/dL
Turbidity	5 %		

QUALITY CONTROL

All Control sera with α-1 Microglobulin values measured by this method may be used. We recommend the Dialab α-1 Microglobulin Control.

CALIBRATION

The assay requires the use of α-1 Microglobulin Calibrators. We recommend the Dialab α-1 Microglobulin Calibrator.

AUTOMATION

Applications for automated systems (with and without sample dilution) are available upon request.

WARNINGS AND PRECATIONS

- The α-1 Microglobulin reagents are intended for in vitro diagnostic use only.
- Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion.
- Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV antibodies, as well as for Hepatitis B surface antigen, using a method approved by the FDA

WASTE MANAGEMENT

Please refer to local requirements.

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