## Antithrombin III

Diagnostic reagent for the quantitative in vitro determination of Antithrombin III (ATIII) in human plasma by turbidimetric assay

### REF

### Content

**A00510** 1x 10 mL ATIII Antibody Reagent  
 5x 25 mL PEG6 Buffer

Additionally offered:

A00580 1x 1 mL Protein Calibrator High  
 A00703 1x 5 mL Protein Calibrator High  
 A00701 1x 1 mL Protein Calibrator Low  
 A00702 1x 5 mL Protein Calibrator Low  
 A00590 1x 1 mL Protein Control  
 A00800 1x 5 mL Protein Control  
 A08591 1x 1 mL Protein Control Low  
 A08823 1x 5 mL Protein Control Low

### GENERAL INFORMATION

**Method** Immunoturbidimetric  
**Reaction** Nonlinear, endpoint  
**Wavelength** 340 nm  
**Assay Temperature** 18 – 37 °C  
**Sample** Plasma  
**Measuring Range** approx. 0 – 65 mg/dL  
**Sensitivity** 8 mg/dL  
**Hook Effect** with sample dilution: > 300 mg/dL  
**Manual Test Procedure** Tests/Kit\*  
 with sample dilution 166

### Automated Test Procedure

Instrument dependent – please ask for applications

\* calculated on amount of antibody reagent; additional buffer on request

### REAGENT COMPOSITION

COMPONENTS	FINAL CONCENTRATION
<b>Antithrombin III Antibody Reagent</b>	
Turbidimetric grade antibody raised in the goat, monospecific for ATIII	variable
Sodium azide	0.095 %
<b>PEG6 Buffer</b>	
Phosphate buffered saline	
PEG	6 %
Sodium azide	0.095 %

### REAGENT PREPARATION

The reagents are liquid and ready to use.

### REAGENT STABILITY AND STORAGE

Conditions: Protect from light. Close immediately after use.  
 Stability: at 2 – 8 °C up to the expiration date  
 at 18 – 25 °C 1 month

Do not freeze!

### SAMPLE STABILITY AND STORAGE

Stability: at 2 – 8 °C 1 week  
 at – 20 °C 3 months

Freeze only once!

### MANUAL TEST PROCEDURE

#### Test Procedure with Sample Dilution:

Sample/ Control: dilute 1:10 in saline 0.9%

Calibration curve: Use Protein Calibrator High to generate a calibration curve by making 1:10, 1:20, 1:40, 1:80 and 1:160 dilutions with 0.9% saline as diluent. Use 0.9% saline as zero point.

Pipette into test tubes	Calibrators	Samples/Controls
Buffer	900 µL	900 µL
Cal./Ctrls/Samples	100 µL	100 µL
Mix. Read A1 of calibrators and samples/controls at 340 nm. Then add:		
Antibody Reagent	60 µL	60 µL
Mix. Incubate 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 340 nm. Calculate: $\Delta A = (A2 - A1)$		

### CALCULATION

Calculate and plot  $\Delta A = (A2 - A1)$  of the calibrators versus assigned concentration values on a linear-linear graph paper. Calculate  $\Delta A$  optical densities of samples and control(s) and read values in mg/dL on the reference curve. Samples yielding absorbances above highest calibrator should be retested after further dilution.

### REFERENCE RANGE

22 – 39 mg/dL

It is recommended that each laboratory establishes its own normal range.

### TEST PRINCIPLE

The assay of Antithrombin III is based on turbidimetric measurement. Turbidity is caused by the formation of antigen-antibody insoluble immuno complexes. The formation of the complexes is accelerated and enhanced by PEG.

### DIAGNOSTIC IMPLICATIONS

ATIII is an inhibitor of thrombin, Factor Xa or Factor VIIa. There is more ATIII in blood than prothrombin; blood is able to clot only because the reaction of the inhibitor with thrombin is much slower than the action of thrombin on fibrinogen. The reactivity of ATIII is regulated by combination with its activator heparin. Decreased levels of ATIII are found in nephrotic syndrome, DIC, deep vein thrombosis, hypercoagulability syndrome, oral contraceptives, pulmonary embolism, direct hepatotoxicity, extra-corporeal circulation and L-asparaginase administration.

### PERFORMANCE CHARACTERISTICS

#### SENSITIVITY

8 mg/dL (Cobas Mira)

#### ACCURACY

Controls are assayed in duplicate on Cobas Mira.

Control	Assigned Value (mg/dL)	Measured Value (mg/dL)
Biorad 1	17 (14 – 21)	20
Biorad 2	47 (38 – 57)	51

### PRECISION

#### Intra-Assay Precision

The Antithrombin III concentration (mg/dl) of 3 unknown serum samples was measured and the variation coefficient determined.

Concentration	n	Mean (mg/dL)	S.D.	C.V
Low	20	17.81	0.51	2.84
Medium	20	34.95	0.71	2.04
High	19	251.74	2.55	4.93

### INTERFERING SUBSTANCES

No interference up to:

Sodium-Citrate	1000 mg/dL	Hemoglobin	1000 mg/dL
Triglycerides	2500 mg/dL	Heparin	50 mg/dL
Bilirubin	20 mg/dL		

### QUALITY CONTROL

All Control sera with Antithrombin III values measured by this method may be used. We recommend the Dialab Protein Control and the Protein Control Low.

### CALIBRATION

The assay requires the use of Antithrombin III serum Calibrators. We recommend the Dialab Protein Calibrator High or the Protein Calibrator Low.

### AUTOMATION

Applications for automated systems (with and without sample dilution) are available upon request.

### WARNINGS AND PRECAUTIONS

- The Antithrombin III reagents are intended for in vitro diagnostic use only.
- Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion.
- Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV antibodies, as well as for Hepatitis B surface antigen, using a method approved by the FDA

### WASTE MANAGEMENT

Please refer to local requirements.

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