β-2 Microglobulin

Diagnostic reagent for the quantitative in vitro determination of β-2 Microglobulin in human serum and urine by turbidimetric assay

REF	Content
A00511	1x 10 mL β-2 Microglobulin Latex Reagent 5x 25 mL β-2 Microglobulin Buffer
Additionally offered:	
A00735	4x 1 mL β-2 Microglobulin Calibrator 4 level series
A00818	1x 1 mL β-2 Microglobulin Control
A00819	1x 5 mL β-2 Microglobulin Control
A00590	1x 1 mL Protein Control
A00800	1x 5 mL Protein Control
A08591	1x 1 mL Protein Control Low
A08823	1x 5 mL Protein Control Low

GENERAL INFORMATION

Method	Immunoturbidimetric
Reaction	Nonlinear, endpoint
Wavelength	600 nm
Assay Temperature	18 – 37 °C
Sample	Serum, centrifuged urine
Measuring Range	approx. 0 – 11 mg/L
Sensitivity	0.15 mg/L (Cobas Mira)
Hook Effect	no risk
Manual Test Procedure	Tests/Kit* 66

Automated Test Procedure

Instrument dependent – please ask for applications

* calculated on amount of antibody reagent; additional buffer on request

REAGENT COMPOSITION

COMPONENTS	FINAL CONCENTRATION
β-2 Microglobulin Latex Reagent	
Suspension of polystyrene latex particles of uniform size coated with β-2 Microglobulin	0.17 %
Sodium azide	0.095 %
β-2 Microglobulin Buffer	
Phosphate buffered saline	
Detergent	0.1 %
Sodium azide	0.095 %

REAGENT PREPARATION

The reagents are liquid and ready to use.

REAGENT STABILITY AND STORAGE

Conditions:	Protect from light. Close immediately after use.	
Stability:	at 2 – 8 °C	up to the expiration date
	at 18 – 25 °C	1 month

Do not freeze!

SAMPLE STABILITY AND STORAGE

Use fresh serum or fresh, centrifuged urine. Freeze only once!

Stability:	at 2 – 8 °C	48 hours (serum and urine)
	at – 20 °C	3 months (serum and urine)

MANUAL TEST PROCEDURE

Sample/ Control: ready to use

Calibration curve: Use β-2 Microglobulin Calibrator Series to generate a calibration curve. Use saline as zero point.

Pipette into test tubes	Calibrators	Samples/Controls
Buffer	990 µL	990 µL
Cal./Ctrls/Samples	6 µL	6 µL
Mix. Read A1 of calibrators and samples/controls at 600 nm. Then add:		
Latex Reagent	150 µL	150 µL
Mix. Incubate 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 600 nm. Calculate: $\Delta A = (A2 - A1)$		

CALCULATION

Calculate and plot $\Delta A = (A2 - A1)$ of the calibrators versus assigned concentration values on a linear-linear graph paper. Calculate ΔA optical densities of samples and control(s) and read values in mg/L on the reference curve. Samples yielding absorbances above highest calibrator should be retested after further dilution.

REFERENCE RANGE

0.8 – 1.8 mg/L (in serum)
 < 0.5 mg/L (in urine)
 It is recommended that each laboratory establishes its own normal range.

TEST PRINCIPLE

The assay of β-2 Microglobulin is based on turbidimetric measurement. Turbidity is caused by the formation of antigen-antibody insoluble immuno complexes. The formation of the complexes is accelerated and enhanced by PEG.

DIAGNOSTIC IMPLICATIONS

β-2 Microglobulin is a low molecular weight (11,000 Daltons) protein found on the membranes of virtually all body cells. Free BMG is a product of cell breakdown. It is secreted by the renal glomeruli, then absorbed and catabolized by the renal tubular cells. Decreased glomerular filtration is associated with high serum levels of BMG, whereas tubular insufficiency is associated with normal serum and high urine levels. Markedly increased cell breakdown, as in acute leukaemia, may also be associated with high serum levels.

PERFORMANCE CHARACTERISTICS

SENSITIVITY

0.15 mg/L (Cobas Mira)

ACCURACY

Controls are assayed in duplicate on Cobas Mira

Control	Assigned Value (mg/L)	Measured Value (mg/L)
DIALAB	2.75 (2.34 – 3.16)	2.86
Dade Behring	1.94 (1.65 – 2.23)	2.07

PRECISION

Intra-Assay Precision

3 Serum Samples were consecutively measured on the Cobas Mira.

Sample	C.V
Low	3.8
Medium	2.95
High	1.52

METHOD COMPARISON

A comparison with RIA gave the following results:
 $y = 1.4725x - 0.242$; $r = 0.0469$

QUALITY CONTROL

All commercially available Control sera with β-2 Microglobulin values measured by this method may be used. We recommend the Dialab β-2 Microglobulin Control, Protein Control and the Protein Control Low.

CALIBRATION

The assay requires the use of β-2 Microglobulin serum Calibrators. We recommend the Dialab β-2 Microglobulin Calibrator 4 level series.

AUTOMATION

Applications for automated systems (with and without sample dilution) are available upon request.

WARNINGS AND PRECAUTIONS

- The β-2 Microglobulin reagents are intended for in vitro diagnostic use only.
- Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion.
- Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV antibodies, as well as for Hepatitis B surface antigen, using a method approved by the FDA

WASTE MANAGEMENT

Please refer to local requirements.

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