Complement C4

Diagnostic reagent for the quantitative in vitro determination of Complement C4 in human serum by turbidimetric assay

REF

Content

A00514 1x 10 mL C4 Antibody Reagent
 5x 25 mL PEG4 Buffer

Additionally offered:

A00704 5x 1 mL Protein Calibrator 5 Level Series
 A00580 1x 1 mL Protein Calibrator High
 A00703 1x 5 mL Protein Calibrator High
 A00701 1x 1 mL Protein Calibrator Low
 A00702 1x 5 mL Protein Calibrator Low
 A00590 1x 1 mL Protein Control
 A00800 1x 5 mL Protein Control
 A08591 1x 1 mL Protein Control Low
 A08823 1x 5 mL Protein Control Low

GENERAL INFORMATION

Method Immunoturbidimetric
Reaction Nonlinear, endpoint
Wavelength 340 nm
Assay Temperature 18 – 37 °C
Sample Serum
Measuring Range approx. 0 – 80 mg/dL
Sensitivity 2 mg/dL
Hook Effect > 600 mg/dL

Manual Test Procedure Tests/Kit*
 without sample dilution 142

Automated Test Procedure

Instrument dependent – please ask for applications

* calculated on amount of antibody reagent; additional buffer on request

REAGENT COMPOSITION

COMPONENTS FINAL CONCENTRATION

Complement C4 Antibody Reagent

Turbidimetric grade antibody raised in the goat,
 monospecific for C4
 Sodium azide 0.095 %

PEG4 Buffer

Phosphate buffered saline
 PEG 4 %
 Sodium azide 0.095 %

REAGENT PREPARATION

The reagents are liquid and ready to use.

REAGENT STABILITY AND STORAGE

Conditions: Protect from light. Close immediately after use.
Stability: at 2 – 8 °C up to the expiration date
 at 18 – 25 °C 1 month

Do not freeze!

SAMPLE STABILITY AND STORAGE

Stability: at 2 – 8 °C 48 hours
 at – 20 °C 3 months

Freeze only once!

MANUAL TEST PROCEDURE

Test Procedure without Sample Dilution:

Samples/Controls: ready to use

Calibration curve: Use Protein Calibrator High to generate a calibration curve by making 1:2 serial dilutions of the calibrator with 0.9% saline as diluent or use the 5 level calibrator series. Use 0.9% saline as zero point.

Pipette into test tubes:	Calibrators	Samples/Controls
Buffer	900 µL	900 µL
Cal./Ctrls/Samples	20 µL	20 µL
Mix. Read A1 of calibrators and samples/controls at 340 nm. Then add:		
Antibody Reagent	70 µL	70 µL
Mix. Incubate 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 340 nm. Calculate: $\Delta A = (A2 - A1)$		

CALCULATION

Calculate and plot $\Delta A = (A2 - A1)$ of the calibrators versus assigned concentration values on a linear-linear graph paper. Calculate ΔA optical densities of samples and control(s) and read values in mg/dL on the reference curve. Samples yielding absorbances above highest calibrator should be retested after further dilution.

REFERENCE RANGE

9 – 36 mg/dL

It is recommended that each laboratory establishes its own normal range.

TEST PRINCIPLE

The assay of C4 is based on turbidimetric measurement. Turbidity is caused by the formation of antigen-antibody insoluble immuno complexes. The formation of the complexes is accelerated and enhanced by PEG.

DIAGNOSTIC IMPLICATIONS

C4 is a constituent of C3 convertase and C5 convertase. Decreased levels are found in hereditary angioneurotic oedema, immune complex diseases and congenital deficiencies.

PERFORMANCE CHARACTERISTICS

SENSITIVITY

2 mg/dL (Cobas Mira)

ACCURACY

Controls were assayed in duplicate.

Control	Assigned Value (mg/dL)	Measured Value (mg/dL)
Liquicheck 1 (BIORAD)	16 (13 - 19)	16.7
Liquicheck 2 (BIORAD)	52 (42 - 62)	47.7
Seronorm L (NYCOMED)	10 (8 - 12)	9.4
Seronorm N (NYCOMED)	34 (27 - 41)	31.9
Seronorm H (NYCOMED)	45 (36 - 54)	43.0
Immunology 1 (CIBA CORNING)	17.7 (14.2 - 21.2)	19.3
Immunology 2 (CIBA CORNING)	35.3 (28.2 - 42.4)	38.2

PRECISION

Intra-Assay Precision

3 Serum Samples were consecutively measured 20 times on the Cobas Mira.

Expected Value	n	Mean	S.D.	C.V.
Low	20	9.9	0.45	4.54
Medium	20	26.3	0.57	2.18
High	20	52.9	2.10	3.97

Inter-Assay Precision

2 Control sera were measured daily on the Mega analyzer after calibration.

Expected Value	n	Mean	S.D.	C.V.
Ortho 1	64	21.4	0.66	3.08
Ortho 2	64	13.5	0.56	4.17

METHOD COMPARISON

A comparison with Nephelometry gave the following results:
 $y = 0.93418x - 0.1518$; $r = 0.9940$

INTERFERING SUBSTANCES

No interference up to:

Triglycerides	2500 mg/dL	Hemoglobin	1000 mg/dL
Bilirubin	20 mg/dL	Sodium Citrate	1000 mg/dL
Heparin	50 mg/dL	Turbidity	5 %

QUALITY CONTROL

All Control sera with C4 values measured by this method may be used. We recommend the Dialab Protein Control and the Protein Control Low.

CALIBRATION

The assay requires the use of C4 serum Calibrators. We recommend the Dialab Protein Calibrator 5 Level Series, the Protein Calibrator High or the Protein Calibrator Low.

AUTOMATION

Applications for automated systems (with and without sample dilution) are available upon request.

WARNINGS AND PRECATIONS

- The C4 reagents are intended for in vitro diagnostic use only.
- Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion.
- Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV antibodies, as well as for Hepatitis B surface antigen, using a method approved by the FDA

WASTE MANAGEMENT

Please refer to local requirements.

REFERENCES

- Dati, F et al., Lab. Med. 13, 87 (1989)
- Müller-Eberhard, H.H., Ann. Rev. Biochem. 44, 697 (1975)
- Lachmann, P.J., Hobart, M.J. and Ashton, W.P. (1973) in Handbook of Experimental Immunology, 2nd Ed., 16, Ed. D.M. Weir, Blackwell Scientific Publications

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