# Transferrin

Diagnostic reagent for the quantitative in vitro determination of Transferrin (TRF) in human serum by turbidimetric assay

REF	Content
<b>A00515M</b>	1x 10 mL Transferrin Antibody Reagent 5x 25 mL PBS Buffer <i>automated Test Procedure without Sample Dilution</i> <i>manual Test Procedure with/without Sample Dilution</i>
<b>A00515A</b>	1x 10 mL Transferrin Antibody Reagent 5x 25 mL PEG4 Buffer <i>automated Test Procedure with Sample Dilution</i>
Additionally offered:	
A00960	5x 1 mL Transferrin Calibrator 5 Level Series
A00580	1x 1 mL Protein Calibrator High
A00703	1x 5 mL Protein Calibrator High
A00701	1x 1 mL Protein Calibrator Low
A00702	1x 5 mL Protein Calibrator Low
A00590	1x 1 mL Protein Control
A00800	1x 5 mL Protein Control
A08591	1x 1 mL Protein Control Low
A08823	1x 5 mL Protein Control Low

## GENERAL INFORMATION

<b>Method</b>	Immunoturbidimetric
<b>Reaction</b>	Nonlinear, endpoint
<b>Wavelength</b>	340 nm
<b>Assay Temperature</b>	18 – 37 °C
<b>Sample</b>	Serum
<b>Measuring Range</b>	approx. 0 – 500 mg/dL
<b>Sensitivity</b>	40 mg/dL (Cobas Mira)
<b>Hook Effect</b>	without sample dilution: > 2,500 mg/dL (Cobas Mira) with sample dilution: > 4,000 mg/dL (Cobas Mira)
<b>Manual Test Procedure</b>	<b>Tests/Kit*</b>
without sample dilution	125
with sample dilution	250
<b>Automated Test Procedure</b>	Instrument dependent – please ask for applications
* calculated on amount of antibody reagent; additional buffer on request	

## REAGENT COMPOSITION

COMPONENTS	FINAL CONCENTRATION
<b>Transferrin Antibody Reagent</b>	
Turbidimetric grade antibody raised in the goat, monospecific for Transferrin	
Sodium azide	variable
	0.095 %
<b>PBS Buffer (A00515M)</b>	
Phosphate buffered saline	
Sodium azide	0.095 %
<b>PEG4 Buffer (A00515A)</b>	
Phosphate buffered saline	
PEG	4 %
Sodium azide	0.095 %

## REAGENT PREPARATION

The reagents are liquid and ready to use.

## REAGENT STABILITY AND STORAGE

<b>Conditions:</b>	Protect from light. Close immediately after use.
<b>Stability:</b>	at 2 – 8 °C up to the expiration date at 18 – 25 °C 1 month

Do not freeze!

## SAMPLE STABILITY AND STORAGE

<b>Stability:</b>	at 2 – 8 °C 48 hours at – 20 °C 3 months
-------------------	---

Freeze only once!

## MANUAL TEST PROCEDURE (FOR REF A00515M)

### Test Procedure without Sample Dilution:

Samples/Controls: ready to use

Calibration curve: Use Protein Calibrator High to generate a calibration curve by making 1:2 serial dilutions of the calibrator with 0.9% saline as diluent or use the 5 level calibrator series. Use 0.9% saline as zero point.

Pipette into test tubes:	Calibrators	Samples/Controls
Buffer	900 µL	900 µL
Cal./Ctrls/Samples	5 µL	5 µL
Mix. Read A1 of calibrators and samples/controls at 340 nm. Then add:		
Antibody Reagent	80 µL	80 µL
Mix. Incubate 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 340 nm. Calculate: $\Delta A = (A2 - A1)$		

### Test Procedure with Sample Dilution:

Sample/ Control: dilute 1:10 in saline 0.9%

Calibration curve: Use Protein Calibrator High to generate a calibration curve by making 1:10, 1:20, 1:40, 1:80 and 1:160 dilutions with 0.9% saline as diluent. Use 0.9% saline as zero point.

Pipette into test tubes	Calibrators	Samples/Controls
Buffer	900 µL	900 µL
Cal./Ctrls/Samples	25 µL	25 µL
Mix. Read A1 of calibrators and samples/controls at 340 nm. Then add:		
Antibody Reagent	40 µL	40 µL
Mix. Incubate 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 340 nm. Calculate: $\Delta A = (A2 - A1)$		

## CALCULATION

Calculate and plot  $\Delta A = (A2 - A1)$  of the calibrators versus assigned concentration values on a linear-linear graph paper. Calculate  $\Delta A$  optical densities of samples and control(s) and read values in mg/dL on the reference curve. Samples yielding absorbances above highest calibrator should be retested after further dilution.

## REFERENCE RANGE

Men & Women: 170 – 340 mg/dL

It is recommended that each laboratory establishes its own normal range.

## TEST PRINCIPLE

The assay of Transferrin is based on turbidimetric measurement. Turbidity is caused by the formation of antigen-antibody insoluble immuno complexes. The formation of the complexes is accelerated and enhanced by PEG.

## DIAGNOSTIC IMPLICATIONS

Iron binding protein, binds ferric ions preventing iron intoxication and loss via kidneys. Increased levels are found in iron deficiency, pregnancy, oestrogen administration and lipoidal nephrosis. Decreased levels may be encountered in hereditary deficiencies, testosterone administration, infection, acute inflammation, some forms of nephrosis, tumors, haemochromatosis, acute malaria and malnutrition.

## PERFORMANCE CHARACTERISTICS

### SENSITIVITY

40 mg/dL (Cobas Mira)

### ACCURACY

Controls were measured on the Cobas Mira.

Control	Assigned Value (mg/dL)	Measured Value (mg/dL)
Liquicheck 1 (BIORAD)	142 (114 - 171)	153
Liquicheck 2 (BIORAD)	328 (262 - 394)	354
Seronorm L (NYCOMED)	196 (167 - 225)	207
Seronorm N (NYCOMED)	319 (271 - 367)	328
Seronorm H (NYCOMED)	441 (375 - 507)	406
Immunol. 1 (CIBA CORNING)	146 (124 - 168)	154
Immunol. 2 (CIBA CORNING)	446 (379 - 513)	440

## PRECISION

### Intra-Assay Precision

3 Serum Samples were consecutively measured 20x on the Cobas Mira.

Expected Value	n	Mean	S.D.	C.V
Low	20	97.6	4.66	4.78
Medium	20	252.7	3.29	1.30
High	20	387.1	3.02	0.78

### Inter-Assay Precision

2 Control sera (ORTHO) were measured daily on the Hitachi 717 analyzer after calibration.

Expected Value	n	Mean	S.D.	C.V
Ortho 1	20	174	4.28	2.46
Ortho 2	20	95	4.41	4.64

## METHOD COMPARISON

A comparison with Nephelometry gave the following results:  
 $y = 0.9878x - 0.2018$ ;  $r = 0.9636$

## INTERFERING SUBSTANCES

No interference up to:

Triglycerides	2500 mg/dL	Hemoglobin	1000 mg/dL
Bilirubin	20 mg/dL	Sodium Citrate	1000 mg/dL
Heparin	50 mg/dL		

## QUALITY CONTROL

All commercially available Control sera with Transferrin values measured by this method may be used. We recommend the Dialab Protein Control and the Protein Control Low.

## CALIBRATION

The assay requires the use of Transferrin serum Calibrators. We recommend the Dialab Transferrin Calibrator 5 Level Series, the Protein Calibrator High or the Protein Calibrator Low.

## AUTOMATION

Applications for automated systems (with and without sample dilution) are available upon request.

## WARNINGS AND PRECATIONS

- The Transferrin reagents are intended for in vitro diagnostic use only.
- Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion.
- Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV antibodies, as well as for Hepatitis B surface antigen, using a method approved by the FDA

## WASTE MANAGEMENT

Please refer to local requirements.

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