Immunoglobulin A (IgA)

Diagnostic reagent for the quantitative in vitro determination of IgA (Immunoglobulin A) in human serum by turbidimetric assay

REF

Content

A00521 1x 10 mL IgA Antibody Reagent
 5x 25 mL IgA Buffer

Additionally offered:

A00704 5x 1 mL Protein Calibrator 5 level series
 A00580 1x 1 mL Protein Calibrator High
 A00703 1x 5 mL Protein Calibrator High
 A00701 1x 1 mL Protein Calibrator Low
 A00702 1x 5 mL Protein Calibrator Low
 A00590 1x 1 mL Protein Control
 A00800 1x 5 mL Protein Control
 A08591 1x 1 mL Protein Control Low
 A08823 1x 5 mL Protein Control Low

GENERAL INFORMATION

Method Immunoturbidimetric
Reaction Nonlinear, endpoint
Wavelength 340 nm
Assay Temperature 18 – 37 °C
Sample Serum
Measuring Range approx. 0 – 650 mg/dL
Sensitivity 1 mg/dL (Hitachi 911)
Hook Effect without sample dilution: > 10,000 mg/dL (Hitachi 911)
 with sample dilution: > 10,000 mg/dL (Hitachi 911)
Manual Test Procedure Tests/Kit*
 without sample dilution 125
 with sample dilution 166

Automated Test Procedure

Instrument dependent – please ask for applications

* calculated on amount of antibody reagent; additional buffer on request

REAGENT COMPOSITION

COMPONENTS FINAL CONCENTRATION

IgA Antibody Reagent
 Turbidimetric grade antibody raised in the goat, monospecific for IgA variable
 Sodium azide 0.095 %
IgA Buffer
 Saline 9 g/L
 Enhancer
 Sodium azide 0.095 %

REAGENT PREPARATION

The reagents are liquid and ready to use.

REAGENT STABILITY AND STORAGE

Conditions: Protect from light. Close immediately after use.
 Stability: at 2 – 8 °C up to the expiration date
 at 18 – 25 °C 1 month

Do not freeze!

SAMPLE STABILITY AND STORAGE

Stability: at 2 – 8 °C 48 hours
 at – 20 °C 3 months

Freeze only once!

MANUAL TEST PROCEDURE

Test Procedure without Sample Dilution:

Samples/Controls: ready to use

Calibration curve: Use Protein Calibrator High to generate a calibration curve by making 1:2 serial dilutions of the calibrator with 0.9% saline as diluent or use the 5 level calibrator series. Use 0.9% saline as zero point.

Pipette into test tubes:	Calibrators	Samples/Controls
Buffer	900 µL	900 µL
Cal./Ctrls/Samples	3 µL	3 µL
Mix. Read A1 of calibrators and samples/controls at 340 nm. Then add:		
Antibody Reagent	80 µL	80 µL
Mix. Incubate 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 340 nm. Calculate: $\Delta A = (A2 - A1)$		

Test Procedure with Sample Dilution:

Sample/ Control: dilute 1:10 in saline 0.9%

Calibration curve: Use Protein Calibrator High to generate a calibration curve by making 1:10, 1:20, 1:40, 1:80 and 1:160 dilutions with 0.9% saline as diluent. Use 0.9% saline as zero point.

Pipette into test tubes	Calibrators	Samples/Controls
Buffer	900 µL	900 µL
Cal./Ctrls/Samples	20 µL	20 µL
Mix. Read A1 of calibrators and samples/controls at 340 nm. Then add:		
Antibody Reagent	60 µL	60 µL
Mix. Incubate 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 340 nm. Calculate: $\Delta A = (A2 - A1)$		

CALCULATION

Calculate and plot $\Delta A = (A2 - A1)$ of the calibrators versus assigned concentration values on a linear-linear graph paper. Calculate ΔA optical densities of samples and control(s) and read values in mg/dL on the reference curve. Samples yielding absorbances above highest calibrator should be retested after further dilution.

REFERENCE RANGE

Men: 83 – 406 mg/dL

Women: 70 – 374 mg/dL

It is recommended that each laboratory establishes its own normal range.

TEST PRINCIPLE

The assay of IgA is based on turbidimetric measurement. Turbidity is caused by the formation of antigen-antibody insoluble immuno complexes. The formation of the complexes is accelerated and enhanced by PEG.

DIAGNOSTIC IMPLICATIONS

The measurement of IgA is important for typing immunodeficiencies and myelomas. Furthermore it plays a role in acute and chronic infections as first line of defence. Increased levels may be found in acute infectious hepatitis, chronic aggressive hepatitis, posthepatic/cryptogenic cirrhosis, active alcoholic cirrhosis, chronic infection, rheumatoid arthritis, polydermatomyositis, mixed connective tissue disease.

PERFORMANCE CHARACTERISTICS

SENSITIVITY

1 mg/dL (Hitachi 911)

ACCURACY

Controls were assayed in duplicate on a Hitachi 911.

Control	Assigned Value (mg/dL)	Measured Value (mg/dL)
Cliniqa Level 1	103 (88 – 118)	111
Cliniqa Level 2	204 (174 – 235)	234
Cliniqa Level 3	312 (265 – 359)	323
Biorad Level 1	135 (108 – 162)	140
Biorad Level 2	251 (201 – 301)	245
DIALAB	218 (185 – 251)	221
Behring	170 (144 – 196)	190
CRM	196 (176 – 216)	208

PRECISION

Intra-Assay Precision

3 Control Sera were measured daily on the Hitachi 911.

Expected Value	n	Mean	S.D.	C.V
Low	20	104.42	0.88	0.84
Medium	20	216.85	2.13	0.98
High	20	646.10	7.74	1.20

Inter-Assay Precision

3 Control Sera were measured daily on the Hitachi 911 after calibration.

Expected Value	n	Mean	S.D.	C.V
Low	20	107.05	1.82	1.70
Medium	20	221.65	3.12	1.41
High	20	649.35	6.93	1.07

METHOD COMPARISON

A comparison with Nephelometry gave the following results:
 $y = 1.2169x - 29.831$; $r = 0.9974$

INTERFERING SUBSTANCES

No interference up to:

Triglycerides 2500 mg/dL Hemoglobin 1000 mg/dL
 Bilirubin 20 mg/dL Sodium Citrate 1000 mg/dL
 Heparin 50 mg/dL

QUALITY CONTROL

All commercially available Control sera with IgA values measured by this method may be used. We recommend the Dialab Protein Control and the Protein Control Low.

CALIBRATION

The assay requires the use of IgA serum Calibrators. We recommend the Dialab Protein Calibrator 5 Level Series, the Protein Calibrator High or the Protein Calibrator Low.

AUTOMATION

Applications for automated systems (with and without sample dilution) are available upon request.

WARNINGS AND PRECAUTIONS

- The IgA reagents are intended for in vitro diagnostic use only.
- Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion.
- Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV antibodies, as well as for Hepatitis B surface antigen, using a method approved by the FDA

WASTE MANAGEMENT

Please refer to local requirements.

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