CRP Control Low

Liquid accuracy control for immunoturbidimetric assays for determination of C-reactive Protein (CRP) in serum

REF	Content
A00520	- 1x 1 mL CRP Control Low
A00803	- 1x 5 mL CRP Control Low

For professional in vitro diagnostic use only.

GENERAL INFORMATION

Method	Immunoturbidimetry
Shelf life	36 months from date of production
Storage	at 2 – 8 °C
Other Info	Liquid, ready to use

INTENDED USE

DIALAB CRP Control Low is used as accuracy control for immunoturbidimetric assays with determination of C-reactive Protein (CRP) in serum.

DIAGNOSTIC SIGNIFICANCE

Please refer to the corresponding section of the used CRP reagent.

TEST PRINCIPLE

Please refer to instructions for use of used CRP reagent

COMPOSITION

Dialab CRP Control Low is a dilution of defibrinated human plasma and recombinant human C-reactive protein with phosphate buffered saline. It is stabilised with 0.095 % sodium azide as preservative

MATERIAL REQUIRED BUT NOT PROVIDED

- DIALAB CRP reagent
- DIALAB Calibrator
- Clinical chemistry analyser or photometer

REAGENT PREPARATION

CRP Control Low is ready to use.

STORAGE AND STABILITY

Storage:	at 2 - 8°C
Stability:	up to expiration date
Stability after opening:	6 weeks at 2 - 8°C if stored tightly closed

The control should not be allowed to stand for repeated long periods of time (up to 6 hours) at room temperature

DO NOT FREEZE!

WARNINGS AND PRECAUTIONS

- For in Vitro Diagnostic use.
- Keep vials tightly capped at all time when not in use to avoid any microbial contamination and evaporation.
- Keep vials refrigerated when not in use.
- Each individual donation intended for use in manufacture of CRP Control Low was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C virus (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required tests and only donations with negative findings were used for its manufacture. Nevertheless, every product obtained from human body fluids should be handled with appropriate care in accordance with recommended procedures for biohazardous materials since absence of infectious agents can never be proven.

- Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

TEST PROCEDURE

1. Allow the control to equilibrate to room temperature (22 – 28 °C) prior to use.
2. Use this control in the same manner as a sample as described in the documentation of the used CRP reagent and instrument.

INTERPRETATION OF RESULTS

Results should be within the specified range.

QUALITY CONTROL AND CALIBRATION

For calibration please refer to the instructions for use of the used CRP reagent, calibrator and instrument. Each laboratory should establish its own quality control scheme and corrective actions if controls do not meet the acceptable tolerances. It is recommended to use quality controls after each calibration, at new measuring sessions, after each change of reagents and every time where a potential increased risk of malfunction can be assumed or when results indicate a potential malfunction.

PERFORMANCE CHARACTERISTICS

Please refer to the section "PERFORMANCE CHARACTERISTICS" of the used CRP reagent

TRACEABILITY

The value of CRP Control Low has been assigned to RPPHS/ERM-DA474/IFCC.

EXPECTED VALUES

For the biological reference range please refer to the instruction for use of the used CRP reagent.

LIMITATIONS

For information about biological limitations please refer to the description for use of the used CRP reagent.

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

Please refer to the corresponding section of the used CRP reagent

USED SYMBOLS

Symbol	Description
	Dispose product, residuals and packaging appropriately

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