RF Control

Liquid human based control serum for use as a control in immunoturbidimetric assays for RF (Rheumatoid Factor)

REF	Cont.	
A00810	1x 1 mL	RF Control
A00811	1x 5 mL	RF Control

COMPOSITION

Dialab RF Control is a liquid dilution of human plasma with high values of RF with PBS.

Sodium azide 0.095 %

PREPARATION

The RF Control is liquid, ready to use.

STABILITY AND STORAGE

Storage: at 2 – 8 °C
 Stability: up to the expiration date
 Stability after opening: 6 weeks at 2 – 8 °C
DO NOT FREEZE!

DIRECTIONS FOR USE

1. Use this control in the same manner as specified with the reagent and/or instrument being used.
2. Allow the control to equilibrate to room temperature (22 – 28 °C) prior to use.
3. The control should not be allowed to stand for repeated long periods of time (up to 6 hours) at room temperature. Keep vials tightly capped at all times when not in use to avoid microbial contamination. Keep vials refrigerated when not in use

WARNINGS AND PRECAUTIONS

1. For in Vitro Diagnostic use.
2. Each individual donation intended for use in manufacture of protein calibrator was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required tests. Since no test method can assure that products derived from human blood do not contain HIV-1/2 and Hepatitis B and Hepatitis C virus, this material and all patient samples should be handled as though capable of transmitting infectious diseases.
3. Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

WASTE MANAGEMENT

Please refer to local legal requirements.

LOT SPECIFIC ASSAY DATA

Values and expiry date are lot specific.

LOT: C1043		EXP: 2017/04
Parameter	Target Value	Range
RF	105 IU/mL	89.0 – 121 IU/mL
Aggregated human IgG Method + Latex Method, based on WHO Standardisation		

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