## Microalbumin Calibrator 5 Level Series

Liquid human based calibrators for use as calibrators  
in immunoturbidimetric assays for Microalbumin

**REF**

**A00955**

**Cont.**

5x 1 mL

Microalbumin Calibrator 5 Level Series

### COMPOSITION

Dialab Microalbumin Calibrators are dilutions of defibrinated human plasma with phosphate buffered saline, liquid stabilized and filtered through 0.2 µ.

Sodium azide 0.095 %

### PREPARATION

The Microalbumin Calibrator Series is liquid, ready to use.

### STABILITY AND STORAGE

Storage: at 2 – 8 °C  
 Stability: up to the expiration date  
 Stability after opening: 6 weeks at 2 – 8 °C  
 DO NOT FREEZE!

### DIRECTIONS FOR USE

- Use these Calibrators in the same manner as specified with the reagent and/or instrument being used.
- Allow the Calibrators to equilibrate to room temperature (22 – 28 °C) prior to use.
- The Calibrators should not be allowed to stand for repeated long periods of time (up to 6 hours) at room temperature. Keep vials tightly capped at all times when not in use to avoid microbial contamination. Keep vials refrigerated when not in use

### WARNINGS AND PRECAUTIONS

- For in Vitro Diagnostic use.
- Each individual donation intended for use in manufacture of protein calibrator was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required tests. Since no test method can assure that products derived from human blood do not contain HIV-1/2 and Hepatitis B and Hepatitis C virus, this material and all patient samples should be handled as though capable of transmitting infectious diseases.
- Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LOT SPECIFIC ASSAY DATA

Values and expiry date are lot specific.

LOT: MAS017		EXP: 2017/04
Parameter	Level	Assigned Value
Microalbumin	Level 1	21.0 mg/L
Microalbumin	Level 2	50.0 mg/L
Microalbumin	Level 3	102.0 mg/L
Microalbumin	Level 4	200.0 mg/L
Microalbumin	Level 5	400.0 mg/L
Values are traceable to RPPHS/ERM-DA470 from IFCC.		

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