## Transferrin Calibrator 5 Level Series

Liquid human based calibrators for use as calibrators  
in immunoturbidimetric assays for Transferrin

**REF**

**A00960**

**Cont.**

5x 1 mL

Transferrin Calibrator 5 Level Series

### COMPOSITION

Dialab Transferrin calibrators are of defibrinated and delipitated human plasma, liquid stabilised and filtered through a 0.2 µ pore-size filter.  
 Sodium azide 0.095 %

### PREPARATION

The Transferrin Calibrator Series is liquid, ready to use.

### STABILITY AND STORAGE

Storage: at 2 – 8 °C  
 Stability: up to the expiration date  
 Stability after opening: 6 weeks at 2 – 8 °C  
**DO NOT FREEZE!**

### DIRECTIONS FOR USE

- Use these Calibrators in the same manner as specified with the reagent and/or instrument being used.
- Allow the Calibrators to equilibrate to room temperature (22 – 28 °C) prior to use.
- The Calibrators should not be allowed to stand for repeated long periods of time (up to 6 hours) at room temperature. Keep vials tightly capped at all times when not in use to avoid microbial contamination. Keep vials refrigerated when not in use

### WARNINGS AND PRECAUTIONS

- For in Vitro Diagnostic use.
- Each individual donation intended for use in manufacture of protein calibrator was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required tests. Since no test method can assure that products derived from human blood do not contain HIV-1/2 and Hepatitis B and Hepatitis C virus, this material and all patient samples should be handled as though capable of transmitting infectious diseases.
- Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LOT SPECIFIC ASSAY DATA

Values and expiry date are lot specific.

LOT: 0616		EXP: 2017/09
Parameter	Level	Assigned Value
Transferrin	Level 1	102.0 mg/dL
Transferrin	Level 2	198.0 mg/dL
Transferrin	Level 3	238.0 mg/dL
Transferrin	Level 4	392.0 mg/dL
Transferrin	Level 5	484.0 mg/dL
Values are traceable to RPPHS/ERM-DA470k/IFCC.		

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