## Pediatric Control

Liquid human based control serum for use as a control in immunoturbidimetric assays

REF	Cont.	
A03808	1x 1 mL	Pediatric Control
A03809	1x 5 mL	Pediatric Control

### COMPOSITION

Dilution of delipitated and defibrinated normal human plasma with PBS.  
 Sodium azide 0.095 %

### PREPARATION

The Control is liquid, ready to use.

### STABILITY AND STORAGE

Storage: at 2 – 8 °C  
 Stability: up to the expiration date  
 Stability after opening: 6 weeks at 2 – 8 °C  
 DO NOT FREEZE!

### DIRECTIONS FOR USE

1. Use this Calibrator in the same manner as specified with the reagent and/or instrument being used.
2. Allow the Calibrator to equilibrate to room temperature (22 – 28 °C) prior to use.
3. The Calibrator should not be allowed to stand for repeated long periods of time (up to 6 hours) at room temperature. Keep vials tightly capped at all times when not in use to avoid microbial contamination. Keep vials refrigerated when not in use

### WARNINGS AND PRECATIONS

1. For in Vitro Diagnostic use.
2. Each individual donation intended for use in manufacture of protein calibrator was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required tests. Since no test method can assure that products derived from human blood do not contain HIV-1/2 and Hepatitis B and Hepatitis C virus, this material and all patient samples should be handled as though capable of transmitting infectious diseases.
3. Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LOT SPECIFIC ASSAY DATA

Values and expiry date are lot specific.

LOT: PED009		EXP: 2016/01
Parameter	Target Value	Range
α-2 Macroglobulin	10.3 mg/dL	8.8 – 11.8 mg/dL
Transferrin	13.0 mg/dL	11.1 – 15.0 mg/dL
IgA	9.9 mg/dL	8.4 – 11.4 mg/dL
IgG	44.9 mg/dL	38.2 – 51.6 mg/dL
IgM	6.0 mg/dL	5.1 – 6.9 mg/dL
Kappa Light Chain **	12.0 mg/dL	10.2 – 13.8 mg/dL
Lambda Light Chain **	5.8 mg/dL	4.9 – 6.7 mg/dL

Values traceable to ERM-DA470k/IFCC.  
 \*\* Values based on Dade Behring Material.

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