

**Liquid Reagents – ready to use**

# GAMMA GT

**(γ-Glutamyl Transferase)**

**SZASZ, standardized to IFCC**

2 Reagents

**Diagnostic reagent for quantitative in vitro determination of gamma-glutamyltransferase (gamma-GT) in human serum or plasma on photometric systems**

REF	Cont.
-----	-------

<b>D95604</b>	<b>5 x 50 ml</b>	4 x 50 ml	Reagent 1
		1 x 50 ml	Reagent 2

Additionally offered:

D98485	5 x 3 ml	Calibrator	Diacal Auto
D98481	12 x 5 ml	Control normal	Diacon N
D98482	12 x 5 ml	Control abnormal	Diacon P

**TEST PARAMETERS**

Method:	Colorimetric, Kinetic, Increasing Reaction SZASZ, standardized to IFCC
Wavelength:	405 nm (400 – 420 nm)
Temperature:	37°C
Sample:	Serum, heparin plasma
Linearity:	up to 284 U/L (Szasz, Substrate start) up to 232 U/L (Szasz, Sample start) Up to 321 U/L (IFCC, Substrate start) Up to 262 U/L (IFCC, Sample start)
Sensitivity:	The lower limit of detection is 2 U/L

**REAGENT COMPOSITION**

COMPONENTS	FINAL CONCENTRATION
<b>Reagent 1:</b>	
Glycylglycine	100 mmol/L
Tris, pH 8.25	100 mmol/L
<b>Reagent 2:</b>	
L-Gamma-Glutamyl-3-Carboxy-4-Nitroanilide	4 mmol/L

**REAGENT PREPARATION**

**Substrate Start:**  
Reagents are ready for use.

**Sample Start:**  
Mix 4 parts of Reagent 1 with 1 part of Reagent 2.  
( = Working Reagent)

**REAGENT STABILITY AND STORAGE**

Conditions: protect from light  
close immediately after use  
do not freeze!

**Substrate Start:**  
Storage: at 2 – 8°C  
Stability: up to the expiration date

**Sample Start (Working Reagent):**  
Storage: at 2-8°C at 15-25°C  
Stability: 4 weeks 5 days

Maximum allowable absorbance of the Working Reagent measured at 405 nm against water as reference is 1.3.

**SAMPLE STABILITY AND STORAGE**

Stability: at -20 - +25°C at least 1 week

**INTERFERING SUBSTANCES**

no interference up to:

ascorbic acid	30 mg/dl
bilirubin	40 mg/dl
triglycerides	2000 mg/dl
hemoglobin	400 mg/dl

**MANUAL TEST PROCEDURE**

Bring reagents and samples to room temperature.

**Substrate start 37°C**

Pipette into test tubes:	Blank	Sample
Reagent 1	1000 µl	1000 µl
Distilled water	100 µl	
Sample		100 µl
Mix. Incubate for approximately 1 minute. Then add:		
Reagent 2	250 µl	250 µl
Mix. Read initial absorbance after 1 minute and start a timer. Read absorbance again after exactly 1, 2 and 3 minutes. Determine $\Delta A/\text{min}$ . during the linear part of the assay.		

**Sample start 37°C**

Pipette into test tubes	Blank	Sample
Working reagent	1000 µl	1000 µl
Distilled water	100 µl	
Sample		100 µl
Mix. Read initial absorbance after 1 minute and start a timer. Read absorbance again after exactly 1, 2 and 3 minutes. Determine $\Delta A/\text{min}$ . during the linear part of the assay.		

**CALCULATION** (light path 1 cm)

$$\Delta A/\text{min} = [\Delta A/\text{min sample}] - [\Delta A/\text{min blank}]$$

$$\text{Gamma-GT (U/L)} = \Delta A/\text{min} \times \text{Factor}$$

**Factors (25/30/37°C):**

	Szasz	IFCC
<b>For Substrate Start:</b>	1421	1606
<b>For Sample start:</b>	1158	1309

**UNIT CONVERSION**

U/L x 0.01667 =  $\mu$ katal/L

## REFERENCE RANGE \*(U/L)

According to Szasz:

<b>Women</b>	< 32
<b>Men</b>	< 49

According to IFCC:

	Female	Male
Adults	< 38 U/L	< 55 U/L
1 day – 6 months	15 – 132 U/L	12 – 122 U/L
6 months – 1 year	1 – 39 U/L	1 – 39 U/L
1 – 12 years	4 – 22 U/L	3 – 22 U/L
13 – 18 years	4 – 24 U/L	2 – 42 U/L

\* It is recommended that each laboratory establishes its own normal range.

## TEST PRINCIPLE

Gamma-GT catalyzes the transfer of glutamic acid to acceptors like glycylglycine. This process releases 5-amino-2-nitrobenzoate, which can be measured at 405 nm. The increase in absorbance at this wavelength is directly related to the activity of gamma-GT.

L- $\gamma$ -glutamyl-3-carboxy-4-nitroanilide+glycylglycine  $\xrightarrow{\gamma\text{-GT}}$   
L- $\gamma$ -glutamyl-glycylglycine + 5-amino-2-nitrobenzoate

## PERFORMANCE CHARACTERISTICS

### LINEARITY

The assay is linear up to a  $\Delta A/\text{min} = 0.20$ . That is 284 U/L (Substr. Start) and 232 U/L (Sample start). Above this concentration, dilute the sample 1+5 with NaCl solution (9 g/L sodium chloride in water) and reassay multiplying the result by 6.

### PRECISION (at 37°C)

Intra-assay n = 20	Mean [U/L]	SD [U/L]	CV [%]
Sample 1	39.9	0.99	2.48
Sample 2	73.6	0.85	1.16
Sample 3	206	1.32	0.64

Inter-assay n = 20	Mean [U/L]	SD [U/L]	CV [%]
Sample 1	41.5	0.62	1.49
Sample 2	72.3	0.61	0.85
Sample 3	205	0.74	0.36

## METHOD COMPARISON

A comparison between Dialab Gamma-GT, standardized to IFCC (y) and the IFCC reference reagent (x) using 51 samples gave following results:  
 $y = 1.005 x - 0.741 \text{ U/L}; r = 0.999$

A comparison between Dialab Gamma-GT, Szasz (y) and a commercially available test (x) using 51 samples gave following results:  $y = 0.996 x + 1.354 \text{ U/L}; r = 1.000$ .

## QUALITY CONTROL

All control sera with Gamma GT values determined by this method can be used.

We recommend:

REF	Cont.		
<b>D98481</b>	12 x 5 ml	<b>DIACON N</b>	Assayed Control Serum Normal
<b>D98482</b>	12 x 5 ml	<b>DIACON P</b>	Assayed Control Serum Abnormal

## CALIBRATION

The use of a Gamma GT Calibrator is optional.

We recommend:

REF	Cont.		
<b>D98485</b>	5 x 3 ml	<b>DIACAL AUTO</b>	Assayed Multi Calibration Serum

## AUTOMATION

Special adaptations for automated analyzers can be made on request.

## WARNINGS AND PRECAUTIONS

1. The reagents contain sodium azide (0.95 g/L) as preservative. Do not swallow! Avoid contact with skin and mucous membranes.
2. Take the necessary precautions for the use of laboratory reagents.

## WASTE MANAGEMENT

Please refer to local legal requirements.

## REFERENCES

1. Thomas L. Clinical Laboratory Diagnostics. 1<sup>st</sup> ed. Frankfurt: TH-Books Verlagsgesellschaft;1998.p.80-6.
2. Persijn JP, van der Silk W. A new method for the determination of gamma-glutamyltransferase in serum. J Clin Chem Clin Biochem 1976; 14:421-7.
3. Szasz G. Gamma-Glutamyltranspeptidase. In: Bergmeyer HU. Methoden der enzymatischen Analyse. Weinheim: Verlag Chemie, 1974. p. 757.
4. Schumann G, Bonora R, Ceriotti F, Féraud G et al. IFCC primary reference procedure for the measurement of catalytic activity concentrations of enzymes at 37°C. Part 5: Reference procedure for the measurement of catalytic concentration of gamma-glutamyltransferase. Clin Chem Lab Med 2002; 40:734-8
5. Fischbach F, Zawta B. Age-dependent reference limits of several enzymes in plasma at different measuring temperatures. Klin Lab 1992;38:555-61.



DIALAB Produktion und Vertrieb von chemisch – technischen Produkten und Laborinstrumenten Gesellschaft m.b.H.  
A – 2351 Wiener Neudorf, Austria  
IZ-NÖ Süd, Hondastrasse, Objekt M55  
Phone: ++43 (0) 2236 660910-0  
Fax: ++43 (0) 2236 660910-30 e-mail: [office@dialab.at](mailto:office@dialab.at)