CRP Uhs Control Set for the universal range

Control Set for use as controls in tests for quantitative in vitro determination of CRP on photometric systems

REF	Cont.		
D99557	2x 2 mL	1x 2 mL Level 1 1x 2 mL Level 2	CRP Uhs Control Set (for the universal range)

COMPOSITION

The CRP Control Set is a set of 2 liquid controls based on human serum intended for the use as controls for CRP in the universal range.

PREPARATION

The Controls are liquid and ready to use.

STABILITY AND STORAGE

Storage: at 2 – 8 °C
 Stability: up to the expiration date
 Stability after opening: at least 2 weeks at 2 – 8 °C
DO NOT FREEZE!

DIRECTIONS FOR USE

- Use these controls in the same manner as specified with the reagent and/or instrument being used.
- Allow the controls to equilibrate to room temperature (22 – 28 °C) prior to use.
- The controls should not be allowed to stand for repeated long periods of time (up to 6 hours) at room temperature. Keep vials tightly capped at all times when not in use to avoid microbial contamination. Keep vials refrigerated when not in use

WARNINGS AND PRECATIONS

- For in Vitro Diagnostic use.
- Each individual donation intended for use in manufacture of protein calibrator was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required tests. Since no test method can assure that products derived from human blood do not contain HIV-1/2 and Hepatitis B and Hepatitis C virus, this material and all patient samples should be handled as though capable of transmitting infectious diseases.
- Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

WASTE MANAGEMENT

Please refer to local legal requirements.


LITERATURE

- Dati F. Reference Materials and guidelines for standardisation of methods in laboratory medicine. In: Thomas L., editor. Clinical laboratory diagnostics. 1st ed. Frankfurt: TH Books Verlagsgesellschaft; 1998. p. 1393-1401
- Röhle G., Siekmann L. Quality Assurance of quantitative determination. In: Thomas L., editor. Clinical laboratory diagnostics. 1st ed. Frankfurt: TH Books Verlagsgesellschaft; 1998. p. 1393-1401
- Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395.

LOT SPECIFIC ASSAY DATA

Values and expiry date are lot specific.

LOT: 16722		EXP: 2013/10	
Level	LOT	Target Value	Range
Level 1	16722	24.80 mg/L	18.10 – 31.50 mg/L
Level 2	16723	62.40 mg/L	45.60 – 79.20 mg/L

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