Instructions for Use

β -HCG ELISA

IVD

REF EIA-1911

96

DRG Instruments GmbH, Germany
Frauenbergstraße 18, 35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

Distributed by:

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

**Please use only the valid version of the Instructions for Use provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
Por favor, use sólo la versión válida de las instrucciones de uso que se suministran con el kit.
Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.**

Table of Contents / Inhaltsverzeichnis / Tabella die Contenuti / Tabla de Contenidos / Table des matières

1	INTENDED USE.....	2	1	FINALIDAD PREVISTA.....	24
2	PRINCIPLE OF THE TEST	2	2	FUNDAMENTO DEL ENSAYO	24
3	WARNINGS AND PRECAUTIONS	3	3	PRECAUCIONES	24
4	REAGENTS.....	4	4	COMPONENTES DEL KIT.....	25
5	SPECIMEN COLLECTION AND PREPARATION.....	5	5	MUESTRAS	26
6	ASSAY PROCEDURE.....	5	6	PROCEDIMIENTO DE ENSAYO	26
7	EXPECTED NORMAL VALUES.....	7	7	VALORES ESPERADOS	28
8	QUALITY CONTROL.....	7	8	CONTROL DE CALIDAD	28
9	PERFORMANCE CHARACTERISTICS.....	8	9	CARACTERÍSTICAS DEL ENSAYO.....	29
10	LIMITATIONS OF USE.....	9	10	LIMITACIONES DE USO	29
11	LEGAL ASPECTS	9	11	ASPECTOS LEGALES	30
12	REFERENCES / LITERATURE.....	9	12	REFERENCIAS / BIBLIOGRAFÍA.....	30
1	ZWECKBESTIMMUNG	10	1	UTILISATION PRÉVUE	31
2	TESTPRINZIP	10	2	PRINCIPE DU TEST	31
3	VORSICHTSMAßNAHMEN.....	11	3	PRECAUTIONS D'UTILISATION	31
4	BESTANDTEILE DES KITS	12	4	RÉACTIFS	32
5	PROBENVORBEREITUNG.....	13	5	ECHANTILLON	33
6	TESTDURCHFÜHRUNG	13	6	RÉALISATION DU TEST	33
7	ERWARTETE WERTE	15	7	VALEURS ATTENDUES.....	35
8	QUALITÄTSKONTROLLE.....	15	8	CONTROLE DE QUALITE	35
9	ASSAY-CHARAKTERISTIKA.....	16	9	CARACTERISTIQUES DU TEST	36
10	GRENZEN DES TESTS	16	10	LIMITES D'UTILISATION.....	36
11	RECHTLICHE GRUNDLAGEN	17	11	ASPECTS LEGAUX.....	37
12	REFERENZEN / LITERATUR	17	12	REFERENCES / BIBLIOGRAPHIE	37
1	DESTINAZIONE D'USO	18		SYMBOLS USED.....	38
2	PRINCIPIO DEL TEST	18			
3	PRECAUZIONI.....	18			
4	COMPONENTI DEL KIT.....	19			
5	CAMPIONI.....	20			
6	ATTUAZIONE DEL TEST.....	20			
7	VALORI NORMALI.....	22			
8	CONTROLLO QUALITÀ.....	22			
9	CARATTERISTICHE DEL TEST	22			
10	LIMITAZIONE DEL TEST	23			
11	ASPETTI LEGALI	23			
12	BIBLIOGRAFIA.....	23			

1 INTENDED USE

The **DRG β -HCG ELISA Kit** is an enzyme immunoassay for the quantitative *in vitro* diagnostic measurement of total human chorionic gonadotropin (hCG and β -hCG) in serum.

This kit is NOT intended to be used for the risk evaluation of trisomy 21.

1.1 Summary and Explanation

Chorionic Gonadotropin (hCG) is a glycoprotein hormone which is normally produced by the placenta during pregnancy. After conception, the hCG concentration increases rapidly to reach a peak near the end of the first trimester. High concentrations are observed throughout pregnancy. After delivery, hCG levels fall rapidly and become undetectable after a few days.

Structurally intact hCG molecules are composed of an alpha and a beta subunit. The alpha subunit is nearly identical to the alpha subunits of other glycoprotein hormones, such as Thyroid Stimulating Hormone (TSH), Luteinizing Hormone (LH), and Follicle Stimulating Hormone (FSH): The differences in the beta subunit of the respective hormones account for their biological specificity and immunochemical distinctiveness.

Monoclonal antibodies recognizing unique sites on the beta chain of the β -hCG/hCG molecule are essential for differentiation between hCG and LH, FSH and TSH.

Specific assays for beta-hCG permit the early detection of pregnancy.

In addition to the elevated hCG levels during pregnancy, high concentrations of β hCG/hCG may be associated with neoplasms of trophoblastic and nontrophoblastic origin such as hydatiform mole, chorionepithelioma, embryonal cell carcinoma, seminoma and many others.

2 PRINCIPLE OF THE TEST

The DRG β -HCG ELISA is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the **sandwich principle**.

The microtiter wells are coated with a monoclonal antibody [mouse] directed towards a unique antigenic site on a β -HCG molecule. An aliquot of patient sample containing endogenous β -HCG and/or HCG is incubated in the coated well with enzyme conjugate, which is an anti- β -HCG antibody conjugated with horseradish peroxidase. After incubation the unbound conjugate is washed off.

The amount of bound peroxidase is proportional to the concentration of β -HCG/HCG in the sample.

Having added the substrate solution, the intensity of colour developed is proportional to the concentration of β -HCG/HCG in the patient sample.

3 WARNINGS AND PRECAUTIONS

1. This kit is for *in vitro* diagnostic use only. For professional use only.
2. All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
3. Before starting the assay, read the instructions completely and carefully. Use the valid version of instructions for use provided with the kit. Be sure that everything is understood.
4. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
5. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
6. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
7. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
8. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
9. Allow the reagents to reach room temperature (21 °C to 26 °C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the patient samples will not be affected.
10. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
11. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
12. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
13. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
14. Do not use reagents beyond expiry date as shown on the kit labels.
15. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
16. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
17. Avoid contact with Stop Solution containing 0.5 M H₂SO₄. It may cause skin irritation and burns.
18. Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
19. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
20. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
21. For information on hazardous substances included in the kit please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from DRG.

4 REAGENTS

4.1 Reagents provided

1. **Microtiterwells**, 12 x 8 (break apart) strips, 96 wells;
Wells coated with anti-β-HCG antibody (monoclonal).
2. **Standard (Standard 0-5)**, 6 vials (lyophilized), 1 mL;
Concentrations: 0, 5; 25; 50; 100; 200 mIU/mL
Conversion: 1 pg/mL = 0.00916 mIU/mL
The standards are calibrated against the following reference material: 5th WHO International Standard Chorionic Gonadotrophin, NIBSC code: 07/364.
See „Preparation of Reagents“;
Contain non-mercury preservative.
3. **Sample Diluent**, 1 vial, 10 mL, ready to use,
Contains non-mercury preservative.
4. **Enzyme Conjugate**, 1 vial, 11 mL, ready to use,
Anti-β-HCG antibody conjugated to horseradish peroxidase;
Contains non-mercury preservative.
5. **Substrate Solution**, 1 vial, 14 mL, ready to use,
Tetramethylbenzidine (TMB).
6. **Stop Solution**, 1 vial, 14 mL, ready to use,
contains 0.5 M H₂SO₄,
Avoid contact with the stop solution. It may cause skin irritations and burns.

Note: Additional *Sample Diluent* for sample dilution is available upon request.

4.2 Materials required but not provided

- A calibrated microtiter plate reader (450 nm, with reference wavelength at 620 nm to 630 nm)
- Calibrated variable precision micropipettes
- Absorbent paper
- Distilled water
- Timer
- Graph paper or software for data reduction

4.3 Storage Conditions

When stored at 2 °C - 8 °C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date.

Opened reagents must be stored at 2 °C - 8 °C. Microtiter wells must be stored at 2 °C - 8 °C. Once the foil bag has been opened, care should be taken to close it tightly again.

Opened kits retain activity for 8 weeks if stored as described above.

4.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature prior to use.

Standards

Reconstitute the lyophilized contents of each vial with 1 ml distilled water and let stand for at least 10 minutes at room temperature. Mix several times before use.

Note: *The reconstituted standards are stable for 2 months at 2 °C - 8 °C.
For longer storage aliquot and freeze at -20 °C.*

4.5 Disposal of the Kit

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Safety Data Sheet, section 13.

4.6 Damaged Test Kits

In case of any severe damage to the test kit or components, DRG has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

5 SPECIMEN COLLECTION AND PREPARATION

Serum should be used in this assay.

Note: Samples containing sodium azide should not be used in the assay.

In general, it should be avoided to use haemolytic, icteric, or lipaemic specimens. For further information refer to chapter “*Interfering Substances*”.

5.1 Specimen Collection

Serum:

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

5.2 Specimen Storage and Preparation

Specimens should be capped and may be stored for up to 7 days at 2 °C - 8 °C prior to assaying.

Specimens held for a longer time (up to 12 months) should be frozen only once at -20 °C prior to assay. Thawed samples should be inverted several times prior to testing.

5.3 Specimen Dilution

If in an initial assay, a specimen is found to contain more than the highest standard, the specimens can be diluted with *Sample Diluent* and reassayed as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

- a) dilution 1:10: 10 µL Serum + 90 µL *Sample Diluent* (mix thoroughly)
- b) dilution 1:100: 10 µL dilution a) 1:10 + 90 µL *Sample Diluent* (mix thoroughly).
- c) dilution 1:1000 10 µL dilution b) 1:100 + 90 µL *Sample Diluent* (mix thoroughly).

NOTE:

Sera of pregnant women must be diluted 1/100 in *Sample Diluent* before starting the assay.

6 ASSAY PROCEDURE

6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination.
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.

6.2 Test Procedure (quantitative method)

Each run must include a standard curve.

1. Secure the desired number of Microtiter wells in the holder.
2. Dispense **25 µL** of each **Standard, Control** and **samples** with new disposable tips into appropriate wells.
3. Dispense **100 µL Enzyme Conjugate** into each well.
Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
4. Incubate for **60 minutes** at room temperature.
5. Briskly shake out the contents of the wells. Rinse the wells 5 times with distilled water (400 µL per well). Strike the wells sharply on absorbent paper to remove residual droplets.
Important note: The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
6. Add **100 µL of Substrate Solution** to each well.
7. Incubate for **15 minutes** at room temperature.
8. Stop the enzymatic reaction by adding **50 µL of Stop Solution** to each well.
9. Determine the optical density of the solution in each well at **450 nm (reading) and at 620 - 630 nm (background subtraction, recommended)** with a microtiter plate reader.
It is recommended that the wells be read **within 10 minutes** after adding the *Stop Solution*.

6.3 Calculation of Results (quantitative)

1. Calculate the average optical density (OD) values for each set of standards, controls and patient samples.
2. Using linear graph paper, construct a standard curve by plotting the mean OD obtained from each standard against its concentration with OD value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean OD value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the Instructions for Use have been calculated automatically using a 4-Parameter curve fit. (4 Parameter Rodbard or 4 Parameter Marquardt are the preferred methods.) Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this standard curve. Samples with concentrations higher than that of the highest standard have to be further diluted or reported as > 200 mIU/mL. For the calculation of the concentrations this dilution factor has to be taken into account.

6.3.1 Example of Typical Standard Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Standard	Optical Density (450 nm)
Standard 0 (0 mIU/mL)	0.04
Standard 1 (5 mIU/mL)	0.15
Standard 2 (25 mIU/mL)	0.28
Standard 3 (50 mIU/mL)	0.53
Standard 4 (100 mIU/mL)	0.94
Standard 5 (200 mIU/mL)	1.50

6.4 Assay Procedure (qualitative method)

This procedure differentiates positive (pregnant) from negative samples by comparing the sample beta hCG levels with the *Standard 0* (0 mIU/mL) and *Standard 3* (50 mIU/mL).

Patient samples are run with the Standard 0 and the 50 mIU/mL Standard. The assay procedure is identical with the Quantitative Method, but step 9 and 10 is omitted.

6.5 Calculation of Results (qualitative)

For a qualitative analysis of the β-hCG level the color development of the specimen is compared with the color of the 0 mIU/mL and 50 mIU/mL standards.

If the blue color is less intense than the color of the 50 mIU/mL standard, the sample is considered as negative.

If the blue color is more intense than or equal to the color of the 50 mIU/mL standard the sample is considered as positive.

7 EXPECTED NORMAL VALUES

It is strongly recommended that each laboratory should determine its own normal and abnormal values.

7.1 Normal healthy adults, non-pregnant

In a study conducted with apparently normal healthy adults, using the DRG β-HCG ELISA the following values are observed:

Population	Age (years)	Valid N	β-hCG [mIU/mL]
men	< 50	40	< 5
	> 50	10	< 5
women	< 50	42	< 5
	> 50	7	< 5

7.2 Pregnant women in the 2nd trimester

Week of pregnancy	Valid samples	5% Percentile [mIU/mL]	95% Percentile [mIU/mL]
14	103	10 303	71 980
15	96	9 246	51 666
16	91	5 266	36 947
17	14	4 632	24 033

CAUTION:

- For the detection of pregnancy in serum, a qualitative assay is used with a cut-off point of 50 mIU/mL. Negative or borderline results should be repeated on a fresh specimen obtained at least 48 hours after the first specimen.
- It has been shown that immunological pregnancy tests may yield false results in cases of several medications and diseases such as rheumatoid arthritis. In such cases, the interpretation of the pregnancy test should be done carefully.

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated to other clinical observations and diagnostic tests.

8 QUALITY CONTROL

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or DRG directly.

9 PERFORMANCE CHARACTERISTICS

9.1 Assay Dynamic Range

The range of the assay is between 0.44 – 200 mIU/mL.

9.2 Specificity of Antibodies (Cross-Reactivity)

The following substances were tested for cross reactivity of the assay:

Protein	Concentration	Produced Colour Intensity Equivalent to β-HCG in serum (mIU/mL)
hLH	300 mIU/mL	6
hLH	200 mIU/mL	< 5
hLH	80 mIU/mL	< 5
TSH	75 µIU/mL	8
TSH	50 µIU/mL	< 5
TSH	25 µIU/mL	< 5
FSH	200 mIU/mL	< 5
FSH	50 mIU/mL	< 5

9.3 Sensitivity

The analytical sensitivity was calculated from the mean plus two standard deviations of twenty (20) replicate analyses of the Standard 0 and was found to be 0.44 mIU/mL.

9.4 Reproducibility

9.4.1 Intra Assay Variation

The within assay variability is shown below:

Sample	n	Mean (mIU/mL)	CV (%)
1	20	5.9	9.4
2	20	18.6	4.0
3	20	147.9	3.5

9.4.2 Inter Assay Variation

The between assay variability is shown below:

Sample	Mean (mIU/mL)	CV (%)
1	5.6	9.9
2	17.1	7.2
3	140.0	6.9

9.5 Recovery

Samples have been spiked by adding β-HCG solutions with known concentrations.

The recovery (%) was calculated by multiplying the ratio of measured and expected values with 100.

	Sample 1	Sample 2	Sample 3
Concentration (mIU/mL)	30.9	69.0	115.2
Average Recovery (%)	94.4	91.3	90.7
Range of Recovery (%)	from	91.3	88.6
	to	99.3	94.2

9.6 Linearity

Samples were measured undiluted and in serial dilutions with standard 0. The recovery (%) was calculated by multiplying the ratio of expected and measured values with 100.

		Sample 1	Sample 2	Sample 3
Concentration (mIU/mL)		39.7	75.6	128.4
Average Recovery (%)		102.0	95.5	96.8
Range of Recovery (%)	from	90.9	86.2	92.9
	to	115.0	109.3	101.3

10 LIMITATIONS OF USE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice.

Any improper handling of samples or modification of this test might influence the results.

10.1 Interfering Substances

Haemoglobin (up to 4 mg/mL), Bilirubin (up to 0.5 mg/mL) and Triglyceride (up to 30 mg/mL) have no influence on the assay results.

10.2 Drug Interferences

Until today no substances (drugs) are known to us, which have an influence to the measurement of β -HCG in a sample.

10.3 High-Dose-Hook Effect

No hook effect was observed in this test up to 15 800 mIU/mL of β -HCG.

11 LEGAL ASPECTS

11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2 are also invalid.

Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

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