## Instructions for Use

# MetCombi Plasma ELISA

## (Metanephrine- Normetanephrine)

*Enzyme Immunoassay for the quantitative determination of free Metanephrine and free Normetanephrine in plasma*

**EIA-4082**

**2 x 96 wells**

**DRG** 

DRG Instruments GmbH, Germany  
Frauenbergstr. 18, D-35039 Marburg  
Telefon: +49 (0)6421-1700 0, Fax: +49-(0)6421-1700 50  
Internet: [www.drg-diagnostics.de](http://www.drg-diagnostics.de)  
E-Mail: [drg@drg-diagnostics.de](mailto:drg@drg-diagnostics.de)

**DRG** 

DRG International, Inc.  
USA

Telephone: (908) 233-2079  
Fax: (908) 233-0758  
E-Mail: [corp@drg-international.com](mailto:corp@drg-international.com)

**Please use only the valid version of the package insert provided with the kit.  
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Arbeitsanleitung.  
Si prega di usare la versione valida dell'inserto del pacco a disposizione con il kit.  
Por favor, se usa solo la version valida de la metodico técnico incluido aqui en el kit.**

**Table of Contents / Inhaltsverzeichnis / Tabella die Contenuti / Tabla de Contenidos**

1	INTRODUCTION.....	2
2	PROCEDURAL CAUTIONS, GUIDELINES, WARNINGS AND LIMITATIONS .....	2
3	STORAGE AND STABILITY .....	3
4	MATERIALS .....	4
5	SAMPLE COLLECTION AND STORAGE .....	5
6	TEST PROCEDURE .....	6
7	CALCULATION OF RESULTS .....	8
8	ASSAY CHARACTERISTICS .....	9
9	REFERENCES/LITERATURE .....	9
1	EINLEITUNG.....	10
2	VERFAHRENSHINWEISE, RICHTLINIEN, WARNUNGEN UND ANWENDUNGSGRENZEN.....	10
3	LAGERUNG UND HALTBARKEIT.....	11
4	MATERIALIEN .....	12
5	PROBENMATERIAL UND LAGERUNG .....	13
6	TESTDURCHFÜHRUNG .....	14
7	BERECHNUNG DER ERGEBNISSE.....	16
8	TESTCHARAKTERISTIKA .....	17
9	REFERENZEN/LITERATUR.....	17
	USED SYMBOLS.....	18

## 1 INTRODUCTION

### 1.1 Intended use and principle of the test

Enzyme Immunoassay for the quantitative determination of free Metanephrine and free Normetanephrine in plasma.

First, the plasma proteins are removed by precipitation. After this Metanephrine (Metadrenaline) and Normetanephrine (Normetadrenaline) are quantitatively acylated.

The subsequent competitive ELISA uses the microtiter plate format. The antigen is bound to the solid phase of the microtiter plate. The acylated standards, controls and samples and the solid phase bound analytes compete for a fixed number of antibody binding sites. After the system is in equilibrium, free antigen and free antigen-antibody complexes are removed by washing. The antibody bound to the solid phase is detected by an anti-rabbit IgG-peroxidase conjugate using TMB as a substrate. The reaction is monitored at 450 nm.

Quantification of unknown samples is achieved by comparing their absorbance with a reference curve prepared with known standards.

*The antibodies used in this test kit only recognise the biologically relevant L-forms of Metanephrines. Commercially available synthetic Normetanephrine or Metanephrine is always a mixture of the D- and L-forms. The ratio between both forms differs widely from lot to lot. This has important implications if synthetic Metanephrines are used to enrich native samples. As only about 50% of the synthetic Metanephrines - the L-portion - will be detected by use of this kit, spiked samples will be underestimated. Therefore native samples containing solely the L-form should be used.*

### 1.2 Clinical application

Metanephrine and Normetanephrine are the metabolites of the catecholamines Epinephrine and Norepinephrine, respectively. Cells derived from neuroendocrine tumors (e.g. pheochromocytoma) are known to produce catecholamines which are secreted episodically via vesicles into the blood stream. But beside this a small portion of the catecholamines is metabolized inside the cells to the corresponding catecholamines metabolites – namely Metanephrine, Normetanephrine and 3-Methoxytyramine – which are secreted at low levels continuously into the blood stream.

Recent studies and publications have shown that the quantification of these plasma free Metanephrine and plasma free Normetanephrine is the most accurate biochemical marker for the clinical diagnosis of pheochromocytoma and follow-up of pheochromocytoma patients.

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as under point "Procedural cautions, guidelines and warnings". Any laboratory result is only a part of the total clinical picture of the patient.

Only in cases where the laboratory results are in an acceptable agreement with the overall clinical picture of the patient it can be used for therapeutic consequences.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

## 2 PROCEDURAL CAUTIONS, GUIDELINES, WARNINGS AND LIMITATIONS

### 2.1 Procedural cautions, guidelines and warnings

1. This kit is intended for professional use only. Users should have a thorough understanding of this protocol for the successful use of this kit. Only the test instruction provided with the kit is valid and has to be used to run the assay. Reliable performance will only be attained by strict and careful adherence to the instructions provided.
2. This assay was validated for a certain type of sample as indicated in *Intended Use* (please refer to Chapter 1). Any off-label use of this kit is in the responsibility of the user and the manufacturer cannot be held liable.
3. Reagents of this kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
4. The principles of Good Laboratory Practice (GLP) have to be followed.
5. In order to reduce exposure to potentially harmful substances, wear lab coats, disposable protective gloves and protective glasses where necessary.
6. All kit reagents and specimens should be brought to room temperature and mixed gently but thoroughly before use. Avoid repeated freezing and thawing of reagents and specimens.
7. For dilution or reconstitution purposes, use deionized, distilled, or ultra-pure water.
8. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch with desiccant and used in the frame provided.

9. Duplicate determination of sample is highly recommended to be able to identify potential pipetting errors.
10. Once the test has been started, all steps should be completed without interruption. Make sure that the required reagents, materials and devices are prepared ready at the appropriate time.
11. Incubation times do influence the results. All wells should be handled in the same order and time intervals.
12. To avoid cross-contamination of reagents, use new disposable pipette tips for dispensing each reagent, sample, standard and control.
13. A standard curve must be established for each run.
14. The controls should be included in each run and fall within established confidence limits. The confidence limits are listed in the QC-Report.
15. Do not mix kit components with different lot numbers within a test and do not use reagents beyond expiry date as shown on the kit labels.
16. Avoid contact with Stop Solution containing 0.25 M H<sub>2</sub>SO<sub>4</sub>. It may cause skin irritation and burns. In case of contact with eyes or skin, rinse off immediately with water.
17. Some reagents contain sodium azide (NaN<sub>3</sub>) as preservative. In case of contact with eyes or skin, rinse off immediately with water. NaN<sub>3</sub> may react with lead and copper plumbing to form explosive metal azides. When disposing reagents, flush with a large volume of water to avoid azide build-up.
18. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them.
19. For information on hazardous substances included in the kit please refer to Safety Data Sheet (SDS). The Safety Data Sheet for this product is made available upon request.
20. The expected reference values reported in this test instruction are only indicative. It is recommended that each laboratory establishes its own reference intervals.
21. The results obtained with this test kit should not be taken as the sole reason for any therapeutic consequence (e.g. medication before a scheduled surgery) but have to be correlated to other diagnostic tests and clinical observations.
22. Kit reagents must be regarded as hazardous waste and disposed according to national regulations.

## 2.2 Limitations

Any inappropriate handling of samples or modification of this test might influence the results.

### 2.2.1 Interfering substances

#### Plasma

Samples containing precipitates or fibrin strands or which are haemolytic or lipemic might cause inaccurate results.

### 2.2.2 Drug interferences

There are no known substances (drugs) which ingestion interferes with the measurement of (Nor-)metanephrine level in the sample.

### 2.2.3 High-Dose-Hook effect

No hook effect was observed in this test.

## 3 STORAGE AND STABILITY

Store the unopened reagents at 2 °C - 8 °C until expiration date.

Do not use components beyond the expiry date indicated on the kit labels.

Once opened the reagents are stable for 1 month when stored at 2 °C - 8 °C.

Once the resealable pouch has been opened, care should be taken to close it tightly with desiccant again.

## 4 MATERIALS

### 4.1 Content of the kit

	<b>REAC-TUBES</b>	<b>Reaction Tubes</b> - Ready to use
Content:	Reaction Tubes in a resealable pouch	
Volume:	2 x 50 tubes	
	<b>FOILS</b>	<b>Adhesive Foil</b> - Ready to use
Content:	Adhesive Foils in a resealable pouch	
Volume:	2 x 4 foils	
	<b>WASH-CONC 50x</b>	<b>Wash Buffer Concentrate</b> - Concentrated 50x
Content:	Buffer with a non-ionic detergent and physiological pH	
Volume:	2 x 20 mL/vial, light purple cap	
	<b>CONJUGATE</b>	<b>Enzyme Conjugate</b> - Ready to use
Content:	Goat anti-rabbit immunoglobulins conjugated with peroxidase	
Volume:	2 x 12 mL/vial, red cap	
	<b>SUBSTRATE</b>	<b>Substrate</b> - Ready to use
Content:	Chromogenic substrate containing tetramethylbenzidine, substrate buffer and hydrogen peroxide	
Volume:	2 x 12 mL/vial, black cap	
	<b>STOP-SOLN</b>	<b>Stop Solution</b> - Ready to use
Content:	0.25 M sulfuric acid	
Volume:	2 x 12 mL/vial, light grey cap	
	<b>ADR MN</b>	<b>Metanephrine Microtiter Strips</b> - Ready to use
Content:	1 x 96 well (12x8) antigen precoated microwell plate in a resealable blue pouch with desiccant	
	<b>NAD NMN</b>	<b>Normetanephrine Microtiter Strips</b> - Ready to use
Content:	1 x 96 well (12x8) antigen precoated microwell plate in a resealable yellow pouch with desiccant.	
	<b>MN-AS</b>	<b>Metanephrine Antiserum</b> - Ready to use
Content:	Rabbit anti- Metanephrine antibody, blue coloured	
Volume:	1 x 6 mL/vial, blue cap	
	<b>NMN-AS</b>	<b>Normetanephrine Antiserum</b> - Ready to use
Content:	Rabbit anti- Normetanephrine antibody, yellow coloured	
Volume:	1 x 6 mL/vial, yellow cap	
	<b>ASSAY-BUFF</b>	<b>Assay Buffer</b> - Ready to use
Content:	TRIS-buffer containing proteins and a non-mercury preservative	
Volume:	1 x 12 mL/vial, orange cap	
	<b>EQUA-REAG</b>	<b>Equalizing Reagent</b> - Lyophilized
Content:	Human serum, negative for HIV I/II, HBsAg and HCV	
Volume:	2 vials, dark green cap	

ACYL-CONC

**Acylation Concentrate** – Concentrated

Content: Acylation reagent in DMSO  
 Volume: 1 x 1.5 mL/vial, dark grey cap  
 Hazards identification:

H302 Harmful if swallowed.

**Standards and Controls** - Ready to use

Component	Colour/Cap	Concentration pg/mL		Concentration pmol/L		Volume/ Vial
		MN	NMN	MN	NMN	
STANDARD A	white	0	0	0	0	12 mL
STANDARD B	light yellow	36	48	183	262	4 mL
STANDARD C	orange	120	160	608	874	4 mL
STANDARD D	dark blue	360	480	1825	2621	4 mL
STANDARD E	light grey	1200	1600	6084	8736	4 mL
STANDARD F	black	3600	4800	18252	26208	4 mL
CONTROL 1	light green	Refer to QC-Report for expected value and				4 mL
CONTROL 2	dark red	acceptable range!				4 mL

Conversion: Metanephrine (pg/mL) x 5.07 = Metanephrine (pmol/L)

Normetanephrine (pg/mL) x 5.46 = Normetanephrine (pmol/L)

Content: Buffer with stabilizer and a precipitating reagent spiked with a defined quantity of Metanephrine and Normetanephrine

Hazards

identification:

H302 Harmful if swallowed.

**4.2 Additional materials and equipment required but not provided in the kit**

- Calibrated precision pipettes to dispense volumes between 25 - 500 µL; 3 mL; 10 mL
- Microtiter plate washing device (manual, semi-automated or automated)
- ELISA reader capable of reading absorbance at 450 nm and if possible 620 - 650 nm
- Centrifuge capable of at least 3.000 x g
- microtiter plate shaker (shaking amplitude 3 mm; approx. 600 rpm)
- Absorbent material (paper towel)
- Water (deionized, distilled, or ultra-pure)
- Vortex mixer

**5 SAMPLE COLLECTION AND STORAGE****EDTA-, Citrate- or Heparin-Plasma**

Whole blood should be collected into centrifuge tubes (Monovette™ or Vacuette™) containing EDTA, heparin or citrate as anti-coagulant and centrifuged (according to manufacturer's instructions) immediately after collection.

Haemolytic and lipemic samples should not be used for the assay.

Storage: up to 6 hours at 2 °C - 8 °C,  
 for longer period (up to 6 months) at -20 °C.

Repeated freezing and thawing should be avoided.

## 6 TEST PROCEDURE

Allow all reagents to reach room temperature and mix thoroughly by gentle inversion before use. Number the Reaction Tubes accordingly. Duplicate determinations are recommended.

The binding of the antibodies and the enzyme conjugates and the activity of the enzyme used are temperature dependent, and the absorption values may vary if a thermostat is not used. The higher the temperature, the higher the absorption values will be. The absorption values also depend on the incubation times. The optimal temperature during the Enzyme Immunoassay is between 20 °C - 25 °C.

*In case of overflow, read the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to 405 nm*

### 6.1 Preparation of reagents

#### Wash Buffer

Dilute the 20 mL Wash Buffer Concentrate with water (deionized, distilled, or ultra-pure) to a final volume of 1000 mL. Storage: 1 month 2 °C - 8 °C

#### Equalizing Reagent

Reconstitute the Equalizing Reagent with 10 mL water (deionized, distilled, or ultra-pure).

Reconstituted Equalizing Reagent which is not used immediately has to be stored in aliquots for max. 1 month at -20 °C and may be thawed only once.

#### Acylation Solution

As the Acylation Solution is only stable for a maximum of 3 minutes it should not be prepared before starting the assay. Therefore its preparation is described in the protocol in chapter 6.3, step 3 and chapter 6.4, step 3.

Discard after use!

### 6.2 Preparation of samples

*The precipitation procedure is the same for Metanephrine and Normetanephrine and has to be done only once.*

#### 6.2.1 Precipitation

1. Pipette **100 µL** of **standards**, **100 µL** of **controls**, and **500 µL** of **plasma samples** into the respective **Reaction Tubes**.
2. Add **500 µL Equalizing Reagent** (*refer to 6.1*) to all tubes containing standards and controls.
3. Add **100 µL Standard A** to all tubes containing plasma samples.
4. Mix **Reaction Tubes** thoroughly (vortex) and centrifuge for **15 minutes** at **3,000 x g**.

Take **75 µL** of the clear supernatant for the Metanephrine ELISA and **25 µL** of the clear supernatant for the Normetanephrine ELISA.

### 6.3 Metanephrine ELISA

1. Pipette **50 µL** of **Assay Buffer** into the appropriate wells of the **Metanephrine Microtiter Strips**.
2. Pipette **75 µL** of the **clear supernatant** from the **standards, controls and samples** into the wells.
3. Preparation of **Acylation Solution**:  
Pipette **80 µL Acylation Reagent Concentrate** to **3 mL water** (deionized, distilled, or ultra-pure) and mix thoroughly.
4. Pipette **25 µL** of the freshly prepared **Acylation Solution** into all wells.
5. Incubate for **15 min** at **RT** (20 °C - 25 °C) on a **shaker** (approx. 600 rpm).
6. Pipette **50 µL** of the **Metanephrine Antiserum** into all wells.
7. Cover the plate with **Adhesive Foil**, shake for **1 min** at **RT** (20 °C - 25 °C) on a **shaker** and incubate for **15 - 20 h** (overnight) at **2 °C - 8 °C** without shaking.
8. Remove the foil. Discard or aspirate the content of the wells. Wash the plate **4 x** by adding **300 µL** of **Wash Buffer**, **discarding** the content and **blotting dry each time** by tapping the inverted plate on absorbent material.
9. Pipette **100 µL** of the **Enzyme Conjugate** into all wells.
10. Incubate for **30 min** at **RT** (20 °C - 25 °C) on a **shaker** (approx. 600 rpm).
11. Discard or aspirate the content of the wells. Wash the plate **4 x** by adding **300 µL** of **Wash Buffer**, **discarding** the content and **blotting dry each time** by tapping the inverted plate on absorbent material.
12. Pipette **100 µL** of the **Substrate** into all wells and incubate for **20 - 30 min** at **RT** (20 °C - 25 °C) on a **shaker** (approx. 600 rpm). **Avoid exposure to direct sunlight!**
13. Add **100 µL** of the **Stop Solution** to each well and shake the microtiter plate to ensure a homogeneous distribution of the solution.
14. **Read** the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to **450 nm** (if available a reference wavelength between 620 nm and 650 nm is recommended).

### 6.4 Normetanephrine ELISA

1. Pipette **50 µL** of **Assay Buffer** into the appropriate wells of the **Normetanephrine Microtiter Strips**.
2. Pipette **25 µL** of the **clear supernatant** from the **standards, controls and samples** into the wells.
3. Preparation of **Acylation Solution**:  
Pipette **80 µL Acylation Reagent Concentrate** to **3 mL water** (deionized, distilled, or ultra-pure) and mix thoroughly.
4. Pipette **25 µL** of the freshly prepared **Acylation Solution** into all wells.
5. Incubate for **15 min** at **RT** (20 °C - 25 °C) on a **shaker** (approx. 600 rpm).
6. Pipette **50 µL** of the **Normetanephrine Antiserum** into all wells.
7. Cover the plate with **Adhesive Foil**, shake for **1 min** at **RT** (20 °C - 25 °C) on a **shaker** and incubate for **15 - 20 h** (overnight) at **2 °C - 8 °C** without shaking.
8. Remove the foil. Discard or aspirate the content of the wells. Wash the plate **4 x** by adding **300 µL** of **Wash Buffer**, **discarding** the content and **blotting dry each time** by tapping the inverted plate on absorbent material.
9. Pipette **100 µL** of the **Enzyme Conjugate** into all wells.
10. Incubate for **30 min** at **RT** (20 °C - 25 °C) on a **shaker** (approx. 600 rpm).
11. Discard or aspirate the content of the wells. Wash the plate **4 x** by adding **300 µL** of **Wash Buffer**, **discarding** the content and **blotting dry each time** by tapping the inverted plate on absorbent material.
12. Pipette **100 µL** of the **Substrate** into all wells and incubate for **20 - 30 min** at **RT** (20 °C - 25 °C) on a **shaker** (approx. 600 rpm). **Avoid exposure to direct sunlight!**
13. Add **100 µL** of the **Stop Solution** to each well and shake the microtiter plate to ensure a homogeneous distribution of the solution.
14. **Read** the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to **450 nm** (if available a reference wavelength between 620 nm and 650 nm is recommended).

## 7 CALCULATION OF RESULTS

Measuring range	Metanephrine	Normetanephrine
		17 – 3600 pg/mL

The standard curve is obtained by plotting the absorbance readings (calculate the mean absorbance) of the standards (linear, y-axis) against the corresponding standard concentrations (logarithmic, x-axis).

Use a non-linear regression for curve fitting (e.g. spline, 4-parameter, akima).

The concentrations of the **samples** and **controls** can be read directly from the standard curve.

Samples found with concentrations higher than the highest standard (Standard F) should be diluted accordingly with Equalizing Reagent and have to be re-assayed.

**Conversion** Metanephrine (pg/mL) x 5.07 = Metanephrine (pmol/L)  
Normetanephrine (pg/mL) x 5.46 = Normetanephrine (pmol/L)

### Expected reference values

It is strongly recommended that each laboratory should determine its own reference values.


Metanephrine	Normetanephrine
< 90 pg/mL	< 180 pg/mL

### 7.1 Quality control

It is recommended to use control samples according to national regulations. Use controls at normal and pathological levels. The kit or other commercial controls should fall within established confidence limits. The confidence limits of the kit controls are indicated in the QC-Report.

### 7.2 Typical standard curves

*Examples, do not use for calculation!*

## 8 ASSAY CHARACTERISTICS

Analytical Sensitivity (Limit of Detection)		Metanephrine	Normetanephrine
	Plasma	17 pg/mL	23 pg/mL

Analytical Specificity (Cross Reactivity)	Substance	Cross Reactivity (%)	
		Metanephrine	Normetanephrine
	Derivatized Metanephrine	100	0.08
	Derivatized Normetanephrine	0.04	100
	3-Methoxytyramin.HCl	< 0.001	1.74
	Adrenaline	< 0.001	< 0.001
	Noradrenaline	< 0.001	< 0.001
	Dopamin.HCl	< 0.001	< 0.001
	VMS	< 0.001	< 0.001
	HMVS	< 0.001	< 0.001
	L-DOPA	< 0.001	< 0.001
	L-Tyrosin	< 0.001	< 0.001
	Tyramine.HCl	< 0.001	< 0.001
	Normetanephrine	< 0.001	< 0.001
	Acetaminophen	< 0.001	< 0.001

Precision							
Intra-Assay				Inter-Assay			
	Sample	Range (pg/mL)	CV (%)		Sample	Range (pg/mL)	CV (%)
Metanephrine	1	155 ± 17	11	Metanephrine	1	133 ± 13	10
	2	245 ± 28	11		2	257 ± 23	8.9
	3	523 ± 38	7.3		3	528 ± 50	9.5
Normetanephrine	1	167 ± 12	7.3	Normetanephrine	1	161 ± 18	11
	2	373 ± 29	7.8		2	370 ± 19	5.1
	3	832 ± 60	7.2		3	844 ± 51	6.0

Linearity			Range (pg/mL)	Serial dilution up to	Mean (%)
	Metanephrine	Plasma	40 – 2100	1: 65	99
	Normetanephrine	Plasma	58 – 5800	1: 129	109

Recovery			Mean (%)	Range (%)	% Recovery after spiking
	Metanephrine	Plasma	97	85 – 111	
	Normetanephrine	Plasma	92	80 – 108	

Method Comparison: ELISA vs. LC-MS/MS	Metanephrine	Plasma	LC-MS/MS = 1.1x – 19.4	r = 0.98; n = 59
	Normetanephrine	Plasma	LC-MS/MS = 1.2x + 10.5	r = 0.99; n = 59

## 9 REFERENCES/LITERATURE

- Jeyaraman et al. The role of urinary fractionated metanephrines in the diagnosis of pheochromocytoma. The Indian Journal of Medical Research, 137(2):316-323 (2013)
- Unger et al. Plasma and Urinary Metanephrines Determined by an Enzyme Immunoassay, but not Serum Chromogranin A for the Diagnosis of Pheochromocytoma in Patients with Adrenal Mass. Exp Clin Endocrinol Diabetes, 120:494-500 (2012)
- Stefanescu et al. Salivary free catecholamines metabolites as possible biochemical markers in pheochromocytoma diagnosis. Acta Endocrinologica (Buc), VII (4): 431-442 (2011)

For updated literature or any other information please contact your local supplier.

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
## 9 REFERENZEN/LITERATUR

- (1) Jeyaraman et al. The role of urinary fractionated metanephrines in the diagnosis of pheochromocytoma. The Indian Journal of Medical Research, 137(2):316-323 (2013)
- (2) Unger et al. Plasma and Urinary Metanephrines Determined by an Enzyme Immunoassay, but not Serum Chromogranin A for the Diagnosis of Pheochromocytoma in Patients with Adrenal Mass. Exp Clin Endocrinol Diabetes, 120:494-500 (2012)
- (3) Stefanescu et al. Salivary free catecholamines metabolites as possible biochemical markers in pheochromocytoma diagnosis. Acta Endocrinologica (Buc), VII (4): 431-442 (2011)

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