Instructions for Use

SARS-CoV-2 (RBD) IgG ELISA

IVD

REF EIA-6150

Σ 96

DRG Instruments GmbH, Germany
Frauenbergstraße 18, 35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

Distributed by:

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

**Please use only the valid version of the Instructions for Use provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
Por favor, use sólo la versión válida de las instrucciones de uso que se suministran con el kit.
Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.**

Table of Contents / Inhaltsverzeichnis / Tabella dei Contenuti / Tabla de Contenidos

1	INTENDED USE.....	2	1	FINALIDAD PREVISTA	26
2	PRINCIPLE OF THE TEST	2	2	FUNDAMENTO DEL ENSAYO	26
3	WARNINGS AND PRECAUTIONS	3	3	PRECAUCIONES	26
4	REAGENTS.....	4	4	COMPONENTES DEL KIT	27
5	SPECIMEN COLLECTION AND PREPARATION	5	5	MUESTRAS	28
6	ASSAY PROCEDURE.....	5	6	PROCEDIMIENTO DE ENSAYO	28
7	RESULTS	7	7	RESULTADOS.....	30
8	QUALITY CONTROL.....	7	8	CONTROL DE CALIDAD.....	30
9	PERFORMANCE CHARACTERISTICS	8	9	CARACTERÍSTICAS DEL ENSAYO	31
10	LIMITATIONS OF USE.....	10	10	LIMITACIONES DE USO.....	31
11	LEGAL ASPECTS	10	11	ASPECTOS LEGALES	32
1	ZWECKBESTIMMUNG	11	12	REFERENCES / LITERATURE.....	33
2	TESTPRINZIP	11			
3	VORSICHTSMAßNAHMEN	12			
4	BESTANDTEILE DES KITS	13			
5	PROBENVORBEREITUNG	14			
6	TESTDURCHFÜHRUNG	14			
7	ERGEBNISSE	16			
8	QUALITÄTSKONTROLLE.....	16			
9	ASSAY-CHARAKTERISTIKA.....	17			
10	GRENZEN DES TESTS.....	17			
11	RECHTLICHE GRUNDLAGEN	18			
				SYMBOLS USED	34
1	DESTINAZIONE D'USO.....	19			
2	PRINCIPIO DEL TEST	19			
3	PRECAUZIONI	19			
4	COMPONENTI DEL KIT.....	20			
5	CAMPIONI.....	21			
6	ATTUAZIONE DEL TEST.....	21			
7	RISULTATI	23			
8	CONTROLLO QUALITÀ.....	23			
9	CARATTERISTICHE DEL TEST.....	23			
10	LIMITAZIONE DEL TEST.....	24			
11	ASPETTI LEGALI	25			

1 INTENDED USE

The **DRG SARS-CoV-2 (RBD) IgG ELISA** is an enzyme immunoassay for the **quantitative *in vitro* diagnostic** measurement of IgG antibodies to SARS-CoV-2 in serum or plasma (EDTA, lithium heparin or citrate plasma)

The DRG SARS-CoV-2 (RBD) IgG ELISA is intended for use as an aid in identifying individuals with an adaptive immune response to SARS-CoV-2, indicating recent or prior infection. It can support the diagnosis of COVID-19 disease and supplement direct pathogen detection via real-time reverse transcription polymerase chain reaction (RT-PCR) or pulmonary CT-Scans.

In addition, quantification of anti-SARS-CoV-2 RBD IgG antibodies can aid to monitor development of neutralizing antibodies after vaccination as well as support epidemiological studies on the disease prevalence.

1.1 Summary and Explanation

Coronavirus disease 2019 (COVID-19) is an infectious disease caused by the newly discovered Severe Acute Respiratory Syndrome Corona Virus 2 (SARS-CoV-2), which was first identified during an outbreak of respiratory illness cases in Wuhan City, China in December 2019 and has been declared a global pandemic by the WHO in March 2020^{1,2}. SARS-CoV-2 is a positive-sense single-stranded RNA virus and belongs to the Betacoronavirus Genus, which also includes SARS-CoV (2003) and MERS (2012), but also other human coronaviruses like strains HKU1, OC43, NL63, and 229E which cause common cold with normally mild symptoms typically throughout winter month³.

SARS-CoV-2 is predominantly transmitted by droplet infection via coughing or sneezing and through close contact with infected patients. Smear infection through contaminated surfaces seem to play a minor role. Infection mainly occurs in the respiratory tract, but may also be caused through conjunctiva of the eyes. The incubation of COVID-19 ranges between 1-14 days, with the majority of cases manifesting within 4 to 6 days. The clinical manifestations of COVID-19 vary widely from asymptomatic courses to severe pneumonia with lung failure and death. However, in around 80% of cases the symptoms are mild to moderate including fever, cough, respiratory problems, fatigue or temporarily loss of taste. Reported fatality rates increase with age, comorbidities (such as diabetes, cardiovascular or pulmonary disease, chronic liver disease, cancer) and genetic predisposition⁴.

The coronavirus spike (S) glycoprotein forms a homo-trimeric class I viral fusion complex on the outer envelope of the virion that mediates the binding of the virus to its host cells. Subunit S1 of the spike protein contains a receptor binding domain (RBD) which interacts with the human ACE-2 receptor on the surface of the host cell and initiates the entry of the virus⁵⁻⁷. Antibodies which are directed against the RBD protein of SARS-CoV-2 are discussed to have neutralizing effects and to protect host cells from viral infection^{8,9}.

Antibody response to SARS-CoV-2 infection can be detected as early as 3 days after symptom onset and reaches a peak after 2-3 weeks. Seroconversion for IgG and IgM occurs simultaneously or sequentially^{10,11}. In consequence, the sensitivity of serological assay increases from the first week of symptom onset (< 40%) to 100% > 15 days after onset¹². Antibody titers have been shown to decrease slowly post-infection¹³.

Diagnosis of COVID-19 mainly relies on real-time reverse transcription polymerase chain reaction (RT-PCR) testing of respiratory specimens. However, PCR can't identify people who went through an infection, recovered, and cleared the virus from their bodies. Furthermore, chest X-rays, pulmonary computed tomography (CT) scans, and lung ultrasounds are important tools in the early diagnosis of pneumonia in patients with COVID-19⁴.

2 PRINCIPLE OF THE TEST

The DRG SARS-CoV-2 (RBD) IgG ELISA is a solid phase enzyme-linked immunosorbent assay (ELISA).

Microtiter wells as a solid phase are coated with recombinant Receptor-Binding-Domain (RBD) of the Spike Protein S1 of SARS-CoV-2.

Diluted patient specimens and **ready-for-use Standards and Controls** are pipetted into these wells. During incubation, SARS-CoV-2 (RBD) specific antibodies of standards/controls and positive specimens are bound to the immobilized RBD antigen.

After a washing step to remove unbound sample and control material, Enzyme conjugate, which is horseradish peroxidase conjugated anti-human IgG antibodies is dispensed into the wells. During a second incubation, this anti-IgG conjugate binds specifically to human IgG antibodies resulting in the formation of enzyme-linked immune complexes.

After a washing step to remove all unbound substances, the solid phase is incubated with the substrate solution. The colorimetric reaction is stopped by addition of stop solution, and optical density (OD) of the resulting yellow product is measured. The intensity of color is proportional to the concentration of the analyte in the sample.

A standard curve is constructed by plotting OD values against concentrations of standards, and concentrations of unknown samples are determined using this standard curve.

3 WARNINGS AND PRECAUTIONS

1. This kit is for in vitro diagnostic use only. For professional use only.
2. All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
3. Before starting the assay, read the instructions completely and carefully. Use the valid version of instructions for use provided with the kit. Be sure that everything is understood.
4. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
5. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
6. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution coloured. Do not pour reagents back into vials as reagent contamination may occur.
7. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
8. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
9. Allow the reagents to reach room temperature (20 °C to 26 °C) before starting the test. Temperature will affect the optical density readings of the assay.
10. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
11. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
12. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
13. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
14. Do not use reagents beyond expiry date as shown on the kit labels.
15. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
16. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
17. Avoid contact with Stop Solution containing 0.5 M H₂SO₄. It may cause skin irritation and burns.
18. Some reagents contain ProClin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
19. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
20. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
21. For information on hazardous substances included in the kit please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from DRG.

4 REAGENTS

4.1 Reagents provided

1. **Microtiterwells**, 12 x 8 (break apart) strips, 96 wells;
Wells coated with SARS-CoV-2 RBD antigen.
2. **Sample Diluent**, 1 vial, 100 mL, ready to use;
colored green
Contains non-mercury preservative.
3. **Standard (Standard 0 - 3)**, 4 vials, 1.0 mL each, ready to use;
Concentrations: 0 – 5 – 35 – 70 DU/mL
Conversion: 1 DU/mL = 5.15 IU/mL
The standards are calibrated against the following reference material: 1st WHO IS for anti-SARS-CoV-2 immunoglobulin (human) NIBSC code: 20/136
Contain non-mercury preservative.
4. **Pos. Control**; 1 vial, 1.0 mL each, ready to use;
For control values and ranges please refer to vial label or *Certificate of Analysis*.
Contain non-mercury preservative.
5. **Neg. Control**; 1 vial, 1.0 mL each, ready to use;
For control values and ranges please refer to vial label or *Certificate of Analysis*.
Contain non-mercury preservative.
6. **Enzyme Conjugate**, 1 vial, 14 mL, ready to use;
colored red,
antibody to human IgG conjugated to horseradish peroxidase;
Contains non-mercury preservative.
7. **Substrate Solution**, 1 vial, 14 mL, ready to use;
Tetramethylbenzidine (TMB).
8. **Stop Solution**, 1 vial, 14 mL, ready to use;
Contains 0.5 M H₂SO₄,
Avoid contact with the stop solution. It may cause skin irritations and burns.
9. **Wash Solution**, 1 vial, 30 mL (40X concentrated);
See "Reagent Preparation".

4.2 Materials required but not provided

- A calibrated microtiter plate reader (450 nm, with reference wavelength at 620 nm to 630 nm)
- Calibrated variable precision micropipettes
- Incubator 37 °C
- Manual or automatic equipment for rinsing microtiter plates
- Absorbent paper
- Distilled water
- Timer
- Graph paper or software for data reduction

4.3 Storage Conditions

When stored at 2 °C to 8 °C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date.

Opened reagents must be stored at 2 °C to 8 °C. Microtiter wells must be stored at 2 °C to 8 °C. Once the foil bag has been opened, care should be taken to close it tightly again.

Opened kits retain activity for 8 weeks if stored as described above.

4.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature (20 °C to 26 °C) prior to use.

Wash Solution

Add distilled water to the 40X concentrated Wash Solution.

Dilute 30 mL of concentrated *Wash Solution* with 1170 mL distilled water to a final volume of 1200 mL.

The diluted Wash Solution is stable for 1 week at room temperature.

4.5 Disposal of the Kit

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet, Section 13.

4.6 Damaged Test Kits

In case of any damage to the test kit or components, DRG must be informed in writing, at the latest one week after receiving the kit. Damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed of according to the official regulations.

5 SPECIMEN COLLECTION AND PREPARATION

Serum or plasma (EDTA, lithium heparin or citrate plasma) can be used in this assay.

Note: Samples containing sodium azide should not be used in the assay.

In general, it should be avoided to use haemolytic, icteric, or lipaemic specimens. For further information refer to chapter "Interfering Substances".

5.1 Specimen Collection

Serum:

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

Plasma:

Whole blood should be collected into centrifuge tubes containing anticoagulant (e.g. Sarstedt Monovette with the appropriate plasma preparation) and centrifuged immediately after collection.

5.2 Specimen Storage and Preparation

Specimens should be capped and may be stored for up to 7 days at 2 °C to 8 °C prior to performing the assay.

Specimens stored for a longer time (up to 6 months) should be frozen only once at -20 °C prior to assay. Thawed samples should be inverted several times prior to testing.

Diluted samples (1+100; see chapter 5.3) can be stored as well for up to 7 days at 2 °C to 8 °C.

5.3 Specimen Dilution

Prior to assaying, dilute each patient specimen **1 + 100** with *Sample Diluent* (mix well).

Example: 10 µL sample + 1 mL of Sample Diluent (mix well)

For patient samples with concentrations greater than Standard 3, a second 1:10 dilution of this 1 + 100 diluted patient sample should be performed.

Example: 20 µL of first sample dilution + 180 µL *Sample Diluent* (mix well).

For the calculation of the concentrations this second dilution factor has to be taken into account.

Note: Standards and Controls are ready to use and do not have to be diluted!

6 ASSAY PROCEDURE

6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination.
- Optical density is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.

6.2 Test Procedure

Each run must include a standard curve.

1. Secure the desired number of Microtiter wells in the frame holder.
2. Dispense **100 µL** of each **Standard, Control** and **pre-diluted sample** with new disposable tips into appropriate wells.
3. Incubate for **60 minutes at 37 °C**.
4. Rinse the wells **3 times** with **400 µL** diluted *Wash Solution* per well, if a plate washer is used.
- OR -
 Briskly shake out the contents of the wells.
 Rinse the wells **3 times** with **300 µL** diluted *Wash Solution* per well for manual washing. Strike the wells sharply on absorbent paper to remove residual droplets.
5. Dispense **100 µL Enzyme Conjugate** into each well.
6. Incubate for **30 minutes** at room temperature (20 °C to 26 °C).
7. Rinse the wells **3 times** with **400 µL** diluted *Wash Solution* per well, if a plate washer is used.
- OR -
 Briskly shake out the contents of the wells.
 Rinse the wells **3 times** with **300 µL** diluted *Wash Solution* per well for manual washing. Strike the wells sharply on absorbent paper to remove residual droplets.
Important note:
 The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
8. Add **100 µL of Substrate Solution** to each well.
9. Incubate for **15 minutes** at room temperature.
10. Stop the enzymatic reaction by adding **50 µL of Stop Solution** to each well.
11. Determine the optical density of the solution in each well at **450 nm (reading) and at 620 nm to 630 nm (background subtraction, recommended)** with a microtiter plate reader.
 It is recommended that the wells be read **within 10 minutes** after adding the *Stop Solution*.

6.3 Calculation of Results

1. Calculate the average optical density (OD) values for each set of standards, controls and patient samples.
2. Automated method:
 The results in the Instructions for Use have been calculated automatically using a 4-Parameter curve fit. (4-Parameter Rodbard or 4-Parameter Marquardt are the preferred methods.) Other data reduction functions may give slightly different results.
3. If no automated method is available:
 Using linear or semi-logarithmic graph paper, construct a standard curve by plotting the mean OD obtained from each standard against its concentration with OD value on the vertical (Y) axis and concentration on the horizontal (X) axis. Using the mean OD value for each sample determine the corresponding concentration from the standard curve.
4. The concentration of the samples can be **read directly** from this standard curve.
 Samples with concentrations greater than that of the highest standard can be further diluted 1:10 or reported as > 70 DU/mL. For the calculation of the concentrations this additional dilution factor (x 10) has to be taken into account.

6.3.1 Example of Typical Standard Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Standard		Optical Density (450 nm)
Standard 0	(0 DU/mL)	0.014
Standard 1	(5 DU/mL)	0.314
Standard 2	(35 DU/mL)	1.460
Standard 3	(70 DU/mL)	2.348

7 RESULTS

7.1 Validation of the Test Run

The test run may be considered valid provided the control values and ranges given on the labels and the *Certificate of Analysis* are met.

7.2 Interpretation of Results

Normal value ranges for this ELISA should be established by each laboratory based on its own patient populations in the geographical areas.

The test results are given in DU/mL (**DRG Units**) and are calibrated against the 1st WHO International Standard for anti-SARS-CoV-2 immunoglobulin (human), NIBSC code 20/136.

The **Cut-off value** is determined with 5.40 DU/mL.

The following values should be considered as a guideline:

Antibody Conc. (DU/mL)	Interpretation	
< 4.86	Negative	
4.86 – 5.94	Grey zone (equivocal)	For a test result in the grey zone, repeat the test after e.g. 2 weeks with a fresh patient sample. If result in the second test is again found in the grey zone, the sample is considered negative .
> 5.94	Positive	

8 QUALITY CONTROL

Good laboratory practice requires that controls be run with each standard curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day-to-day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or DRG directly.

9 PERFORMANCE CHARACTERISTICS

9.1 Assay Dynamic Range

The range of the assay is between 0.161 DU/mL - 70.0 DU/mL.

9.2 Specificity of Antigen (Cross-Reactivity)

The following substances were tested for cross-reactivity of the assay:

Sample positive to	n
Human Coronavirus NL63 Nucleoprotein*	5
Human Coronavirus 229E Nucleoprotein IgG*	9
MERS IgG*	3
Coronavirus HKU1 IgG*	13
Coronavirus OC43 IgG*	12
Human Coronavirus 229E/NL63 Nucleoprotein IgG*	5
Human Coronavirus NL63 Nucleoprotein and Parainfluenza 1/2/3 IgG*	1
Human Coronavirus 229E/NL63 Nucleoprotein IgG and Parainfluenza 1/2/3 IgG*	1
Adenovirus IgG	20
Bordetella pertussis IgG	16

Sample positive to	n
Chlamydia pneumonia IgG	9
Chlamydia trachomatis IgG	9
EBV IgG	16
Enterovirus IgA	5
Enterovirus IgG	6
Enterovirus IgM	2
Influenza A IgG	7
Influenza B IgG	7
Mycoplasma pneumonia IgG	15
Parainfluenza 1/2/3 IgG	16
RSV IgG	19

*marked samples were used as Coronavirus panel for evaluation of diagnostic specificity.

Sample positive to	n
Rheumatoid factor IgM	6
Rheumatoid factor IgA	4
different ANA	19

Sample	n
Pregnant	18

The antigen used for the DRG SARS-CoV-2 (RBD) IgG ELISA shows no cross-reactivity to the evaluated antibodies.

9.3 Sensitivity

The analytical sensitivity of the DRG ELISA was calculated by adding 2 standard deviations to the mean OD of 20 replicate analyses of *Standard 0* and was found to be 0.009 DU/mL.

The Limit of Blank (LoB) is 0.008 DU/mL.

The Limit of Detection (LoD) is 0.046 DU/mL.

The Limit of Quantification (LoQ) is 0.161 DU/mL.

9.4 Diagnostic Specificity, Sensitivity and Cut-Off

ROC-Analysis is used for the determination.

Using the Cut-Off 5.40 DU/mL, the following specificities were found.

Cohort	n	Non-reactive	Reactive	Specificity (%)
Diagnostic Routine panel	194	194	0	100.0
Blood Donors panel	128	128	0	100.0
Coronavirus panel*	49	49	0	100.0
Total	371	371	0	100.0

* See chapter 9.2 Specificity of Antigen

Using the Cut-Off 5.40 DU/mL, the following sensitivity was found.

Days post PCR confirmation	n	Non-reactive	Reactive	Sensitivity (%)
> 14	84	4	80	95.2

Please note that the sensitivity decreases for samples collected earlier than 14 days after infection with the coronavirus (9-12).

9.5 Reproducibility

9.5.1 Intra-Assay

The within-assay variability was determined by measuring each sample 10 times in one run (n = 10):

Sample	n	Mean (DU/mL)	CV (%)
1	10	3.6	7.4
2	10	13.9	5.7
3	10	23.2	5.4
4	10	40.7	5.6

9.5.2 Inter-Assay

The between-assay variability was determined by measuring each sample 10 times per run on 3 days (n = 30):

Sample	n	Mean (DU/mL)	CV (%)
1	30	3.9	13.4
2	30	13.6	7.4
3	30	24.0	7.0
4	30	42.2	11.2

9.5.3 Inter-Lot

The inter-assay (between-lots) variation was determined by measuring each sample 6 times with 3 different kit lots:

Sample	n	Mean (DU/mL)	CV (%)
1	18	7.9	8.0
2	18	9.6	3.9
3	18	24.9	8.1
4	18	60.9	6.9

9.6 Recovery

Samples have been spiked by adding solutions with known concentrations of analyte.

The recovery (%) was calculated by multiplying the ratio of measured and expected values with 100.

		Sample 1	Sample 2	Sample 3	Sample 4
Concentration (DU/mL)		3.7	14.6	23.7	44.8
Average Recovery (%)		103.5	100.0	94.7	99.2
Range of Recovery (%)	from	98.4	88.0	87.6	92.8
	to	108.1	109.5	98.1	106.3

9.7 Linearity

Samples were measured undiluted and in serial dilutions with *Sample Diluent*. The recovery (%) was calculated by multiplying the ratio of expected and measured values with 100.

		Sample 1	Sample 2	Sample 3	Sample 4
Concentration (DU/mL)		30.5	31.2	32.2	51.1
Average Recovery (%)		105.6	103.4	101.4	100.0
Range of Recovery (%)	from	101.8	102.6	98.7	93.8
	to	109.6	104.8	105.4	104.7

10 LIMITATIONS OF USE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice.

Any improper handling of samples or modification of this test might influence the results.

Bacterial or fungal contamination of patient samples or reagents, or cross-contamination between reagents may cause erroneous results.

This test is for quantitative detection of anti-SARS-CoV-2 IgG. Test results should not be the sole basis for clinical diagnosis and treatment.

The confirmation of infection with SARS-CoV-2 must be combined with the patient's clinical signs in conjunction to other tests.

In the first week of the onset of the infection with SARS-CoV-2, patient results may be negative for antibodies. The sensitivity increases over the first two weeks and reaches a plateau after two to three weeks.

In addition, false negative results for anti-SARS-CoV-2 antibodies may be obtained for patients with low immunity or other diseases that affect immune function, for patients with failure of important systemic organs, and for patients who use of drugs that suppress immune function.

10.1 Interfering Substances

Haemoglobin (up to 4 mg/mL), Bilirubin (up to 0.5 mg/mL) and Triglyceride (up to 7.5 mg/mL) have no influence on the assay results.

11 LEGAL ASPECTS

11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover, the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2 are also invalid.

Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

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

12 REFERENCES / LITERATURE

1. Ji T. et al. Detection of Covid-19: A review of the current literature and future perspectives. *Biosensors and Bioelectronics*. 2020 166; 112455.
2. World Health Organization. Coronavirus disease 2019 (COVID-19) situation report-85, Data as received by WHO from national authorities by 10:00 CET, 14 April 2020. Accessed April 14, 2020. <https://www.who.int/docs/default-source/coronaviruse/situation-reports/20200414-sitrep-85-covid-19>.
3. Zhong N et al. Epidemiology and cause of severe acute respiratory syndrome (SARS) in Guangdong People's Republic of China in February 2003. *Lancet*. 2003 362:1353-1358
4. da Silva SJR et al. Clinical and laboratory diagnosis of SARS-CoV-2, the virus causing COVID-19. *ACS Inf Dis*. August 2020; <https://dx.doi.org/10.1021/acsinfectdis.Oc00274>.
5. Brielle ES, Schneidman-Duhovny D, and Linial M. The SARS-CoV-2 exerts a distinctive strategy for interacting with the ACE2 human receptor. *Viruses*. 2020 12(5):497.
6. Gui M et al. Cryo-electron microscopy structures of the SARS-CoV spike glycoprotein reveal a prerequisite conformational state for receptor binding. *Cell Research*. 2017 27:119-129.
7. Turonova B et al. In situ structural analysis of SARS-CoV-2 spike reveals flexibility mediated by three hinges. *Science*. 2020 10.1126/science.abd.5223.
8. Premkumar L et al. The receptor binding domain of the viral spike protein is an immunodominant and highly specific target of antibodies in SARS-CoV-2 patients. *Science Immunology* 2020 10.1126/sciimmunol.abc8413.
9. Robbiani DF et al. Convergent antibody responses to SARS-CoV-2 in convalescent individuals. *Nature*. 2020 584:437-456.
10. Kontou PI et al. Antibody tests in detecting SARS-CoV-2 infections: a meta-analysis. *Diagnostics*. 2020 10:319-34.
11. Long QX et al. Antibody responses to SARS-CoV-2 in patients with COVID-19. *Nature Medicine*. 2020 26:845-848.
12. Zhao J et al. Antibody responses to SARS-CoV-2 in patients of novel coronavirus disease 2019. *Clin Inf Dis*. 2020 Accepted manuscript doi:10.1093/cid/ciaa344.
13. Ward, H., Cooke, G., Atchison, C. J., Whitaker, M., Elliott, J., Moshe, M., ... & Ashby, D. (2020). Declining prevalence of antibody positivity to SARS-CoV-2: a community study of 365,000 adults. medRxiv.

SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français
	European Conformity	CE-Konformitäts-kennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation
	<i>In vitro</i> diagnostic medical device *	<i>In-vitro</i> -Diagnostikum *	Dispositivo medico-diagnostico in vitro	Producto sanitario para diagnóstico in vitro	Dispositif médical de diagnostic in vitro
	Catalogue number *	Artikelnummer *	No. di Cat.	No de catálogo	Référence
	Batch code *	Fertigungslosnummer, Charge *	Lotto no	Número de lote	No. de lot
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen *	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests
	Temperature limit *	Temperaturbegrenzung *	Temperatura di conservazione	Temperatura de conservación	Température de conservation
	Use-by date *	Verwendbar bis *	Data di scadenza	Fecha de caducidad	Date limite d'utilisation
	Manufacturer *	Hersteller *	Fabbricante	Fabricante	Fabricant
	Biological risks *	Biologische Risiken *	Rischi biologici	Riesgos biológicos	Risques biologiques
	Caution *	Achtung *	Attenzione	Precaución	Attention
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches
<i>Distributed by</i>	Distributed by	Vertreiber	Distributore	Distribuidor	Distributeur
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Conditionnement
<i>Volume/No.</i>	Volume / No.	Volumen / Anzahl	Volume / Quantità	Volumen / Número	Volume / Quantité
<i>Microtiterwells</i>	Microtiterwells	Mikrotiterwells	Micropozzetti	Placas multipocillo	Microplaques
<i>Antiserum</i>	Antiserum	Antiserum	Antisiero	Antisuero	Antisérum
<i>Enzyme Conjugate</i>	Enzyme Conjugate	Enzymkonjugat	Tracciante enzimatico	Conjugado enzimático	Conjugué enzymatique
<i>Enzyme Complex</i>	Enzyme Complex	Enzymkomplex	Complesso enzimatico	Complejo enzimático	Complexe enzymatique
<i>Substrate Solution</i>	Substrate Solution	Substratlösung	Soluzione di substrato	Solución de sustrato	Solution substrat
<i>Stop Solution</i>	Stop Solution	Stopplösung	Soluzione d'arresto	Solución de parada	Solution d'arrêt
<i>Zero Standard</i>	Zero Standard	Nullstandard	Standard zero	Estándar cero	Zero Standard
<i>Standard</i>	Standard	Standard	Standard	Estándar	Standard
<i>Control</i>	Control	Kontrolle	Controllo	Control	Contrôle
<i>Assay Buffer</i>	Assay Buffer	Assaypuffer	Tampone del test	Tampón de ensayo	Tampon d'essai
<i>Wash Solution</i>	Wash Solution	Waschlösung	Soluzione di lavaggio	Solución de lavado	Solution de lavage
<i>1N NaOH</i>	1N NaOH	1N NaOH	1N NaOH (idrossido di sodio 1N)	1N NaOH	1N NaOH
<i>1 N HCl</i>	1 N HCl	1 N HCl	1 N HCl	1 N HCl	1N HCl
<i>Sample Diluent</i>	Sample Diluent	Probenverdünnungs-medium	Diluyente dei campioni	Solución para dilución de la muestra	Solution pour dilution de l'échantillon
<i>Conjugate Diluent</i>	Conjugate Diluent	Konjugatverdünnungs-medium	Diluyente del tracciante	Solución para dilución del conjugado	Solution pour dilution du conjugué
<i>Pos. Control</i>	Positive control	Positive Kontrolle	Controllo positivo	Control positivo	Contrôle positif
<i>Neg. Control</i>	Negative control	Negative Kontrolle	Controllo negativo	Control negativo	Contrôle négatif