HbA1c Direct Calibrator 5 Level Series

(en) English

REF

Content

Y04603 - 1x 1 mL (Level 1) + 4x 1 mL (Level 2 – 5)

For professional in vitro diagnostic use only.

INTENDED USE

HbA1c Direct Calibrator 5 Level Series is a set consisting of lyophilized hemolysates and buffer used for the calibration in course of the quantitative in vitro determination of Haemoglobin A1c in blood on photometric systems. For professional laboratory use.

DIAGNOSTIC SIGNIFICANCE

Please refer to the corresponding section of the used HbA1c Direct reagent.

TEST PRINCIPLE

Please refer to the corresponding section of the used Dialab HbA1c Direct reagent.

REAGENT COMPOSITION

The HbA1c Calibrator Level 1 is a liquid buffer with preservative: 0.095% sodium azide. The HbA1c Calibrators Level 2 – 5 are lyophilized, stabilized hemolysates prepared from packed human erythrocytes.

MATERIAL REQUIRED BUT NOT PROVIDED

- Deionized water for the preparation of calibrator Levels 2 – 5 and the haemolysis of patient samples.

- HbA1c Direct reagent, e.g.:

REF	Name	Content
Y04602	HbA1c Direct	4x 22.5 mL R1 + 1x 30 mL R2
Y04606	HbA1c Direct	4x 7.5 mL R1 + 1x 10 mL R2
Y04601	HbA1c Direct	4x 45 mL R1 + 1x 60 mL R2

- Controls:

REF	Name	Content
Y04604	HbA1c Direct Control Set	4 x 1 mL

- Clinical Laboratory Analyser or Photometer
- General laboratory equipment

REAGENT PREPARATION

Level 1 is ready for use.

Level 2 – 5:

- Open the vial carefully, avoiding any loss of the lyophilized material.
- Add exactly 1 mL of deionized water (inaccurate reconstitution of the calibrators and error in assay technique can cause erroneous results).
- Close the vial carefully and gently mix for 10 minutes, or until all material has dissolved, avoiding the formation of foam. Do not shake!
- Reconstituted calibrators are now ready for use. Do not lyse calibrators.

STORAGE AND STABILITY

Storage:	at 2 – 8 °C
Stability:	up to the expiration date if sealed tightly
Stability after reconstitution:	30 days at 2 – 8 °C 3 months at -20 °C

CLOSE IMMEDIATELY AFTER USE. PROTECT FROM LIGHT AND HEAT.

Do not use the product in case of suspect of microbial growth in the vial. Improper handling and/or storage of the calibrators and controls can affect results.

WARNINGS AND PRECAUTIONS

- For in Vitro Diagnostic use.
- Although this product has been tested and found non-reactive for HBsAg, HCV and anti-HIV 1&2, no known test can offer assurance that products derived from human blood will not transmit disease. Therefore, all human serum products and patient specimens should be handled in the same manner as an infectious agent. Do not pipette by mouth. Avoid contact with skin and mucous membranes.
- Please refer to the safety data sheets (SDS) and take the necessary precautions for the use of laboratory reagents.
- In case of product malfunction or altered appearance that could affect the performance, contact the manufacturer.
- In the event of an incident related to the device, report it to the manufacturer and your competent authority as required.

SPECIMEN COLLECTION AND STORAGE

Not applicable

TEST PROCEDURE

Measure the calibrators after preparation as described in the section for reagent preparation in the same manner, like haemolysed human samples in accordance with the instruction for use of the HbA1c Direct reagent and with your instrument manual. In case of using a Dialab application sheet for your instrument please note eventual recommendations regarding the calibration.

INTERPRETATION OF RESULTS

The results of the Levels should show increasing ΔA results with the increase of concentration. Use controls to verify your calibration.

QUALITY CONTROL AND CALIBRATION

Please refer to the instructions for use of the HbA1c reagent, used controls and The manual of your instrument. It is recommended to use quality controls after each calibration, at new measuring sessions, after each change of reagents and every time where a potential increased risk of malfunction can be assumed or when results indicate a potential malfunction and to repeat the calibration where required. Each laboratory should establish its own quality control scheme and corrective actions if controls do not meet the acceptable tolerances.

PERFORMANCE CHARACTERISTICS

Please refer to the corresponding section of HbA1c Direct reagent package insert.

TRACEABILITY

Values are based on reference material from NGSP/DCCT.

EXPECTED VALUES

Please refer to the corresponding section of HbA1c Direct reagent package insert.

LIMITATIONS

The product should not be used after the expiry date.

WASTE MANAGEMENT

Please refer to local legal requirements.

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Kontrollen:

REF	Name	Inhalt
Y04604	HbA1c Direct Control Set	4 x 1 mL

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