Liquid Calibrator/Control Set

ETHANOL CALIBRATOR/CONTROL SET

For the calibration and validation of enzymatic tests for the quantitative in vitro determination of ethanol in human serum, plasma or urine.

REF

Cont.

- Z05880** 1 x 1 ml Negative Calibrator
 1 x 1 ml Calibrator High:
 contains approx. 100 mg/dL ethyl alcohol
 1 x 1 ml Control Level 1:
 contains approx. 50 (45 – 55) mg/dL ethyl alcohol
 1 x 1 ml Control Level 2:
 contains approx. 300 (270 – 330) mg/dL ethyl alcohol

COMPOSITION

Contains negative human urine spiked with ethyl alcohol with sodium azide as preservative.

CALIBRATOR/CONTROL PREPARATION

- The calibrators and controls are ready to use. No reconstitution is required.
- Use calibrators/controls in the same manner as specified with the reagent and/or instrument being used.
- Allow the calibrators/controls to equilibrate to room temperature (18°C to 25°C) prior to use.
- Cap vials immediately after use and keep them tightly capped and refrigerated when not in use.

CALIBRATOR/CONTROL STABILITY AND STORAGE

Storage: 2 - 8°C
 Stability: up to the indicated expiration date
 Do not freeze!

Once the calibrators and controls are put into the sample cups and exposed to air, alcohol evaporation will take place rapidly.

Within 15 minutes concentration may decrease up to 10%, and within 30 minutes, the alcohol concentration may decrease by 20%.

Never use any old calibrators or controls that have been left on the instrument; always use fresh solutions from the bottles.

WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use only.
- The calibrators and controls contain sodium azide (0.95 g/l) as preservative. Harmful if swallowed. Avoid contact with skin and mucous membranes.
- Handle calibrators/controls as if they are potentially infectious.
- Take the necessary precautions for the use of laboratory reagents.

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