DIAQUICK DOA Multi-6 Cup

for the qualitative detection of multiple drugs of abuse in human urine

DOA Multi-6 Cup AMP, MET, BZO, COC, MOP, THC
 - REF Z06667CE Cont: 1 test, individually packed

All products contain a package insert!

For in vitro diagnostic use only. For use by medical professionals only.
 For diagnosis and therapeutic monitoring only.

INTENDED USE

Multi-test panel with integrated urine collection cup and temperature strip for the simultaneous, qualitative detection of multiple drugs and metabolites in human urine. Urine based tests for multiple drugs of abuse range from simple immunoassay tests to complex analytical procedures. The speed and sensitivity of immunoassays have made them the most widely accepted method to test urine for multiple drugs of abuse. The DIAQUICK DOA Multi-6 Cup is a lateral flow chromatographic immunoassay for the qualitative detection of following drugs of abuse without the need of additional instruments:

Parameter	Short	Calibrator Substance	Cut-off
Amphetamines	AMP	D-Amphetamine	1 000 ng/mL
Benzodiazepines	BZO	Oxazepam	300 ng/mL
Cocaine	COC	Benzoylcegonine	300 ng/mL
Methamphetamines	MET	D-Methamphetamine	1 000 ng/mL
Opiates, Heroine	MOP	Morphine	300 ng/mL
Marihuana/Cannabis	THC	11-nor-Δ9-THC-9-COOH	50 ng/mL

This test will detect other related compounds, please refer to the Analytical Specificity table in this insert. This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are obtained. For in vitro diagnostic use only

TEST PRINCIPLE

The DIAQUICK DOA Multi-6 Cup is an immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against their respective drug conjugate for binding sites on their specific antibody. During testing, a urine specimen migrates upward by capillary action. A drug, if present in the urine specimen below its cut-off concentration, will not saturate the binding sites of its specific antibody coated on the particles. The antibody coated particles will then be captured by the immobilized drug conjugate and a visible colored line will show up in the test line region of the specific drug strip. The colored line will not form in the test line region if the drug level is above its cut-off concentration because it will saturate all the binding sites of the antibody coated on the particles. A drug-positive urine specimen will not generate a colored line in the specific test line region of the strip because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

WARNINGS AND PRECAUTIONS

- For medical and other in vitro diagnostic use only. Do not use after the expiration date.
- The test should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used tests should be discarded according to federal, state and local regulations.

REAGENTS

Each test line contains anti-drug mouse monoclonal antibody and corresponding drug-protein conjugates. The control line contains goat anti-rabbit IgG polyclonal antibodies and rabbit IgG.

STORAGE

The test can be stored refrigerated or at room temperature (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

SAMPLE COLLECTION AND PREPARATION


The urine must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitations should be centrifuged, filtered or allowed to settle to obtain a clear specimen for testing. Urine specimens may be stored at 2-8°C for up to 48 h prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

ASSAY PROCEDURE



Allow the test cup, urine specimen, and/or controls to equilibrate to room temperature (15-30°C) prior to testing.

- Bring the pouch to room temperature before opening it. Remove the cup from the sealed pouch and use it as soon as possible.
- Pull tab to remove cap, collect specimen in the cup and secure cap by pressing down on all three corners.
- Check the temperature label up to 4 minutes after specimen collection. A green colour will appear to indicate the temperature of the urine specimen. The proper range for an unadulterated specimen is 32 – 38 °C.
- Check the cap for a tight seal, date and initial the security seal label, then place it over the cap.
- Remove the key from the kit, place the cup on a flat surface, and push the key into the socket of the cup to begin the test. Start timer.
- Remove the peel off label covering the test results and wait for the coloured line(s) to appear. Read results at 5 minutes. Do not interpret results after 10 minutes.

- Place cup on a flat surface, insert key and push in.
- Peel off label to view results.

- Wait 5 minutes to read drug results.
- Interpret drug results at 5 minutes.

INTERPRETATION OF RESULTS

NEGATIVE: A colored line in the control region (C) and a colored line in the test line region (T) for a specific drug indicate a negative result. This indicates that the drug concentration in the urine specimen is below the designated cut-off level for that specific drug.

*NOTE: The shade of color in the test region (T) may vary, but it should be considered negative whenever there is even a faint pink line.

POSITIVE: A colored line in the control line region (C) but no line in the test line region (T) for a specific drug indicates a positive result. This indicates that the drug concentration in the urine specimen exceeds the designated cut-off level.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test panel. If the problem persists, discontinue using the lot immediately and contact your local distributor

QUALITY CONTROL

A procedural control is included in the test. A red line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

- The DIAQUICK DOA Multi-6 Cup (urine) provide only a preliminary analytical result. A more specific chemical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.
- It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
- Adulterants, such as bleaching agents in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
- A positive result indicates presence of the drug or its metabolites but does not indicate the level of intoxication, administration route or concentration in urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained if a drug is present but below the cut-off level of the test.
- The DIAQUICK DOA Multi-6 Cup (urine) do not distinguish between drugs of abuse and certain medications.
- A positive result might be obtained from certain foods or food supplements.

PERFORMANCE CHARACTERISTICS

ACCURACY

A side-by-side comparison was conducted using the DIAQUICK DOA Multi-6 Cup (urine) and commercially available drug rapid tests. Testing was performed on approximately 250 specimens. Presumptive positive results were confirmed by GC/MS. The following results were tabulated.

% Agreement with Commercial Kit

	Positive Agreement	Negative Agreement	Total Results
AMP	> 99.9 %	> 99.9 %	> 99.9 %
BZO	> 99.9 %	> 99.9 %	> 99.9 %
COC	> 99.9 %	> 99.9 %	> 99.9 %
THC	> 99.9 %	> 99.9 %	> 99.9 %
MET	> 99.9 %	> 99.9 %	> 99.9 %
MOP	> 99.9 %	> 99.9 %	> 99.9 %

% Agreement with GC/MS

	Positive Agreement	Negative Agreement	Total Results
AMP	98.1 %	97.9 %	98.0 %
BZO	98.4 %	99.2 %	98.8 %
COC	98.2 %	97.8 %	98.0 %
THC	97.9 %	98.1 %	98.0 %
MET	96.2 %	97.1 %	96.8 %
MOP	95.0 %	95.3 %	95.2 %

ANALYTICAL SPECIFICITY

The following tables lists the concentration of compounds (ng/mL) that are detected positive in urine by the DIAQUICK DOA Multi-6 Cup (urine) at 5 minutes.

AMPHETAMINES	AMP	COCAINE	COC
D,L-Amphetamine sulfate	300	Benzoyllecgonine	300
L-Amphetamine	25 000	Cocaine HCl	200
(±) 3,4-Methylenedioxyamphetamine	500	Cocaehtylene	20 000
Phentermine	800	Ecgonine HCl	30 000
Maprotiline	50 000	OPIATES	MOP
Methoxyphenamine	6 000	Codeine	200
D-Amphetamine	1 000	Ethylmorphine	6 000
BENZODIAZEPINES	BZO	Hydrocodone	50 000
Alprazolam	100	Hydromorphone	3 000
a-hydroxyvalprazolam	1 500	Levorphanol	1 500
Bromazepam	900	β-Monoacetylmorphine	300
Chlordiazepoxide	900	Morphine 3-β-D-glucuronide	800
Clobazam	200	Morphine	300
Clonazepam	500	Norcodeine	6 000
Clorazepate dipotassium	500	Normorphone	50 000
Delorazepam	900	Oxycodone	30 000
Desalkylflurazepam	200	Oxymorphone	50,000
Diazepam	300	Procaine	15,000
Estazolam	6 000	Thebaine	6,000
Flunitrazepam	200	MARIJUANA	THC
(±) Lorazepam	3 000	Cannabinol	35 000
RS-Lorazepam glucuronide	200	11-nor-Δ ⁸ -THC-9 COOH	50
Midazolam	6 000	11-nor-Δ⁸-THC-9 COOH	50
Nitrazepam	200	Δ ⁸ -THC	17 000
Norchlordiazepoxide	100	Δ ⁹ -THC	17 000
Nordiazepam	900	METHAMPHETAMINES	MET
Oxazepam	300	p-Hydroxymethamphetamine	25 000
Temazepam	100	D-Methamphetamine	1 000
Triazolam	3 000	L-Methamphetamine	20 000
		(±)-3,4-Methylenedioxymethamphetamine	12 500
		Mephentermine	50 000

CROSS-REACTIVITY


A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or drug-positive urine. The following compounds did not show a cross-reactivity when tested with the DIAQUICK Multi-6 Cup (urine) at a concentration of 100 µg/mL.

Non Cross-Reacting Compounds:

Acetophenetidin	Cortisone	Zomepirac	d-Pseudoephedrine
N-Acetylprocainamide	Creatinine	Ketoprofen	Quinidine
Acetylsalicylic acid	Deoxycorticosterone	Labeltalol	Quinine
Aminopyrine	Dextromethorphan	Loperamide	Salicylic acid
Amoxicillin	Diclofenac	Meprobamate	Serotonin
Ampicillin	Diffunisal	Methoxyphenamine	Sulfamethazine
l-Ascorbic acid	Digoxin	Methylphenidate	Sulindac
Apomorphine	Diphenhydramine	Nalidixic acid	Tetracycline
Aspartame	Ethyl-p-aminobenzoate	Naproxen	Tetrahydrocortisone, 3-acetate
Atropine	β-Estradiol	Niacinamide	Tetrahydrocortisone
Benzoic acid	Estrone-3-sulfate	Nifedipine	Tetrahydrozoline
Benzoic acid	Erythromycin	Norethindrone	Thiamine
Bilirubin	Fenoprofen	Noscapine	Thioridazine
d,l-Brompheniramine	Furosemide	d,l-Octopamine	d,l-Tyrosine
Caffeine	Gentisic acid	Oxalic acid	Tolbutamide
Cannabidiol	Hemoglobin	Oxolinic acid	Triamterene
Chloral hydrate	Hydralazine	Oxymetazoline	Trifluoperazine
Chloramphenicol	Hydrochlorothiazide	Papaverine	Trimethoprim
Chlorothiazide	Hydrocortisone	Penicillin-G	d,l-Tryptophan
d,l-Chlorpheniramine	o-Hydroxyhippuric acid	Perphenazine	Uric acid
Chlorpromazine	3-Hydroxytyramine	Phenelzine	Verapamil
Cholesterol	d,l-Isoproterenol	Prednisone	
Clonidine	Isoxsuprine	d,l-Propranolol	

REFERENCES

- Baselt, RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488
- Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company. 1986; 1735
- Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986.

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