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

- Oesterling J E. Prostate specific antigen: a critical assessment of the most useful tumor marker for adenocarcinoma of the prostate[J]. The Journal of Urology, 1991, 145:907-923.
- Balk S P, Ko Y J, Bubley G J. Prostate-Specific Antigen[J]. Journal of Clinical Oncology, 2003, 21:383-391.
- Mccormack R T, Wang T J, Ritterhouse H G, et al. Molecular forms of prostate-specific antigen and the human kallikrein gene family: a new era[J]. Urology, 1995, 45: 729-744.
- Benson M C, Wang I S, Olsson C A, et al. The use of prostate specific antigen density to enhance the predictive value of intermediate levels of serum prostate specific antigen[J]. The Journal of Urology, 1992, 147:817-821.
- Becker C, Noldus J, Diamandis E, et al. The role of molecular forms of prostate-specific antigen (PSA or hK3) and of human glandular kallikrein 2 (hK2) in the diagnosis and monitoring of prostate cancer and in extra-prostatic diseases[J]. Critical Reviews in Clinical Laboratory Sciences, 2001, 38(5):357–399.
- Abrahamson P A, Lijja H, Oesterling J E. Molecular forms of serum prostate-specific antigen: the clinical value of percent free prostate-specific antigen[J]. Urologic Clinics of North America, 1997, 24: 353-365.
- Catalona W J, Smith D S, Ratliff T L, et al. Detection of organ-confined prostate cancer is increased through prostate-specific based Screening[J]. JAMA, 1993, 270:948-954.
- Oesterling J E, Jacobsen S J, Chute C G, et al. Serum prostate-specific antigen in a community-based population of healthy men: establishment of age-specific reference ranges[J]. JAMA, 1993, 270:860-864.
- Carter H B, Pearson J D, Metter J, et al. Longitudinal evaluation of prostate-specific antigen levels in men with and without prostate disease[J]. JAMA, 1992, 267:2215-2220.
- Berry S J, Coffey D S, Walsh P C, et al. The development of human benign prostatic hyperplasia with age [J]. The Journal of Urology, 2000, 132:474–479.
- Richie J E, Ratliff T L, Catalona W J, et al. Effect of patient age on early detection of prostate cancer with serum prostate-specific antigen and digital rectal examination[J]. Urology, 1993, 42(4):365-374.
- Sturgeon C M, Duffy M J, Stenman U H, et al. National academy of clinical biochemistry laboratory medicine practice guidelines for use of tumor markers in testicular, prostate, colorectal, breast, and ovarian cancers [J]. Clinical Chemistry, 2008, 54 (12):e11–e79.
- CLSI. Statistical Quality Control for Quantitative Measurement Procedures: Principles and Definitions. 4th ed. CLSI guideline C24. Wayne, PA: Clinical and Laboratory Standards Institute; 2016.
- Collins G N, Martin P J, Wynn-Davies A, et al. The effect of digital rectal examination, flexible cystoscopy and prostatic biopsy on free and total prostate specific antigen, and the free-to-total prostate specific antigen ratio in clinical practice[J]. The Journal of Urology, 1997, 157(5):1744-1747.
- Yuan J J, Coplen D E, Petros J A, et al. Effects of rectal examination, prostatic massage, ultrasonography and needle biopsy on serum prostate specific antigen levels[J]. Journal of Urology, 1992, 147(3 Part 2):810-814.
- Tchetgen M B, Song J T, Strawderman M, et al. Ejaculation increases the serum prostate-specific antigen concentration[J]. Urology, 1996, 47(4):511-516.
- Morgan W R, Zincke H, Rainwater L M, et al. Prostate specific antigen values after radical retropubic prostatectomy for adenocarcinoma of the prostate: impact of adjuvant treatment (hormonal and radiation)[J]. The Journal of Urology, 1991, 145(2):319-323.
- Robert W, Schroff, Kenneth A, Foon, Shannon M, Beatty, et al. Human Anti-Murine Immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy [J]. Cancer Research, 1985, 45(2):879-885.
- Primus F J, Kelley E A, Hansen H J, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy [J]. Clinical Chemistry, 1988, 34(2):261-264.
- Boscolo L M, Stuart M C. Heterophilic antibodies: a problem for all immunoassays [J]. Clinical Chemistry, 1988,34(1):27-33.

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**Shanghai International Holding Corp. GmbH (Europe)**